

Product datasheet for **MG208883**

Asns (NM_012055) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Asns (NM_012055) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Asns
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MG208883 representing NM_012055
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGTGTGGCATCTGGCCCTCTTCGGCAGCGATGACTGCCTTTCCGTGCAGTGTCTGAGTGCATGAAGA
 TCGCGCACAGGGGCCAGATGCATTTTCGCTTTGAGAAATGTCAATGGATACACCAACTGCTGCTTTGGCTT
 TCACCGCTTGGCTGTGGTTGACCCGCTGTTTGAATGCAGCCGATAAGAGTGAGGAAATACCCCTATTTG
 TGGCTCTGTTACAATGGTAAATCTACAACCACAAGGCGCTACAGCAACGTTTTGAATTTGAATATCAGA
 CCAATGTGGATGGTGAATTATCCTCCACCTCTATGACAAAGGAGGCATCGAGAAAACCTCTGTATGCT
 GGACGGGTGTTTGCATTCATCTACTGGACACTGCCAATAAGAAAGTATTTCTGGGCAGAGACACCTAT
 GGAGTCAGGCCCTGTTTAAAGCCATGACAGAAGATGGGTTTCTGGCTGTGTGTTCAGAAGCTAAAGGCC
 TTGTTTCCTTAAACACTCCACCACTCCCTTCTTAAAAGTGGAGCCCTTCTTCTGGACACTATGAAGT
 TTTGGATTTAAACCAAATGGCAAAGTTGCGTCTGTGAAATGGTCAAATACCATCACTGTACGGATGAA
 CCATTGCATGCCATCTATGACAGCGTGGAGAACTCTTCCAGGCTTTGACCTAGAGACCGTGAAGAACA
 ATCTGCGTATCCTTTTTGACAACGCTATCAAGAAACGCTTGATGACAGACCGGAGGATTGGCTGCCTTTT
 ATCAGGGGGCCTGGACTCGAGCTTGGTTGCTGCCTCTCTGCTGAAGCAACTCAAGGAGGCCCAAGTTCAG
 TATCCTCTCCAGACATTTGCTATTGGCATGGAGGACAGCCCCGATCTCCTGGCCGCTAGAAAAGTGGCAA
 ATTATATTGGAAGCGAGCATCATGAAGTCTTTTTAACTCTGAAGAAGGCATTACGGCCCTGGATGAAGT
 CATATTTTCCTTGGAACTTATGATATTACGACAGTTCGGGCATCTGTGGGCATGATTTAATTTCCAAG
 TATATTCGGAAGAACACAGACAGCGTGGTATCTTCTCCGAGAGGGGTGAGATGAACTTACACAGGGCT
 ATATATTTCCACAAGGCTCCTTCCCTGAGAAGGCCGAGGAGAGAGTGAAGACTGCTGAAGGAACT
 CTACCTGTTTGTGTTCTCCGGGCCGACCGCACTACTGCTGCACACGGTCTCGAACTGAGAGTCCCTTTT
 CTGGATCATCGGTTTTCTTCTATTACCTGTCTCTGCCGCCAGATATGAGAATTCAAAAAATGGCATAG
 AAAACATCTCCTGAGAGAGACTTTTGAGGACTGCAACCTGCTACCCAAAGAGATTCTCTGGCGACCCAA
 AGAAGCCTTCAGTGATGGGATCACCTCAGTCAAGAAGTCTGTTCAAGATTTTGCAGGACTATGTTGAA
 CATCAGGTTGATGATGAGATGATGTCTGCAGCCTCCAGAAGTTTCCCTTCAATACTCCAAAACTAAGG
 AAGGCTACTTCTACCGTCAGATCTTGAACGCCATTACCCAGGCCGGCTGATTGGCTGACTCATTATTG
 GATGCCTAAGTGGATCAATGCTACTGACCCTTCTGCCCGCACTCTGACCCATTATAAGTCAGCTGCCAAA
 GCT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence:

>MG208883 representing NM_012055
 Red=Cloning site Green=Tags(s)

MCGIWAIFGSDCLSVQCLSAMKIAHRGPDAFRFENVNGYTNCCFGFHRLAVVDPLFGMQPIRVKYPYL
 WLCYNGEIYNHKAQQRFEFYQTNVDGEIILHLYDKGGIEKTIKMLDGVFAFILLDTANKKVFLGRDTY
 GVRPLFKAMTEDGFLAVCSEAKGLVSLKHSTTFFLKVPEFLPGHYEVLDLKPNGKVASVEMVKYHHCTDE
 PLHAIYDSVEKLFPGFDLETVKNNLRILFDNAIKKRLMDDRRIGCLLGGLLDSSLVAASLLKQLKEAQQV
 YPLQTFAIMEDSPDLLAARKVANYIGSEHHEVLFNSEEGIQALDEVIFSLETYDITTVRASVGMYLISK
 YIRKNTDSVVIFSGEGSDELQGYIYFHKAPSPEKAEESERLLKELYLFDVLRADRTTAAHGLELRVFP
 LDHRFSSYYLSLPPDMRIPKNGIEKHLRETFCNLLPKEILWRPKEAFSDGITSVKNWFKILQDYVE
 HQVDDEMSAASQKFPFNTPKTKEGYFYRQIFERHYPGRADWLTHYWMPKWINATDPSARTLTHYKSAK
 A


TRTRPLE - GFP Tag - V

Restriction Sites:


Sgfl-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:



ACCN:

NM_012055

ORF Size:	1683 bp
OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_012055.4
RefSeq Size:	1928 bp
RefSeq ORF:	1686 bp
Locus ID:	27053
UniProt ID:	Q61024
Cytogenetics:	6 A1