

## Product datasheet for **MG208963**

### Zfp768 (NM\_146202) Mouse Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Zfp768 (NM_146202) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Zfp768
Synonyms:	BC026432; Znf768
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide Sequence:**

>MG208963 representing NM\_146202  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGGAGCGGGAGCGCTCGTCGTGGGGCCTGGAGTCCCAGGATGTGCACAGTCCCAACGCAGTGGGCAGCC  
 CCGAAGGGTCTCTCAAAGATCCTGCAGGCAACACAAGTGAGAATGAAGAAGGGGAAATTTCTCAGCGAGA  
 AGGCAATGGGACTATGAAGTTGAAGAAATACCTTTTGGGCTTGAACCCCAAAGCCAGAGTTTGAAGCCA  
 CAAAGCCAGAATTTGAATCCCAGAGCCCCAGGTTTGGAGCCTGAAAGCCAGGGTTTGAATCCCAGAGCC  
 CTGGGTTTGTGCCCCAAGCCCTGAGTTTGCACCCAGAAGTCTGAATCAGATCCTCAGAGCCCGGAGTT  
 TGAATCCCAGAGTCTAAGTATGAGCCCCGAAGCCCTGGCTGCCATCCCAGGAGCCCTGGCTGTGAACCT  
 GGGAGTCTAGGTATGAACCCAAAAGCCCTGGGTACGGTTCAAAGAGCCCCGAATTTGAATCTCAAAGTC  
 CAGGGTATGAATCTCAGAGCCCCGGGTATGAACCCGAGAAGTCTGGAGATGGCGTTCAGAATTCTGAGTT  
 CAAAACCCATAGCCCTGAGTTTGAAGTCAAAGTTCAAATTCAGGAAGGTGCAGAGATGCCTCTGAGC  
 CCCGAGGAAAAGAATCCCTTGAAGTCTCCCTAGGAGTCCACCCCTGGATTCTTTACTCAAGGGTTTGG  
 GGGAGCAGCCACAGGAGCGCTGCCTCCATTTGACATGCCTTCGGGGGCTCTGCTGGCTGCACCCAGTT  
 TGAGATGCTCCAGAACTCACTGAATCTGACAGGAACCCTGCGTGGGCCAGGCCGTAGGGGTGGCCGGCC  
 AGGGGTGGGAAGGCCCTCGACCTAACATCTGTGGCATCTGCGGGAAGAGCTTTGGACGAGGCTCGACCC  
 TGATCCAGCACCAGCGCATCCATACCGCGGAGAAGCCTTACAAATGTGAGGTTTGCAGCAAGGCCTTTTC  
 CCAGAGCTCTGACCTCATCAAACACCAGCGCACGCATACAGGCGAGCGGCCCTACAAGTGTCCAGGTGC  
 GGCAAGGCCTTCGCTGATAGCTTTACCTACTCCGCCACCAGCGTACCCACTCTGGCCAGAAGCCCTACA  
 AATGTCCACACTGTGGGAAGGCCTTTGGCGACAGCTCTACCTCTGCGACACCAGCGCACCCACAGCCAA  
 CGAGAGGCCCTACAGCTGCCCGAGTGTGGCAAGTGTATAGCCAGAAGTCTCCCTGCGCAGTCAACCAA  
 AGGGTGCACACTGGGCAACGGCCTTTGAGTGTGGCATCTGTGGCAAAAGCTTCTCCCAACGGTCAAGC  
 TTATCCCCACGCCCGCAGCCACGCCCGTAAAAGCCCTTCAAGTGCCAGAGTGCAGGCAAGCGCTTCGG  
 CCAGAGCTCAGTTCTGGCCATCCATGCCCGCACCCACCTGCCAGGCCGACGTACAGTGCCTCCGATTGT  
 GGCAAGACCTTCAACCGCTCCTTACACTATTACAGCACCAGCGTTCCACACAGGCGAGAGGCCGTACC  
 GCTGCGTGTGTGGCAAGGGCTTCTGTCGCTCCTCCACGCTGTACAGCACCATCGCGTGCATAGTGG  
 GGAGCGACCTTAAGTGTGATGACTGTGGAAGGCCTTCTCGCAGAGCTCCGACCTCATCCGCCACCAG  
 CGGACCCACGAGCCGGCCCGCT

**ACGCGT**ACGCGGCCGCTCGAG – GFP Tag – GTTTAA

**Protein Sequence:**

>MG208963 representing NM\_146202  
 Red=Cloning site Green=Tags(s)


MEREASSWGLSRDVHSPNAVGSPEGSLKDPAGNTSENEEGEISQREGNGDYEVVEEIPFGLPEQSPFEFEP  
 QSPEFESQSPRFEPESPGFESRSPGFVPPSPFEFAPRSPESDPQSPEFESQSPKYEPRSPGCHPRSPGCEP  
 GSPRYEPKSPGYGSKSPEFESQSPGYESQSPGYEPQNSGDGVQNSEFKTHSPEFETQSSKFQEGAEMPLS  
 PEEKNPLSISLGVHPLDSFTQGFGEQPTGALPPFDMPSGALLAAPQFEMLQNLNLTGLTRGPGRRGGRA  
 RGGQGPRPNICGICGKSFGRGSTLIQHQRHTGEKPYKCEVCSKAFSSDLIKHQRHTGERPYKCPRC  
 GKAFADSSYLLRHQRTHSGQKPYKCPHCGKAFGDSSYLLRHQRTHSHERPYSCPECGKCYSQNSSLRSHQ  
 RVHTGQRPFSQICGKSFQRSALIPHARSHAREKPFKCECGKRFQSSVLAIHARTHLPGRTYSPDC  
 GKTFNRSSTLIQHQRSHTERPYRCAVCGKGFCSSTLLQHHRVHSGERPYKDDCGKAFSSDLIRHQ  
 RTHAAGR

**TRTRPLE** – GFP Tag – V


**Restriction Sites:**

SgfI-MluI

## Cloning Scheme:



## Plasmid Map:



ACCN: NM\_146202

ORF Size: 1704 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_146202.1</a> , <a href="#">NP_666314.1</a>
<b>RefSeq Size:</b>	2268 bp
<b>RefSeq ORF:</b>	1707 bp
<b>Locus ID:</b>	233890
<b>UniProt ID:</b>	<a href="#">Q8R0T2</a>
<b>Cytogenetics:</b>	7 F3
<b>Gene Summary:</b>	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]