

Product datasheet for **MG209000**

Zfp128 (NM_153802) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Zfp128 (NM_153802) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Zfp128
Synonyms:	9630016P15Rik; mZnf8; Znf8
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MG209000 representing NM_153802
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGACCACCAGGACAAAGCAGCAACAGTGGCAATGGCCTCGCGCCGCAAGCAACAACACTTCAGGAGC
 CAGTGACTTTCCGAGATGTGGCTGTGGACTTTACCCAAGAGGAGTGGGGGACAGCTGGACCTACCCAGAG
 GACTCTGTATCGTGATGTGATGCTGGAGACCTTTGGGCACCTTCTCTGTGGGTCTGATCTTCCCAAA
 CCTGCAGTTATCTCCAGCTGGAACAAGGTGCTGAGCTGTGGTGGCAGACAGAGGAGGCACCGGGGCT
 GCCATCCAGGCTGGATTCTGAACCTGAAGACCACACATTGCTCAAGGACCAGGGCTCCCAAAAATGGA
 ACCATCCCCTATTACAGAAAAGGATGGTTTGC AAAAGGCTGTTCCCTGTCGTTCTATGATAGGGATAGAT
 CAGGAAAGTGATGGCCAGAGGCAGGCACTCAAGAAGGACCAGAGTAACCTGAATGACCCAAAGAGATAC
 CACTTCAAAGTCAAAGCCACAAAAGTCTTGGACTTGTGGAAGCCTGTGTTCTTGGCTTAAACACTTATCT
 ATTACCAGATATTTCTGGAAGAGAGTATGGCTGTACTTATGACTCACAAGTTAAAAATTCAGAACATAAT
 CCAAGCTTAGTTAGACAACGGACAGACTCCCCTGCAACACAGTCCTTTGATGACAATGGCAGTCAAAAAG
 CCTTTGATCAAATTATGCCCATACAGAACTCACAAAAAGTCAGGTGCAAGATAAACCTATAAGTGTAC
 TGATTGTGGGAAGTCATTTAACCATAACGCACATCTCACAGTGCACAAGAGAATTCATACGGGAGAAAGA
 CCTTACATGTGCAAAGAGTGTGGGAAAGCCTTCAGCCAGAACTCCTCCCTGTTCAACATGAGCGAATCC
 ACACTGGAGACAAGCCTTATAAGTGTGATGAATGTGGGAAATCTTTTGCACAGTACACACCTTACAGT
 CCATAGGAGAATTCACACTGGAGAGAAACCCTATGAGTGTGAGGACTGTGGGAGGCTTCAATCAGAAT
 TCATCCCTAGGTGCACACAAGAGGACACACACTGGGAGAAAGCCTTATACCTGCAGTGTGTTGGGAAAT
 CTTTCTCCGGACCCTTGCCTCTTCTGCACCTGCGGACTCACACTGAGGAAAGCCATACGAGTGATAA
 TCACTGTGGAAAAGGCTTCAGGCACAGTTCCTCCCTGGCCAGCATCAGCGGAAGCACGCTGGTGAGAAA
 CCCTATGAGTGCCGACAGAGGCTGATCTTTGAGCAGGCGCCGGCTTAAATAAAGTACGAATGGACAGAAC
 CCCTCGGCTGTGACTCACCGTTGAGTCAAGGTGAAAGGACTCAACGGAGTGATAGGCCCTTAAAGTGTA
 TCAGTGTGAAAAGTGTTCCTCAGAGTCTCACCTCATCCGACATCAGCTCACCCACAGCAGAGAGGAA
 GAACCACTACGTGGACGTAGCAGAAGGCAAGAACGCCGTGTAGACGGGCTCAAGACTCATTGAGAATA
 CAAACTCAAATTCAGGGAGCTTCTGTTGCCAGCCAAAGGCAGGACAAGCCAGCAGAACCTGGCCCT
 GTTTGACCTTCGTGAAATATGCAGGAGCAAAACCCTGTGCATGTTATTGGGGTAGAAGAACCCTCAGTG
 GGCAACTCCATGCTGTTTACACACAGAGAATCTAGA

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence:

>MG209000 representing NM_153802
 Red=Cloning site Green=Tags(s)


MDHQDKAATVAMASRPQATQLQEPVTFRDVAVDFTQEEWGQLDPTQRTL YRDVMLETFGHLLSVGPDLPK
 PAVISQLEQGAELWYADRGGTGACHPGWILEPEDHTLLKDQGLPKMEPSPITEKDFAKAVPCRSMIGID
 QESDQQRQALKKQSNLNDPKEIPLQSQSHKSLGLVEACVLGLNTYLLPDISGREYGYCTYDSQVKNSEHN
 PSLVLRQRTDSPATQSFDDNGSQKAFDQIMPITELTKSQVQDKPYKCTDCGKSFHNHAHLTVHKRIHTGER
 PYMCKEKGKAFSQNSLVQHERIHTGDKPYKCECGKSFCHSTHLTVHRRIHTGEKPYECQDCGRSFNQN
 SSLGRHKRTHTEKPYTCSVCGKSFRTTCLFLHLRTHTEERPPECNHCCKGFRHSSSLAQHQKHAGEK
 PYECRQRLIFEQAPALIKYEWTEPLGCDSPLSQGERTQRSRPFKCNQCGKCFQTQSSHLIRHQLTHSREE
 EPLRGRSRRQEQPCRRGSRLIQNTNSNSREL PVAQPKAGQASRTLALFDLREIMQEONPVHVIQVEEPSV
 GNSMLFDTRESR

TRTRPLE – GFP Tag – V


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_153802
 ORF Size: 1716 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_153802.3 , NP_722497.1
RefSeq Size:	3850 bp
RefSeq ORF:	1719 bp
Locus ID:	243833
UniProt ID:	Q8BGV5
Cytogenetics:	7 A1
Gene Summary:	Transcriptional repressor. May modulate BMP and TGF-beta signal transduction, through its interaction with SMAD proteins.[UniProtKB/Swiss-Prot Function]