

Product datasheet for **MG209244**

HnrpII (BC012849) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	HnrpII (BC012849) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	HnrpII
Synonyms:	PHGPx
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MG209244 representing BC012849
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGCATCGCC

ATGTCCTCCACCTCCTCCTCCTCCCCAAGGAGGAGACGTACGAGGAGGACCGGGAGTTCGAGAGCCAGG
 CCAAGCGCCTCAAGACCGAGGAGGAGAGATCGTCTACTCGGCCGAGGAGAGCGAGAATCGCCAGGAAGC
 GACGCCTCAGGCGGGGAGCGATAGCGACAGCGCGCGCGCGACGAGGCGACGGAGACGGAGGCAGCGGC
 GGCGGCGCGACGGCGAGGAAGGCGAAGGTGGAGAAGAAGGCGACGAAGGCGACGGCGATGAGGGAGGGA
 GCGGCGGCGACGAAGGCGGACGGCGGTGGCGGCGCGGCCCTCGGAGCATGCCCTGTCCACGGAGGG
 AGGTGGAAGTCATCATAAAGTTTCTGTTCCCTGTTGTCCATGTTGAGGGCTCTGTGAATCTGTTGTG
 GAAGCTGACCTTGTGGAGGCACTGGAGAAATTTGGGACAATATGCTATGTGATGATGATGCCATTTAAAC
 GACAGGCTCTCGTGGAGTTTAAAAATAGATAGTCCAAAGAATGTGTGACATTTGCTGCAGATGTGCC
 TGTGTACATTGCTGGCCAGCAGGCTTTCTCACTATTCTACAAGTAAAGAATCACTCGGCCAGGAAAT
 ACCGACGATCCATCAGGGGGCAACAAAGTCTCTTGTCTCTATTAGAAATCCTCTTTACCAATTACGG
 TGGATGTTCTGTATACTGTATGCAACCTGTTGGAAAAGTACAGCGTATTGTTATCTTCAAGAGAAATGG
 GATACAAGCAATGGTTGAGTTTGTGAGTCAGTCTTTGCGCCAGAAAGGCTAAAGCAGCACTCAATGGAGCT
 GATATATATGCTGGATGTTGCACACTAAAAATGAATATGCAAGGCCAACTCGTCTAAATGTTATTAGGA
 ATGACAATGACAGTTGGGACTACACTAAACCTTATTTGGGAAGACGAGATAGAGGAAAGGGTCGCCAGAG
 ACAAGCCATTTTGGGAGACCCTCTTCTGTTTACAGATGATGGCTATGGATCGCATGGTCCGCTACTG
 CCTTTACCAAGTCGTTATAGAATGGGCTCTAGAGATACACCTGAACCTGTTGCTTATCCGTTACCACAAG
 CGTCTTCTCTTACATGCATGGAGGCACTCCGCTGGCTCCGTTGTAATGGTTAGTGGATTACATCAGCT
 CAAAATGAATTGTTCAAGAGTCTTCAACCTGTTCTGCTTGTATGGAAATATTGAAAAGGTAAAGTTTATG
 AAGACCATTCTGGCACAGCACTTGTAGAAATGGGTGACGAGTATGCTGTGGAAAGAGCTGTCAACCCACC
 TTAATAATGTCAAACCTTTTGGGAAAAGACTTAATGTTTGTGTCTAAACAGCACTCAGTTGTTCCAG
 TCAAATATTTGAGCTAGAGGATGGTACGAGCAGCTACAAAGATTTTGAATGAGCAAAAATAATCGCTTT
 ACAAGTGTGGACAAGCATCTAAGAATATAATCCAGCCACCGTCATGTGTGCTGCATTATTACAATGTTT
 CGCTGTGTGTACAGAAGAGACCTTCAAAAGTTGTGTAATGACCACGAAGTTCTCCGTTCAACAAATA
 TAAAGTGTGTGACGCGAAAGCTTCTGCGAAAACCTGTCTGGGCTATTGGAGTGGAAAGTCAAAACCTGAT
 GCCGTGGAAGCTCTCACTGCGCTCAATCACTACCAGATACGAGTCCCGAATGGCTCAACCCCTATACGT
 TGAAGCTTTGCTTTTCTACATCATCCCATTTA

ACGCGTACGCGGCGGCTCGAG – GFP Tag – GTTTAA

Protein Sequence:

>MG209244 representing BC012849
 Red=Cloning site Green=Tags(s)


MSSTSSSSPKEETYEEDREFESQAKRLKTEEGEIVYSAESENQAEATPQAGSDSDSGGDDGGDGGSG
 GGGDGEEGEGGEEGDEGDGDEGGSGGDEGGSGGGGGPRSMPLSTEGGSHHKVSVSPVVHVRGLCESV
 EADLVEALEKFGTICYMMMPFKRQALVEFENIDSAKECVTFADVPVYIAGQQAFFNYSTSKRITRPGN
 TDDPSGGNKVLLLSIQNPLYPIITVDVLYTVCNPVGKVQRIVIFKRNGIQAMVEFESVLCQAQAKAALNGA
 DIYAGCCTLKIEYARPTRLNVIRNDNDSDWYTKPYLGRDRGKGRQRQAILGDHPSSFRHDGYGSHGPLL
 PLPSRYRMGSRDTPELVAYPLQASSSYMHGGSPSGSVVMVSGLHQLKMNCSTVFNLFCLYGNIEKVKFM
 KTIPGTALVEMGDEYAVERAETHLNNVKLFGKRLNVCVSKQHSVVPSTQIFELDGTSSYKDFAMSKNNRF
 TSAGQASKNIIQPPSCVLHYNNVPLCVTEETFTKLCNDHEVLPFIKYKVFDAKASAKTSLGLEWKCKTD
 AVEALTALNHYQIRVPNGSNPYTLKLCFSTSSHL

TRTRPLE – GFP Tag – V

Restriction Sites:

SgfI-MluI

Cloning Scheme:




ACCN:	BC012849
ORF Size:	1782 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	BC012849 , AAH12849
RefSeq Size:	3017 bp
RefSeq ORF:	1784 bp
Locus ID:	72692

Cytogenetics: 17 E3

Gene Summary: RNA-binding protein that functions as regulator of alternative splicing for multiple target mRNAs, including PTPRC/CD45 and STAT5A. Required for alternative splicing of PTPRC. [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MG209244