

Product datasheet for **MG209543**

Slc13a4 (NM_172892) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Slc13a4 (NM_172892) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Slc13a4
Synonyms:	9630060C05Rik; SUT-1; SUT1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MG209543 representing NM_172892
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGCATCGCC

ATGGGCTTGCTGCAGGGCCTTCTCCAAGCGAGGAAGCTGCTATTGGTCATCTGTGTGCCGCTGCTGCTGC
 TGCTCTGCCCCACCATCTACCCACCAAGCAAGCTGCGTGTGCTTACGTGTTGCTTGTGACCGCTGTGTA
 CTGGGTGTCTGAGGCCGTGCCTCTCGGAGCTGCAGCCCTCGTGCCTGCCTTCTCTACCCGTTCTTCGGT
 GTTCTCCGATCCAGTGAGGTGGCGGCTGAGTACTTCAAGAACACCAAGCTGCTGTTGGTGGCGCTCATCT
 GTGTGGCAGCAGCTGTGGAGAAGTGAATTTACACAAGCGCATCGCCCTGCGCATGGTCTTGATGGCGGG
 AGCCAAGCCTGGCATGCTCCTGCTGTGTTTCATGTGCTGTACCACCATGCTGTCCATGTGGCTCTCCAAC
 ACCTCCACCACCGCCATGGTGATGCCATCGTGGAGGCGGTGCTGCAGGAGCTCATCAACGCTGAGGAGG
 AACAGCTGGCAGCCGCGACTGAGGAGGCTGAGCTCATGGGTCTTGATGTAACAACAGACAGACGTCCTAT
 GGAATCATCTTTGTCAATGAAGACACGTCGGCTGCAGACGTCCTTCTCTGATGCAGAGCAAGAATCTG
 AATGGTGTGCCATGGTACCAAATCCATCAACACAGCAAAACCAAGCAACAGAAGAAGCAGCAGCCGT
 CACAGGAAAAGCCAGGAGTCCCAACCCCGAGCTCCAAGACCCAAGAACTGAACAAGAAGAAGTACAGATC
 CCACCATGACCAAATGATCTGCAAGTGCCTTTCTCTGAGCATTTCTATGCTGCCACCATTTGGTGGCCTG
 ACCACCATCATAGGCACCTCCACCAGCCTCATTTTCTTAGAGCATTTCAACAACCAAGTACCCGGCGGCAG
 AGGTGGTGAATTCGGAACCTGGTTCCTCTTACGCTTCCCATCTCGTCTCATGCTGGTGGTCAAGTGT
 GTTCTGGATCCACTGGCTCTTCTGGGCTGCAACTTCAAAGAGACCTGCTCTCTGAGCAAGAAGAAGAAG
 ACCAGGAGGGAGGAGCTATCAGAGAAGAGGATCCGGGAGGAATATGAGAAGCTAGGAGACATTAGCTACC
 CTGAGATGGTGACCGCCTTTTCTTCTGATGACGGTGTGTTTACCCGGGAGCCCGCCTTTGT
 CCTGGCTGGGATTCTTTCTTTGAAAAGAAAGGCTACCGCACTGATGCCACAGTCTCTGTATTCCTTGGC
 TTCTGCTTTTCTCATTCCAGCAAGAAGCCCTGCTTTGGGAAAAAGAGTGATGAACAGGCCAGGAAG
 CCTCAAAGGGATAGAGCCCATCATCACGTGGAAGGACTTCCAGAAAATATGCCCTGGGAGATTGTCAT
 TCTGGCGGGAGGAGGCTACGCCCTGGCTTCTGGCAGCAAGAGCTCTGGCCTCTCCACATGGATCGGGCAC
 CAGATGCTGTCCCTGAGCAGCCTCCCGCATGGGCTATCACCTGCTGGCGTGTGCTCTGGTGTCCATTG
 TCACGGAGTTTGTGAGCAATCCAGCCACCATCACCATCTTCTGCCCATCCTTTGATCCCTGTCCGAAAC
 GCTGCACATCAATCCACTGTACACCTGGTCCCGTTACCATGTCAATCTCCTTCGCTGTGATGCTGCC
 GTGGGCAACCCCAACGCCATCGTCTTACGCTACGGGCACTGCCAGATCAAAGACATGGTGAAAGCCG
 GTCTGGGAGTCAATGTATCGGTCTGGTATCGTGATGGTGCCATCAATACTTGGGAGTTAGCCTCTT
 CCACTTGATGCTTTCCAGCCTGGGCTAAGGTACGAACATCACTGATCAGACC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence:

>MG209543 representing NM_172892
 Red=Cloning site Green=Tags(s)


MGLLQGLLQARKLLLVICVPLLLLPLPTIYPTSEAACAYVLLVTAVYVWSEAVPLGAAALVPAFLYPFFG
 VLRSEVAAEYFKNTTLLLVGVICVAAAVEKWNLHKRIALRMVLMAGAKPGMLLLCFMCCTTMSMWLSN
 TSTTAMVMPIVEAVLQELINAEELQAAGTEEAELMGLDVNNRQTSMEIFVNEDTSAADVTSMLMQSKNL
 NGVPMVTKSINTANQQQKKQPSQEKPGVPTPSSKTQELNKKKYRSHDQMICKCLSLISYAATIGGL
 TTIIGTSTSLIFLEHFNQYPAAEVNFGTWFLFSFPISLIMLVSWFWIHWLFLGCNFKETCSLSKSKK
 TRREELSEKRIREEYEKLDISYPEMVTAFFFILMTVLWFTREPGFVPGWDSFFEKKGYRTDATVSVFLG
 FLLFLIPAKKPCFGKSDGTGQEASKGIEPIITWKDFQKTMPEIVILAGGGYALASGKSSGLSTWIGH
 QMLSLSSLPWAITLLACVLVSIIVTEFVSNPATITIFLPILCTLSETLHINPLYTLVPVTMSISFAVMLP
 VGNPPNAIVFSYGHCIQKDMVKAGLGVNVIGLVIMVAINTWGVSLFHLDAFPAWAKVSNITDQT

TRTRPLE - GFP Tag - V

Restriction Sites:

Sgfl-MluI

Cloning Scheme:



ACCN: NM_172892

ORF Size: 1875 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_172892.1](#), [NP_766480.1](#)


RefSeq Size: 3426 bp

RefSeq ORF: 1878 bp

Locus ID: 243755

Cytogenetics: 6 B1

Product images:



Circular map for MG209543