

Product datasheet for **MG209573**

Ttc27 (BC021912) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ttc27 (BC021912) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Ttc27
Synonyms:	2610511O17Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MG209573 representing BC021912
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGAAACAGGAGAATTTATTTGAAGGTGATTAGGTCGACTTTTGGCTATTCAATTTTCATCTGGAATGTG
 CACATGTATTTTTATATTATTATGAATATAAAGAAGCAAAGGACCAGTTCAGTACTGCTAAAGACATCAG
 CAAGTTAGAGATTGACCTGACAGGTGCCTTGGGAAAAAGGACACGATTCCAGGAAAATTACGTGGCGCAG
 CTGATTGTAGATGTCAGGAGAAAAGAGGCTGTCCCTTTTCAGTTGTGAGTTCAGCCAGCACCCACGCCCC
 AGGAGTGTGGCCAAGAATCTTGAGCTCAATGATGATACTGTTCTGAATGAAATAAAGTTAGCAGATAG
 TGAACGGTTTCAGATGCCAGACCTGTGTGCGGAGGAGCTTGCCGTGGTCTTGGAGTCTGCACAAATTTTC
 CAGAAGAACAACCCAGTTCATAAGCTAACGGAGGAGGAGCTGCTGGCCTTCACATCGTGTCTTACTTTTAC
 AACCAAAGTTCTGGGCCATTGATGTCAGCCTTGATCCTCCGACAAAACCTTGAGAGAGGAAGCACACG
 GCGAGTGGAGCGGCCATGAGGCAGACCCAGGCTCTTGACAGCAATTCGAAGATAAAGCCACATCTGTA
 CTGGAACGTCTGAAGATTTTCTATTGTTGTCAAGTACCCCTCACTGGGCTGTTACGCGCCAGCTCGCAG
 GTTTGCTGTTTGAAGTTGGGATGCACCTCCTCAGCCCTGCAGATATTTGAGAAGCTGGAGATGTGGGAGGA
 TGTGGTCACTGCCCACGAGAGGGCTGGGCGGCACGGGAAGGCAGAAGAACTACTAGGCAAGAACTGGAG
 AAAAAGGAAACACCAGGTTTATATTGCTTACTTGGAGATGTCTCGAGGACCATTCTGTACGACAAGG
 CCTGGGAGCTGTCCCGCATCGCAGTGCACGTGCCAGCGCTCCAAAGCTCTCTCCATCTTCGGAACAA
 GGAGTTTCGGGAATGTGTGGAGTGCTTCGAGCGCTCAGTGAAGATCAATCCCATGCAGCTCGGGGTCTGG
 TTCTCTCTAGGCTGCGCCTACTTGGCCTTGAAGACTATGGTGGATCAGCAAAGGCGTTTCAGCGTTGTG
 TGACCCTAGAACCCGATAATGCTGAGGCTTGAACAATTTATCTACTTCTATATTCGATTAAAGCAAAA
 AGTCAAAGCTTTTGAACGTTACAAGAAGCCCTCAAGTGTAACACGAACACTGGCAAATCTGGGAAAAAC
 TACATTCTCACCAGCACCAGTGTGCGGGAGTTTGGGAAGCCATCAAAGCTTACCACCGGCTCTGGACC
 TCCGAGACAAATACAAAGATATTCAGGTCCTTAAATTTCTGGTCCAGGCAGTGGTTAATGATATGACCGA
 CCGAAGTGGAGCTGTGCTCCAGCCTCAAAGGCAAGCTGCAGGAGTTGTTTGGCAGAATAACTTCAAGG
 GTGACCAATGACGGAGAAGTCTGGAGACTGTATGCCCAAGTGCATGGGAATGGGCAGAGTGAAGGCCAG
 ATGAGAACGACAAGGCATTCCAGTGTCTCTGAAGGCGTACAAGTGTGACACACAGTCCAGTTGTTGGGA
 GAAAGATGCTACAGCATTAAAGGAAGTTGTTCAAAGAGCCATAGGACTCGCTCATGTGGCTATGAAATGC
 GCTGAGAGCAAGTCCATTCCCAGGAAGCTGTACAAACACTTTTTCAGTTCGACTCAATTTACGTGGCT
 TGTATCCAAAGCAAAGAATTTACAGATGTGGTAAGTGGAGAAGTCTCTGGAGAGCTCGCCAGTGAAAT
 AGCAGCTGTGAATGCTCTGGAGGCAGAGCTGCAAGACCTCAGCAACCAGCTTCGAAACCGGTAC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence:


>MG209573 representing BC021912
 Red=Cloning site Green=Tags(s)

MKQENLFEGDSGRLLAIQFHLECAHVFLYYEYKEAKDQFSTAKDISKLEIDL TGALGKTRFQENYVAQ
 LIYDVRRKEAVPFSCEFSAPTQEQECLAKNLELNDDTVLNEIKLADSERFQMPDLCAEELAVVLGVCTNF
 QKNNPVHKLTEEELLAFTSCLLSQPKFWAIQMSALILRTKLERGSTRRVERAMRQTALADQFEDKATSV
 LERLKIYFCCQVPPHWAVQRQLAGLLFELGCTSSALQIFEKLEMWEDVVICHERAGRHGKAEELRQELE
 KKETPGLYCCLLDVLDQHSCYDKAWELSRHRSARAQRSKALLHLRNKEFRECECFERSVKINPMQLGVW
 FSLGCAYLALEDYGGSAKAFQRCVTLEPDNAEAWNNLSTSYIRLKQKVKAFTLQEAALKCNIEHWQIWEN
 YILTSTDVGEFGAIAKAYHRLLDLTDKYKDIQVLKILVQAVVNDMTDRSGAVASSLKGLQELFGRITSR
 VTNDGEVWRLYAQVHNGQSEKPDENDKAFQCLSKAYKCDTQSSWEKDATAFKEVVQRAIGLAHVAMKC
 AESKSIPQEAQVTLSSVRLNLRGLLSKAKNFTDVSSEVSGELASEIAAVNALEAELQDLSNQLRNRY

TRTRPLE - GFP Tag - V

Restriction Sites:

SgfI-MluI

Cloning Scheme:


ACCN: BC021912

ORF Size: 1884 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [BC021912](#), [AAH21912](#)


RefSeq Size: 2688 bp

RefSeq ORF: 1886 bp

Locus ID: 74196

Cytogenetics: 17 E2

Product images:



Circular map for MG209573