

Product datasheet for **MG209689**

Fubp1 (BC023780) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Fubp1 (BC023780) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Fubp1
Synonyms:	9530027K12Rik; D3Ertd330e; FBP; Fubp; Fubp4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MG209689 representing BC023780
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGCCGACTACTCCACAGTGCCTCCGCCGTCGTCTGGCTCGGCTGGCGCGCGCGCGCGGAGTCGTTA
 ACGACGCTTTCAAAGATGCCCTGCAGAGAGCGCGCAGATTGCAGCAAAAATTGGGGGTGATGCTGGTAC
 ATCATTGAATTTAAATGACTATGGTTATGGGGGACAAAAAGACCCTTAGAAGATGGAGATCAGCCAGAT
 GCTAAGAAAGTACCTCCCCAAAATGACTCTTTGGAGCACAGTTACCTCCAATGCATCAGCAGCAAAGCA
 GATCTGTAAATGACAGAAGAATACAAAGTCCAGATGGGATGGTCGGATTATAATTGGCAGAGGAGGTGA
 ACAGATCTCACGAATACAGCAGGAATCTGGATGCAAGATACAGATAGCACCTGATAGTGGTGGCCTACCA
 GAAAGGTCTTGTATGCTAACTGGAACACCTGAATCTGTCCAATCAGCAAAAAGATTATTGGACCAGATTG
 TTGAAAAGGGAAGACCAGCCCTGGCTTTCATCATGGTGTGGACCTGGAATGCAGTTCAGGAAATCAT
 GATTCCAGCCAGCAAAGCAGGACTAGTTATTGAAAGGGGGCGAGACTATTAACAACCTCAGGAAACGG
 GCTGGTGTAAAAATGGTAAATGATTCAAGATGGGCCTCAAAACACTGGTGTGATAAACCTCTTAGGATTA
 CGGGTGACCCATACAAAGTTCAGCAAGCAAGGAAATGGTATTAGAGTTAATTCGTGATCAAGGTGGTTT
 CAGAGAAGTGGGAATGAGTATGGCTCAAGAATAGGAGGCAATGAAGGGATAGATGTCCCAATTCGAAGA
 TTTGCTGTTGGCATTGTAATAGGAAGAAATGGAGAAATGATTAACAAAATACAAAATGATGCTGGTGTTC
 GAATTCAGTTTAAAGCCAGATGATGGAACAACACCTGATAGGATAGCACAGATAACAGGACCTCCAGACAG
 GTGTCAGCACGCTGCAGAAATAACACAGACCTTCTACGAAGTGTTCAGGCTGGCAATCCTGGTGGACCG
 GGACCTGGTGGTCGAGGACGAGGTAGAGGTCAAGGAACTGGAATATGGGGCCCCGGGTGGACTCCAGG
 AGTTTAAATTTCAATGTGCCAACTGGGAAAACCTGGACTGATCATTGGAAAAGGAGGTGAAACCATAAAAA
 CATAAGCCAACAGTCTGGTGAAGAATAGAAGTGCAGAGAAGCCCTCCACCTAATGCAGATCCCAATATG
 AAGTTATTTACAATTCGGGGCACTCCACAGCAAAATAGACTATGCTCGACAACCTCATAGAAGAGAAGATTG
 GGGGCCAGTAAATCCTTTAGGGCCACCTGTACCCCATGGGCCCATGGGGTTCCAGGTCCTCATGGGCC
 TCCTGGACCTCCAGGGCCGGAACTCCAATGGGACCATAACAACCTGCACCTTACAATCCAGGACCACCT
 GGCCAGCTCCTCAGGTCCTCCAGCCCCATATGCTCCCAGGGATGGGAAATGCGTATCCACATTGGC
 AGCAACAGGCTCCTCCTGACCCAGCCAAGGCAGGAGCAGATCCAACTCGGCAGCTTGGGCTGCTTATTA
 TGCTCACTATTACCAGCAGCAGGACACACCCACCTGCAGCTCCTGCCGGTGCACCAGCTACAACCCAA
 ACGAACGGACAAGGAGATCAGCAGGCTCCAGCTCCAGCTGGACAGTTGATTATACAAGGCTTGGGAAG
 AATACTACAAGAAAATGGGTCAAGCAGTTCCTGCTCCTGCTGGGGCCCCACCAGTGGTCCAGCCGATTA
 TAGTGACGCTGGGCTGAGTACTATAGACAGCAAGCAGCGTATTATGCCAGACAAGTCCCAGGGGATG
 CCGCAGCATCCTCCAGCACCTCAGGGCCAA

ACGCGTACGCGGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence:

>MG209689 representing BC023780
 Red=Cloning site Green=Tags(s)


MADYSTVPPSSGSAGGGGGVVNDAFKDALQRARQIAAKIGGDAGTSLNSNDYGYGGQKRPLEDGDQPD
 AKKVPPQNSFQAQLPPMHQQSRSVMTEEYKVPDGMVGFIIIGRGEQISRIQQESGCKIQIAPDSGGLP
 ERSCLMTGTPESVQSAKRLLDQIVEKGRPAPGFHHGDGPGNAVQEIMIPASKAGLVIGKGETIKQLQER
 AGVKMVMIQDGPQNTGADKPLRITGDPYKVVQAKEMVLELIRDQGGFREVRNEYSGRIGGNEGIDVPIPR
 FAVGIVIGRNGEMIKKIQNDAGVRIQFKPDDGTPDRIAQITGPPDRCQHAAEIIITDLLRSVQAGNPGGP
 PGGRGRGRGQGNWMMGPPGGLQEFNFIVPTGKTGLIIGKGETIKSISQSGARIELQRSPPPNADPNM
 KLFTIRGTPQQIDYARQLIEEKIGGPVNPLGPPVPHGPHGVPVPHGPPGPPGPTMPGPNPAPYNP
 GPAPHGPPAPYAPQGWGNAYPHWQQAPPDPAKAGADPNSAAWAAYYAHYYQQQAQPPAAPAGAPATTQ
 TNGQGDQAPAPAGQVDYTKAWEEYKMGQAVPAPAGAPPGGQPDYSAAWAAYYRQQAAYYAQTSPPQGM
 PQHPPAPQGG

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



Plasmid Map:



ACCN: BC023780

ORF Size: 1920 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	BC023780.1
RefSeq Size:	2425 bp
RefSeq ORF:	1922 bp
Locus ID:	51886
Cytogenetics:	3 77.05 cM
Gene Summary:	Regulates MYC expression by binding to a single-stranded far-upstream element (FUSE) upstream of the MYC promoter. May act both as activator and repressor of transcription (By similarity).[UniProtKB/Swiss-Prot Function]