

Product datasheet for **MG209693**

Lrrc4c (NM_178725) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Lrrc4c (NM_178725) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Lrrc4c
Synonyms:	6430556C10Rik; NGL-1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MG209693 representing NM_178725
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGC**C

ATGTTGAACAAGATGACCTTACATCCACAGCAGATAATGATAGGTCCTAGGTTTAAACAGGGCCCTATTTG
 ACCCCCTGCTTGTGGTGCTGTTGGCTCTTCAACTTCTTGTGGTGCTGGTTTGGTTCGTGCTCAACCTG
 CCCTTCAGTGTGCTCTTGTAGCAACCAGTTCAGCAAGGTGATTTGTGTCGGAACAACTTCGTGAAGTT
 CCGGATGGCATCTCCACCAACACAAGGCTGCTGAACCTCCATGAGAACCAATCCAGATCATTAAGTGA
 ACAGCTTCAAGCATTTGAGGCACCTAGAAATCCTCCAGTTGAGCAGGAACCATATTCGAACCATGAAAT
 CGGGGCCCTTCAATGGTCTGGCAACCTCAACACCTGGAACCTTTGACAATCGTCTTACTACCATACCG
 AATGGAGCTTTTGTATATTTGTCTAACTGAAGGAGCTCTGGTTGAGGAACAACCTATTGAAAGCATCC
 CTTCTATGCTTTTAAACAGAATCCCTTCTTTCGCCCGCTGGACTTAGGGGAATTGAAAAGGCTTTCATA
 CATCTCAGAAGGTGCCTTTGAAGGTCTGTCCAACCTGAGGTATTTGAACCTTGCCATGTGAACCTTCGG
 GAAATCCCTAACCTCACACCACTCATCAAACTGACGAGCTAGATCTTTCTGGGAACCATTTGTCTGCAA
 TCAGGCCTGGCTCTTTTCAAGGGTTGATGCACCTTCAAAAACCTGTGGATGATACAGTCTCAGATTCAAGT
 GATTGAACGGAATGCCTTTGATAACCTTCAGTCACTAGTGGAGATAAACTTGGCACACAACATCTAACA
 TTATTGCCTCATGACCTCTTACACCCCTTGCATCATCTAGAGAGGATACACCTCCATCACAACCCGTGGA
 ACTGTAACGTGATATCTGTGGCTCAGCTGGTGGATAAGAGACATGGCCCCCTCCAACACAGCTTGCTG
 TGCCAGGTGTAACTCCCCCAACCTGAAAGGGAGGTACATCGGAGAGCTGGACCAGAATTACTTTACA
 TGCTATGCTCCCGTAATTGTGGAGCCCCCTGCAGACCTCAATGTCACTGAAGGCATGGCAGCTGAGCTGA
 AATGTCGGGCTTCTACGTCCTGACTCCGTATCTTGGATTACTCCAAATGGAACAGTCATGACCCACGG
 GGCATACAAAGTGCGGATAGCTGTGCTCAGCGACGGTACTTTAAATTTCAAAATGTAACCTGTGCAAGAC
 ACAGGCATGTACACATGTATGGTGAATAATTCTGTTGGCAACACCACTGCTTCTGCCACCTTGAATGTTA
 CTGCGGCAACCACTACTCCGTTCTCGTACTTTTCAACTGTAACAGTAGAGACTATGGAACCTTCTCAGGA
 TGAGGCACGGACCACAGATAACAATGTGGGCCCACTCCAGTGATCGATTGGGAGACCACCAATGTAACC
 ACATCTCTTACGCCACAGAGCACAAGGTCGACAGAAAAACATTACCATCCCAGTAACCTGACATCAACA
 GCGGAATCCCAGGAATTGATGAGGTCATGAAAACAATAAAATCATTATTGGGTGTTTTGTGGCCATCAC
 ACTCATGGCTGCCGTGATGCTGGTCATTTTCTACAAGATGAGGAAACAGCACCATCGGCAAAACCACCAT
 GCTCCAACAAGGACTGTTGAAATCATTAACTGGATGATGAGTCACTGGGGACACGCCCATGGAAGCC
 ACCTGCCCATGCCTGCGATTGAGCATGAGCACCTAAACCACTATAACTCTTACAAATCTCCCTTCAACCA
 CACAACAACAGTAAACACAATAAATTCAATACACAGTTCACTGATGAACCGTTATTGATCCGAATGAAC
 TCTAAAGACAATGTACAAGAGACTCAGATA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence:


>MG209693 representing NM_178725
 Red=Cloning site Green=Tags(s)

MLNKMTLHPQQIMIGPRFNRLFDPLLVLALLQLLVAGLVRAQTCPSCVCSNQFSKVICVRKNLREV
 PDGISTNTRLLNLHENQIQIIVNSFKHLRHLEILQLSRNHRTIEIGAFNGLANLNTLELFDNRLTTIP
 NGAFVYLSKLELWLRNNPIESIPSYAFNRIPSLRRLDLGELKRLSYISEGAFGLSNLRYNLAMCNLR
 EIPNLTPLIKLELDLGNHLAIRPGSFQGLMHLQKLWMIQSQIQVIERNAFDNLQSLVEINLAHNNLT
 LLPHDLFTPLHHLERIHLLHNPWNCNDILWLSWWIRDMAPSNTACCARCNTPPNLKGRYIGELDQNYFT
 CYAPVIVEPPADLVNTEGMAAELKCRASLTSSWITPNGTVMTHGAYKVRIAVLSDGTLNFTNVTVDQ
 TGMVTCMVSNSVGNNTASATLNVTAATTPFSYFSTVTVETMEPSQDEARTDNNVGPPTVIDWETTNVT
 TSLTPQSTRSTEKFTTIPVTDINSIGIPGIDEVMKTKIIIGCFVAITLMAAVMLVIFYKMRKQHRQNH
 APTRTVEIINVDEITGDTMPESHLPMPAIEHEHLNHNYSYKSPFNHTTTVNTINSIHSSVHEPLLRNM
 SKDNVQETQI

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_178725

ORF Size: 1920 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.


RefSeq: [NM_178725.6](#), [NP_848840.3](#)

RefSeq Size: 3475 bp

RefSeq ORF: 1923 bp

Locus ID: 241568
UniProt ID: [Q8C031](#)
Cytogenetics: 2 E1
Gene Summary: May promote neurite outgrowth of developing thalamic neurons.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MG209693