

Product datasheet for **MG209821**

Larp5 (BC078640) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Larp5 (BC078640) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Larp5
Synonyms:	A630096F19; AI256361; D13Wsu64e; Larp5
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>MG209821 representing BC078640
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGC**

ATGACTTCTGATCAGGACGCTAAAGTTGTGGCTGAACCGCAGGCACAGCGAGTCCAGGAGGGCAAGGACA
 GCTCTCATCTGATGAATGGTCCTATATCGCAAACCACTTCTCAAACAAGGTCCCTCCAGCTTTGACTCA
 GGTCCCAACAATAAGGTTTCAGAGCTAAACCCTAATGCAAAGGTATGGGGGACTCATATGTTACATTTG
 GAAGCGAGTAGTGCAGCTGTTGGTGTGAATGCTGCGTGGGAAGAGGCTCCCGGCCATCTACAGACTGTG
 ACCAGCAGGTTCTAGGATTGGATGCCAATGGTGATGGCGACAAAAGTCGAGAGAATGCAGCGCTCCCGGA
 CGCACAGGAGGCAGAGCAGACAGACATGAGCACTCTGGCTCTGGACCATTGAGAAATGAACCTCTGCCT
 GAGAACAACGATACGGGAGGAAATGAGTCTCAACCAGAAAGCCAGGAAGACCTCGGGAAGTACTTAAAA
 AAACATTGGAATTCTGCTTATCTAGGGAGAATCTTGCTAGCGACATGTACCTTATATCCAGATGGACAG
 TGATCAGTATGTGCCAATTACAACGGTAGCTAACCTTGACCATATCAAGAACTCAGTACTGATGTGGAT
 TTGATTGTTGAAGTGCTAAGATCTTTACCTTTAGTCCAAGTGGATGAGAAAGGAGAGAAAGTAAGGCCAA
 ATCAGAACCGGTGCATAGTAATACTGCGTGAGATATCGGAGTCTACCCCGTGGAGAAGTAGAAGCACT
 ATTTAAAGGAGATAATTTACAAAATTTATAAACTGTGAATTTGCCTATAATGATAATTGGTTTATTACC
 TTTGAAACAGAAGCTGATGCACAGCAGGCTTACAAATACCTTCGAGAAGAAGTCAGAACTTTCCAAGGAA
 AGCCAATTAAGGCGCGGATAAAAGCAAAGGCAATAGCGATAAACACATTTTTGCCAAAGAATGGATTTAG
 ACCTCTGGACATGAACCTTTATACCCAGCAGCGCTATGCTACATCGTTTTACTTACCTCCAGTGTACAGT
 CCCCAGCAGCAGTTCCTCTGTACAGCCTGATCACCCCCAGACGTGGTCAACAACACATAGCTATCTTG
 ATCCACCCTTGGTAACTCCTTTTCCAAGTACTGGATTTATAAATGGGTTTACATCTCCGACCTTCAAACC
 TGCAACGCTCCTCTGACATCACTCAGACAGTATCCTCCTAGAAGCAGGAAAAATTCCTTTGGCTATCGG
 AAGAAAAGAGAAGAGAAGTTTACAAGCAGTCAGACGAGTCTCCAACACCACCAAGCCTCCATCACCGA
 GCTTCGAACTGGGACTGTCCAATTCCCCCATTACCTGGAGCTGCAGGCAATTTGAAGACAGAGGACTT
 GTTTGAAAACAGGCTTTCCAGCTTGATCATAGGATCTTCAAAGAAAGAAACCTCAGTACAGATGCAAGC
 ACAAACTGTTCCTGTAGTGGGACCCAGAGAGCCCTCCGTGCCAGCGCCCTGTGCTGTGTCGGCAGCAT
 TTGAACGTTCCCTTCCCGGTCCATTTACCTGAGGATCCCAAGGTGGCAGAAAAGCAGAGGGAAACTCA
 GAGTGTGGACAGACTCCCTTCTACCCCTACTACAACGTCATGTAATCCGTGCAGGTGAATGGAGCAGCA
 ACAGAATTGAGAAAGCCAGTTATGCAGAGATTTGTGAGAGAACGAGTAAAGACCCTTCTCCTCCTCTC
 CATTGCAACCCCCAAAAGAACAAAAGCCAAGCACTGTGCGCTGTGGGAAGGAGGAGAAGCAGCTGTCTGA
 GCCTGTGAGAGACACCGAGAGCCCCCTGCCCTTAAGTCCACGCCTGGAGTCCCCAAGGATCAGAGGAGG
 CAGCCTGGGCGCGAGCCTCTCCTCCAGCTGCCGGAAACGCCTGAGCAAGGAACAGAACACGCCTCCCA
 AGTCTCCTCAG

ACGCGTACGCGGCGGCTCGAG – GFP Tag – GTTTAA

Protein Sequence: >MG209821 representing BC078640
 Red=Cloning site Green=Tags(s)


MTSDQDAKVVAEPQAQRVQEGKDSSHLMNGPISQTTSTQTRSLPALTQVPTTKVSELNPNKAVWGTHMLHL
 EASSAAGVNAWEEAPGHPTDCDQQLGLDANGDGDKSRENAALPDAQEAEQTMSTLALDHSEYEPLP
 ENNDTGGNESQPESQEDPREVLKKTLEFCLSRENLASDMLISQMSDQYVPITTVANLDHIKKLSTDVD
 LIVEVLRSLPLVQVDEKGEKVRPNQNCIVILREISESTPVEEVEALFKGDNLPKFINCEFAYNDFWIT
 FETEADAQQAAYKYLRREEVRTFQGKPIKARIKAKAIAINTFLPKNGFRPLDMNLYTQORYATSFYLPVYS
 PQQQFPLYSLITPQTWSTHSDYLDPLVTPFPSTGFINGFTSPTFKPATSPSLTSLRQYPPRSRKNSFGYR
 KKREEKFTSSQTQSPTPPKPPSPFELGLSNFPPLPGAAGNLKTEDLFENRSSLIGSSKERNLSTDA
 TNTVPVVGPREPSVPAPCAVSAAFERSPPVHLPEDPKVAEKQRETQSVDRPLSTPTTTACKSVQVNGAA
 TELRKPSYAEICQRTSKDPSSSSPLQPPKEQKPSTVACGKEEKQLSEPVERHREPPALKSTPGVPKQRR
 QPGRRASPPAAGKRLSKEQNTPPKSPQ

TRTRPLE – GFP Tag – V

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: BC078640

ORF Size: 1971 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [BC078640](#), [AAH78640](#)

RefSeq Size: 4382 bp


RefSeq ORF: 1973 bp

Locus ID: 217980

Cytogenetics: 13 3.66 cM

Gene Summary: Stimulates mRNA translation.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MG209821