

## Product datasheet for **MG209824**

### **Ddx3y (NM\_012008) Mouse Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	Ddx3y (NM_012008) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Ddx3y
Synonyms:	8030469F12Rik; D1Pas1-rs1; Dby
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide  
Sequence:


>MG209824 representing NM\_012008  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGCATCGCC**

ATGAGTCAAGTGGCAGCGGAAAGTACGGCCGGCCTGGACCAGCAATTTGTTGGTCTTGACCTGAAATCAT  
CAGATAATCAGAATGGAGGAGGAAATACAGAGAGCAAAGGGCGCTATATACCTCCTCACTTAAGGAACAG  
AGAAACATCTAAAGGGTCTGTGATAAGGACAGTTCAGGATGGAGCTGTAGTAAAGATAAAGATGCCTAC  
AGCAGTTTTGGATCTCGTGATCCAGAGGGAAGCCCAATTATTTTCAGTGATCGTGGAAAGTGGATCCAGGG  
GAAGGTTTCGATGATCATGGTTCGAAATGACTATGATGGTATTGGTGGTCTGACCGAACTGGATTTGGCAA  
ATTTGAACGGAGTGGTCATAGTCGTTGGAGTGACAGATCAGATGAAGATGACTGGTCAAAACCACTTCCA  
CCAAGTGAACGCTTAGAGCAAGAACTTTTTCTGGAGGAAACACAGGGATTAACCTTTGAGAAATATGATG  
ATATACCAAGTAGAAGCAACTGAAATAATTGCCTCCTCACATTGAAAATTTTCAGTGATATTGAGATGGG  
AGAAATTATCATGGGGAACATTGAACCTACTCGTTACACTCGTCTACTCCAGTGCAAAAACATGCTATT  
CCAATTATTAAAGAAAAAGGGACTTAATGGCTTGTGCTCAAACAGGGTCTGGAAAAACTGCTGCATTTTC  
TACTACCCATTCTCAGTCAGATTTATACAGATGGTCCAGGAGAGGCTTTGAAGGCTATGAAGGAAAATGG  
AAGATATGGTCGCCGTAACAATACCCAATCTCCTTGGTGTAGCCCCAACAAAGAGAATTGGCTGTGCAG  
ATCTATAGGGAAGCCAGAAAATTTTCATACCGATCTAGAGTCCGTCCTTGTAGTGTATGGTGGTGGTCTG  
ATACTGTTTCAGCAGATTCGGGACTTAGAACGTGGATGCCACTTGTAGTTGCCACACCAGGACGTCTAGT  
GGACATGATGGAAAGAGGAAAGATTGGATTAGATTTCTGCAAACTTAGTATTGGATGAAGCTGATAGG  
ATGCTGGATATGGGATTTGAACCTCAAATACGTCGTATAGTTGAGCAGGACACAATGCCACCAAAGGGGG  
TTCGTCACACTATGATGTTTAGTGTACTTTTTCTAAGGAGATACAGATGCTTGTGCTGACTTTTTTAGA  
TGAATATATCTTTCTGGCTGTAGGCAGAGTAGGATCTACTTCTGAGAACATCACACAGAAAAGTAGTATGG  
GTTGAAGAACTAGACAAACGGTCATTTTTACTTGATCTCTTAAATGCAACAGGGAAGGATTCATTGACTT  
TAGTATTTGTGGAGACTAAAAGGGTGCAGATTCGCTGGAGAACTTCTTATTTCAAGAAAGATATGCGTG  
TACTAGTATTCATGGAGACAGATCACAGAAAGATCGAGAGGAGGCCCTTACCAGTTTCGCTCAGGGCGA  
AAGCCAATTTTAGTGGCTACAGCTGTGGCAGCAAGAGGTCTAGACATTTCAAATGTGAAACATGTTATCA  
ATTTTGATTTGCCAAGCGATATTGAAGAATATGTGCATCGAATTGGCCGCACAGGACGTGTAGGAAACCT  
TGGACTTGCCACCTCTTTTTCAATGAGAGGAATTTGAACATCACAAAAGATTTGTTAGATCTCCTTGT  
GAGGCTAAACAAGAAGTGCCTTCTGGTTGAAAGTATGGCTTATGAACACCACTACAAGGGAAGTAGCC  
GTGGACGTTCTAAAAGCAGATTCAGTGGAGGATTTGGAGCCAGAGACTATCGACAGAGCAGTGGTTCTGC  
AAATGCTGGATTTAATAGTAACCGTGCCAATAGCAGCCGAAGTAGTGGTAGTAGCCACAACAGAGGATTT  
GGTGGAGGTGGCTATGGAGGTTTCTACAATAATGATGGTTATGGAGGAAATTATAACTCCCAGGCAGTTG  
ACTGGTGGGGCAAT

**ACGCGT**ACGCGGGCCGCTCGAG – GFP Tag – GTTTAA



**Plasmid Map:**


**ACCN:** NM\_012008

**ORF Size:** 1974 bp

**OTI Disclaimer:** Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in *E. coli* are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_012008.2</a> , <a href="#">NP_036138.1</a>
<b>RefSeq Size:</b>	4640 bp
<b>RefSeq ORF:</b>	1977 bp
<b>Locus ID:</b>	26900
<b>UniProt ID:</b>	<a href="#">Q62095</a>
<b>Cytogenetics:</b>	Ypter
<b>Gene Summary:</b>	Probable ATP-dependent RNA helicase. May play a role in spermatogenesis.[UniProtKB/Swiss-Prot Function]