

Product datasheet for **MG209845**

Pus7 (BC008544) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Pus7 (BC008544) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Pus7
Synonyms:	C330017115Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>MG209845 representing BC008544
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGAATGACAAGCACGTCCCTGAAACGTGGGTGCCTTGTGTGGAAGATAATGACAGTGTGACCCAC
ACGACGAAACAAAGAAACAAGAAGGTGTCCGAGGGCTGTCTGACCAGCAGTCAGGATGGGGTAGAGAATGA
CGGTCTACACCGCAGTGAGAACGAGCCTGGGCCTCCAGAAGCTGAAAGCACCGTAAAGATGACGAGAAT
TCCAGCGCTCAGGTGCAAGAGGAAGAGGAGGAGGAGGAGGAAGACGGCCTTTCAGAAGCAGGGGAAG
AGGAAGAAGCAGAGAGTTTTGCAGACATGATGAAACATGGGCTCACAGAACTTGATGTGGGCATCTGCAA
GTTCTGTGAGTTCCCATCATGGCTTCTCAGGGATCTTAAAAGAAAGATACTCTGACTTTGTCGTTACAGAA
ATCGGGAAAGATGGACGGATCAGCCATTTGGATGACTTATCTGTTCCGGTGGACGAAGAGGACCCTCCAG
AAGACGCTCTCACTGTCTTGACAGCTGAGGACAGGCAGCAGTTGGAGGAGCTCCAGCTTTTCAAAAATAA
AGAAACCAAGTGTGCCATCGAGGTCAATGAGGACACCAAAAGAGAAGCGAACAGTCATCCATCAAGCCATT
AAATCTCTCTCCAGGACTAGAGACAAAGACAGAGGATAGAGAGGGGCGGAAATACATCGTAGCCTACC
ATGCAGCTGGCAAAAAGGCTTTGGCAAATCCAAGAAAACATTCTTGCCAAAATCCCGGGGAAGTACTG
CCATTTTGTCTCTACAAGGAAAACAAAGACACCATGGATGCCATCAATGTTCTATCCAAGTACCTGCGG
GTAAAGCCAAACATATTCTCCTATATGGGAACCAAGACAAAAGGGCAATCACAGTCCAGGAGATTGCTG
TTCTCAAATAAGTGCACAAAGACTTGCCACCTCAATAAGTGCCTGATGAACCTGAAGCTCGGGAATTT
CAGCTATCAGAAAACCCCTCTGAAGTTAGGAGCACTCCAAGGAAACCACTTCACTGTTGTTCTCAGGAAT
ATAACAGGAAGTATGAGCAGGTGCAGCAAGCCATGCAGTCTCTGAGGGAGACCGGCTTCATTAAGTACT
ATGGCATGCAAAGGTTCCGGACCAAGCTGTTCCAACGTACCAGGTTGGAAGAGCAATTCACAGAATTC
CTGGACAGAAGTCATGGATTTAATACTGAAACCCGTTCTGGAGCTGAGAAGGGTTACTTGGTGAATGC
AGAGAAGAATGGGCCAAGACCAAGACCCAGCCTCTGCTCTCAAAAAGCTGCCCGTCAAAAAGGTGTGTGG
AAGGGCAGCTGCTCCGAGGACTTCCAGATATGGAATGAAGAATATAGTCTCTGCATTTGGCATAATTCC
AAGAAACAATCGTTAATGTACATCCACAGTTATCAGAGCTATGTGTGGAACACCATGGTGAAGCAGAAGA
ATTGAGGAGTACGGGCTGAGGCCTGTGCCCGGGACCTTGTCTCAAAGGAGCCACGCCACCTACATTG
AGGAAGATGACGTTGACAACACTCCATCCACGATGTGGTATGCCTTTGCCTGGGTTTGTATTTATCTA
CCCCAAGCATAAAATCAGTGAGGCCTACAGGGAAATGCTTGTCTGCAGACAACCTTGATATTGACAACATG
AGACACAGATTGAGATTACTCCTTATCGGGGCTACAGAAGGATCATCATCAGGCCTCAGAGCGTCA
GCTGGGAAGTTGTTGCATATGATGATCCTAAAATTCACCTTTTCAACACGGATGTGGACAACCTCGAGGG
GAAACCACCACAGTGTGCTTCTGAAGGCAAAACAGGGCTCTGAAGATGGACTTTTCTCTTCTCCT
TCTACTTACGCCACCATGGCCATTTCGAGAAGTGCTCAAAATGGACACCAGCATCAAGAACCAGACCCAGC
TAAATACGAGCTGGCTGCGA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG209845 representing BC008544
 Red=Cloning site Green=Tags(s)


MEMTSTSLKRGCLVVEDNDSVTPHDETKKQKVSEGCLTSSQDGVENDGLHRSENEPGPPEAESTVKDDEN
 SSAQVQEEEEEEEEEDGLSEAGEEEAESEFADMMKHGL TELDVGICKFVSSHGFSGILKERYSDFFVHE
 IGKDGRISHLDDL SVPVDEEDPPEDAL TVLTAEDRQQLLEELQLFKNKETSVAIEVIEDTKEKRTVIHQAI
 KSLFPGLETKTEDREGRKYIVAYHAAGKKALANPRKHSWPKSRGSYCHFVLYKENKDTMDAINVLSKYLR
 VKPNIFSYMGTDKRAITVQEIAVLKISAQRLAHLNKCLMNLKLGNF SYQKTPLKLGALQGNHFTVVLRN
 ITGTDEQVQAMQSLRETGF INYYGMQRFGT TAVPTYQVGRAILQNSWTEVMDLILKPRSGAEKGYLVKC
 REEWAKTKDPASALKKLPVKRCVEGQLLRGL SRYGMKNI VSAFGIIPRNNRLMYIHSYQSYVWNTMVSRR
 IEEYGLRPVPGDLVLKGATPTYIEEDVDNYSIHDVVMPLPGFDVIYPKHKISEAYREMLAADNLDIDNM
 RHTIRDYSLSGAYRRIIIRPQSVSWEVVAYDDPKIPLFNTDVDNLEGGKPPPVFASEGKYRALKMDFSLPP
 STYATMAIREVLKMDTSIKNQTLNTSWLR

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



Plasmid Map:



ACCN: BC008544

ORF Size: 1982 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [BC008544](#), [AAH08544](#)

RefSeq Size: 3137 bp

RefSeq ORF: 1982 bp

Locus ID: 78697

Cytogenetics: 5 A3

Gene Summary: Pseudouridylate synthase that catalyzes pseudouridylation of RNAs. Acts as a regulator of protein synthesis in embryonic stem cells by mediating pseudouridylation of RNA fragments derived from tRNAs (tRFs): pseudouridylated tRFs inhibit translation by targeting the translation initiation complex. Also catalyzes pseudouridylation of mRNAs: mediates pseudouridylation of mRNAs with the consensus sequence 5'-UGUAG-3'. In addition to mRNAs and tRNAs, binds other types of RNAs, such as snRNAs, Y RNAs and vault RNAs, suggesting that it can catalyze pseudouridylation of many RNA types.[UniProtKB/Swiss-Prot Function]