

Product datasheet for **MG209860**

Enox1 (BC052062) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Enox1 (BC052062) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Enox1
Synonyms:	B230207J08; D230005D02Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>MG209860 representing BC052062
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGTAGATGCGGCTGGATTTGAGAGCGTTGCCAGTGTCCACGGGAGCTTCATCAGATGATGGCTGCAG
 CAGCTGATGGTTTGGGGAGCATAGCCCTAGATACCAACCAGCTCAACATGTCCGTACCGATCCACAGC
 CTGGGCCACAGCTATGAATAACTTGGGCATGGTCCCGTGGGGTTGCCTGGACAGCAGCTCGTGTCCGAC
 TCAATCTGCGTGCCGGGCTTCGATCCAGGCCTCAACATGATGACTGGAATCACCCCATTAACCCGATGA
 TACCCGGTCTGGGGCTGGTACCACCCCAACCGAAGTGGCTGTCGTCAAGGAAATAATCCACTG
 CAAAAGCTGTACTCTTTCCCAAAATCCAAATCTTCCACCTCCTTCCACAAGAGAACGACCTCCAGGG
 TGTAAGACTGTGTTTGTGGAGGATTACCAGAAAATGCAACGGAAGAAATTATTCAAGAAGTCTTTGAGC
 AGTGTGGAGATATCACAGCCATCCGAAAGAGCAAGAAGAAATTTCTGTACATTGCTTTGCCGAGGAATT
 CATGGTTGACAAAGCCATTTACCTCTCTGTTACAGAATGCGATTAGGGTCGAGCACGGATAAAAGGAT
 TCCGGCCGCTGCACGTGGATTTTGCTCAGGCCAGGGATGACTTCTATGAGTGGGAATGCAAGCAAAGGA
 TGCGCGCCCGGGAGGAGCGGCATCGACGCAAACTGGAAGAGGACAGGCTGCGGCCCTCCCTCGCCGCCGGC
 CATCATGCACTACTCGGAGCACGAAGCTGCCCTGCTGGCTGACAACTGAAAGATGACAGCAAGTTCTCT
 GAGGCTATCACAGTTCTGCTGTCTGGATCGAGCGTGGGAGGTGAACCGCGCTCGGCAAAACGATTCT
 ACTCCATGGTTTCACTCGGCCAACAGCCAGCTGCGCAGGCTCATGAATGAAAGGCTACCCACGAGCAGGA
 AATGGAGGAAGCAAAGGAGAAGTTCAAAATGCCTTAACGGGGATCCTCACTCAATTTGAGCAGATTGTG
 GCCGTTTTCAACGCTTCTACCAGACAAAAGCTTGGGACCATTCTCGAAAGCCAGCGCAAGAACATAG
 ACATTTGGCGAAAGCACTCTGAGGAGCTCCGGAACGCTCAAAGCGAGCAGCTCATGGGCATCCGCCGTGA
 GGAGGAGATGGAAATGTCTGATGATGAGAACTGTGACAGCCCGACTAAAAAATGAGAGTCGATGAATCA
 GAGGCAGACATTGAAGACTGTGACGAGACCCGCGGATCGTGGCTTGCAATTCTGCTGCCCTGGCTGCTC
 AGGCCTACGCCCTCAAAGAGGAGAATGACAGTCTCCGGTGGCAACTGGATGCCTACAGGAACGAGGTGGA
 GCTGCTGAAGCAAGAGAAGGAACAGCTTCCGGACAGAAGAAACCTCACCAAGGACCAACAGCTACAG
 TTCCTGCAGCAAAACCATGCAAGGCATGCAGCAGCAATTGCTGGCCATCCAGGAGGAGCTGAACAACAAA
 AGTCGGAATTGGAGCAGGCAAAGGAGGAGCAGTCTCACACGCAAGCGTTGCTGAAAGTCTCCAGGAGCA
 GTTAAAGGTACCAAGGACTTGGTGGAGACGAATGGCCACAGCCACGAGGATGCAATGAAATCAATGTG
 TTGACGGTCGCGTTGGTTAACCAGGACCGAGAGAAACAACACTGAAAAGAGAAGCCAAGGTTTGAATCTG
 AGAAGGAGGCCCTGCTGATAGGCATCATATCAACGTTTCTTACGTTTACCCCTTTGGAGCTAATATAGA
 GTATCTGTGGTCATACATGCAGCAGCTGGACTCTAAGATATCCGCAATGAAATTGAAATGCTTCTGATG
 AGGCTGCCACGGATGTTCAAACAGGAATCACTGGTGTGGGAGCGACGCTAGAGAAGAGATGGAAGTTGT
 GTGCCTTTGAAGGAATTAACGACC

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence: >MG209860 representing BC052062
 Red=Cloning site Green=Tags(s)


```

MVDAAGFESVAQCPRELHQMAAADGLGSIALDTTQLNMSVTDPTAWATAMNNLGMVPVGLPGQQLVSD
SICVPGFDPGLNMMTGITPINMIPGLGLVPPPPTEVAVVKEIIHCKSCTLPQNPPLPPSTRERPPG
CKTVFVGGLPENATEEIIQEVFEQCGDITAIRSKKNFCHIRFAEEFMDKAIYLSGYRMRLGSSTDKKD
SGRLHVDFAQARDDFYWECKQRMREERHRRKLEEDRLRPPSPPAIMHYSEHEAALLADKLKDDSKFS
EAITVLLSWIERGEVNRSSANQFYSMVQSANSHVRRMLMNEKATHEQEMEEAKENFKNALTGILTQFEQIV
AVFNASTRQKAWDHFSAQRKNIDIWRKHSEELRNAQSEQLMGIREEEMEMSDDENCDSPTKMRVDES
EADIEDCAGDPADRGHSAALAAQAYALKEENDSLRWQLDAYRNEVELLKQEKEQLFRTEENLTKDQQLQ
FLQQTMQGMQQLLAIQEELNNKKSELEQAKKEQSHTQALLKVLQEQLKGTDLVETNGHSHEDANEINV
LTVALVNQDRENNTEKRSQGLKSEKEALLIGIISTFLHVHPFGANIEYLWSYMQQLDSKISANEIEMLLM
RLPRMFKEFTGVGATLEKRWKLCAFEGIKT
  
```

TRTRPLE – GFP Tag – V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: BC052062

ORF Size: 1986 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [BC052062.1](#)

RefSeq Size: 3117 bp


RefSeq ORF: 1988 bp

Locus ID: 239188

Cytogenetics: 14 D3

Gene Summary: Probably acts as a terminal oxidase of plasma electron transport from cytosolic NAD(P)H via hydroquinones to acceptors at the cell surface. Hydroquinone oxidase activity alternates with a protein disulfide-thiol interchange/oxidoreductase activity which may control physical membrane displacements associated with vesicle budding or cell enlargement. The activities oscillate with a period length of 24 minutes and play a role in control of the ultradian cellular biological clock (By similarity).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MG209860