

## Product datasheet for **MG209879**

### **Cd6 (NM\_009852) Mouse Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	Cd6 (NM_009852) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Cd6
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide  
Sequence:

>MG209879 representing NM\_009852  
Red=Cloning site Blue=ORF Green=Tags(s)


TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCCGCGATCGCC

ATGTGGCTCTTCCTTGGGATCGCTGGATTGCTGACAGCAGTTCTCTCAGGTCTCCCATCTCCAGCCCCCT  
CTGGCCAGCACAAAAATGGCACCATCCCGAACATGACCTTGGACCTAGAGGAGCGGCTGGGATTGCGCT  
GGTGAACGGGAGCAGTCGCTGCAGCGGATCTGTGAAGTTCTACTGGAGTCTGGGAGCCTGTGTGTGCG  
GCACACTGGAACAGGGCTGCTACGGAGGCTGTGTGCAAAGCGCTGAACTGTGGTATTCGGGAAAGGTTA  
CCTACCTGATGCCACCCACCTCTGAGCTTCCACAGGAGCCACTTCTGGGAACACAAGCAGCGCCGGGAA  
TACGACGTGGGCGCGGCCACCACTGAGAGATGCAGGGGAGCCAATTGGCAGTTCTGCAAGGTGCAGGAC  
CAGGAGTGTAGCAGCGACAGGCGGCTGGTGTGGGTACCTGTGCAGAGAACCAAGCGGTGCGGTTAGTGG  
ATGGTAGCAGCCGCTGTGCTGGCCGAGTGGAGATGCTGGAGCATGGTGAATGGGGTACCGTGTGCGACGA  
CACGTGGGACCTGCAGGATGCACACGTAGTGTGAAGCAACTGAAGTGGCGCTGGGCAGTGAAGGCCTG  
GCCGGCTTGCACCTCACCCAGGCCAAGGACCCATCCACCGCATCAGGTGAACTGCTCTGGGACTGAGG  
CTTACCTGTGGGACTGCCCGGGACGGCCAGGAGACCAGTACTGTGGTCAACAAGGAGGACGACAGGAGTGGT  
GTGCTCAGAGCACAATCCTGGCGCCTGACTGGAGGCATTGACTCCTGCGAAGGGCAAGTGGAAAGTGTAC  
TTCCGGGTGTCTGGAGCACGGTGTGTGACAGCGAGTGGTACCCATCTGAGGCCAAAGTGTCTGTGCGT  
CCTTGGGTTGCGGATCAGCAGTGGCTCGGCCAGGGGGCTGCCACACTCCCTGGACGGTAGAATGTACTA  
CTCCTGCAAAGGACAGGAGCCCGCCCTCAGCACCTGTTCTGGCGGTTCAACAACCTCAACCTGTGCAGC  
CAGTCGCGGGCAGCCAGGGTTGTCTGCTCAGGCTCCCAGAGACATCTCAATCTGTCTACTTCTGAAGTGC  
CTTCCAGAGTTCAGTCACGATAGAATCTTCTGTGCCGGTGAAGGATAAGGACTCACAAGGGCT  
GACACTTCTCATTCTGTATTGTCCTGGGAATTCTCCTCCTCGTCTCCACCATCTCATAGTTATTCTT  
CTCTTGAGAGCCAAAGGACAATACGCCCTCCCTGCTTCAGTGAACCACCAGCAACTGTCTACTGCCAACC  
AAGCAGGGATCAATAACTACCATCCTGTCCCTATCACCATTGCCAAAGAAGCTCCTATGCTGTTTATCCA  
GCCCGGGTCCCTGCAGACTCGGACTCCAGTTCTGACTCCGACTATGAGCACTATGACTTCAGCTCCAG  
CCGCTGTGGCCCTGACCACCTTCTACAATTCCCAGCGACACCGGGTACAGAGGAGGAAGCTCAGCAGA  
ACAGGTTTTCAGATGCCACCCTTGGAGGAAGGACTTGAAGAGCTGCATGTTTCCACATCCCAGCGGCTGA  
CCCCAGACCCTGTGTGGCAGATGTCCCCTCTCGGGGCTCTCAGTACCACGTGAGGAACAACAGTATTCC  
AGCACCTCCTCTGAGGAGGGCTACTGCAATGATCCCAGCAGCAAGCCACCTCCATGGAATCCCAGGCGT  
TTTATTCTGAGAAGAGTCCCCTCACAGAGCAGCCCCCAACTTGGAACTGGCTGGCTCTCCGGCAGTGT  
TTCAGGGCCCTCAGCCGATGACAGCTCCAGCACCTCATCCGGAGAGTGGTACCAGAACTCCAGCCACCA  
CCCCAGCATCCTCCAGCAGAACAGTTTGAATGTCCAGGACCCCTGGCCCCAAACAGACTCCATTGATG  
AGGAGGAGGAGGATTACGATGACATCGGAGCAGCC

ACGCGTACGCGGGCCGCTCGAG – GFP Tag – GTTTAA



## Plasmid Map:



**ACCN:** NM\_009852

**ORF Size:** 1995 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_009852.2](#)

RefSeq Size: 3034 bp

RefSeq ORF: 1998 bp

Locus ID: 12511

UniProt ID: [Q61003](#)

Cytogenetics: 19 7.16 cM

**Gene Summary:** Cell adhesion molecule that mediates cell-cell contacts and regulates T-cell responses via its interaction with ALCAM/CD166. Contributes to signaling cascades triggered by activation of the TCR/CD3 complex (PubMed:24584089). Functions as costimulatory molecule; promotes T-cell activation and proliferation. Contributes to the formation and maturation of the immunological synapse. Functions as calcium-dependent pattern receptor that binds and aggregates both Gram-positive and Gram-negative bacteria. Binds both lipopolysaccharide (LPS) from Gram-negative bacteria and lipoteichoic acid from Gram-positive bacteria. LPS binding leads to the activation of signaling cascades and down-stream MAP kinases. Mediates activation of the inflammatory response and the secretion of pro-inflammatory cytokines in response to LPS.[UniProtKB/Swiss-Prot Function]