

Product datasheet for **MG209963**

Uroc1 (NM_144940) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Uroc1 (NM_144940) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Uroc1
Synonyms:	A1265605; BC022133
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MG209963 representing NM_144940
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGTCCAGCCTCCAGGAAGTGTGCTCTGGTCTGCCCTACGACCCCTTCTGAGAACCAGGGTTCGCTGGG
 CTGGGGTCCACATGCTCCGGTTAGAAGTCCGAACCTCAGCCCTGAAGAAGAGCAGCTGGCCCTGAGGAA
 TGCTTTCGCTATTTCCACCTGATGTCCAGAACTGTTGGCCCTGGAGTTTGCCAGGAGCTAAGACAG
 TTTGGACACATCTACATGTACCGTTCTGCCCCAGCATTGAGATGAGAGCCTACCCAATTGATCAGTACC
 CATGCCGGACACGGGAGCTGCAGCCATCATGCACATGATCATGAACAATTGGACCCTGCTGTGGCCCA
 GTTTCCACAGGAGCTGGTACCTACGGAGGCAACGGGCAAGTGTTCAGCAACTGGGCCAGTTCCGGCTG
 ACCATGTCCTACCTGTGAAGATGACTGAGGAACAGACTCTGGTCATGTACAGTGGACACCCCTCTGGGCC
 TCTTCCCAGCAGCCCCGTGCCACGGCTGGTCACTAACGGGATGGTCATTCCCAACTACTCCTC
 CAGGACAGAGTATGAGAAGCTGTTTGCCTTGGGGGTGACCATGTATGGCCAGATGACAGCGGGCAGCTAC
 TGCTACATTGGTCCCCAGGGAATCGTCCATGGCACAGTGTCTACTGTGCTGAATGCTGGTCTGGTACC
 TGGGCATCGAGAAGTGGCAGGGAAGGCTTTGTACCTCTGGGCTTGGCGGGATGAGTGGGGCTCAAGC
 CAAGGCAGCAGCCATCGTTGGATGCATTGGGGTATTGCGGAGGTGGATAAAGCTGCACTTGTAAAACGC
 CACCGGAAGGCTGGCTGATGGAGGTGACAGACAGCTTGGACCCTGCATTGCAAGACTCAGGGAAGCAA
 GGAAAAAGAAGGAGGTTCTCAGCCTCGGCTACCATGGCAACGTGGTGGATCTTTGGGAGCGCCTTGCCA
 TGAAGTGGACACAACGGGGAAGTCTTGGTGGACCTCGGGTCAGACCAGACATCCTGTACAACCCGTTT
 AATGGCGGCTACTACCTGTGCAGCTCAGCTTCTCAGAGGCCAGAGCCTCATGCTCCTCAACCCGCTG
 CCTTCAAGCACCTGGTGCAGAAAGCCTGAGGAGGCATGTCGACCCATCAACAGTTGGCCCAAGAGAA
 GTTCTTCTTGGGACTATGGTAATGCCTTCTTGGAGGCTCAGAGAGCAGGAGCAGATGTAGAGAAG
 AAAGGAGCCAACAAGATGGAGTTCGCTACCCCTCCTATGTCCAGCACATTATGGGGATATATTCTCCC
 AGGGCTTTGGGCCCTTCCGCTGGGTATGCACATCAGGGGACCCCAAGGACCTGGCTGCTACTGATCATCT
 GGCCACATCTGTACTGGAGAAAGCCATTGCTGATGGAGTGAAGGCATCTGTGAAGCTACAGTATATGGAT
 AATATACGCTGGATCCGAGAGGCTGCCAAGCATCAGCTGGTGGTGGGCTCCAGGCAAGGATCCTGTACT
 CAGACCAGAAAGCCGTGTGGCCATTGCTGTGGCCATTAACCAGGCCATCGCCAGTGGGAAGATCAAGGC
 ACCAGTGGTCTGAGCCGTGACCACCATGATGTGAGTGGCACTGACAGCCCTTCAGGGAGACTTCGAAT
 ATTTATGACGGTTCTGCCTTCTGCGCAGACATGGCCGTGCAGAACTTCGTGGGAGATGCCGTGCTGGTG
 CCACCTGGGTGCACTTACAATGGAGGGGGTGTAGGCTGGGGTGAAGTCATCAATGGGGGATTTGGCCT
 TGTCTAGATGGAAGTGCAGAGGCTGAGCAGAAAGCCAGGATGATGCTCAGCTGGGATGTCTCCAATGGT
 GTGGCACGACGCTGCTGGTCTGGAACCCCAAGGCCTATGAGATCATCTGCCAGACAATGCAGGAGAACA
 GTGGCCTGGTGGTACGCTTCCCATGAGGTGGCAGATGAACAAGTGTACAGCAGGCCCTCGCGCCC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG209963 representing NM_144940
Red=Cloning site Green=Tags(s)

```
MSSLQELCSGLPLRPLPENRGRWAGVPHAPVRTPNLSPEEEQLALRNALRYFPPDVQKLLALEFAQELRQ
FGHIYMYRFCPSIEMRAYPIDQYPCRTRAAAAIMHMIMNNDPAVAQFPQELVTYGGNGQVFSNWAQFRL
TMSYLSKMTTEEQTLVMYSGHPLGLFPSSPRAPRLVITNGMVIIPNYSSRTEYEKLFALGVTMYGQMTAGSY
CYIGPQGI VHGT VLVNAGRRYLGIENLAGKVFVTSGLGGMSGQAQAKAAAI VGCIGVIAEVDKAAALVKR
HRQGWLMEVTDSDLDRCIARLREARKKKEVLSLGYHGNVVDLWERLVHELDTTGELLVDLGDSDQTSCHNPF
NGGYYPVQLSFSEAQSLMSSNPAAFKHLVQESLRRHVAAINRLAQEKFFWDYGN AFLLEAQRAGADVEK
KGANKMEFRYPSYVQHIMGDIFSQGF GPF RWVCTSGDPQDLAVTDHLATSVLEKAIADGVKASVKLQYMD
NIRWIREAAKHQLVVG SQARILYSDQKGRVAIAVA INQAIASGKIKAPVVL SRDHHDVSGTDSPFRETSN
IYDGS AFCADMAVQNFVGDACRGATWVALHNGGGVGVGEVINGG FGLVLDGTAEAEQKARMMMLSWDVNSG
VARRCWSGNPKAYEII CQTMQENSGLVVTL PHEVADEQVLQQALRP
```

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_144940

ORF Size: 2028 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_144940.2](#), [NP_659189.2](#)

RefSeq Size: 3365 bp


RefSeq ORF: 2031 bp

Locus ID: 243537

UniProt ID: [Q8VC12](#)

Cytogenetics: 6 D1

Product images:



Circular map for MG209963