

Product datasheet for **MG209982**

Tbc1d14 (BC042515) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Tbc1d14 (BC042515) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Tbc1d14
Synonyms:	2810413P16Rik; AU043625; C86258; D5Ertd110e; mKIAA1322
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>MG209982 representing BC042515
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCCGCATCGCC

ATGGGCATCCTAGATGGTCGGCAAGGAGACTCCCTTCAGGACCTACAACACCTTAGTATCAAGGCGGCTC
CCAGATCCCTTTTCAAGTGCCTGAATACGGACCTTCACTAAAACCTGGTGCTTTGGAAGATCGACACAGCCT
TCAATCAGTGGACTCGGGCATTCTACCTGGAGATTGGCAACCCAGAGCCTGTTCCCTGCAGTGTGGTC
CATGTGAAGAGAAAGCAGTCTGAGTCAGAGATCGTCCCGAGCGGGCCTTCCAGAGTGCATGCCCTCGC
CGTCTGCACACCCTCAGCTCCACCTGCAGCGAGCGGGAGCAGGTTGTGCGGAAGTCTTCCACGTTTCC
TAGGACAGGCTATGACTCAGTCAAACCTACAGCCCCACCTCAAAGCCTTGAGCCGAAGTGACAATGTC
TCTGTCTGCAGTGTGTCTAGTCTTGGCACAGAGCTGTCAACTACGTTGTGCGGTGAGCAATGAGGACATCT
TGGACCTCATGGTCACGAGCAATTCAGCGCCATTGTGACCCTGGAGAATGATGATGACCCACAGTTTAC
TGACGTACCCCTGAGCTCCATCAATGAAACCAGCGACTTACACCAGCAGGACTGTGTTGCTGAGACTGAG
GAGGGGAGGAAATTTGAAACTATTGCAACCATTCAGTCACTTCTTTACAAGGAATTTGCTTGTAGAAAAC
AAAATGCAAGGCTTGATAGACAGAGAGACTTGGGATGGAAGCTGTTTGGAAAAGTCCCACTCCGAGAGAC
TGCTCAGAAAAGATTCAAAGAAGACCCAGAAGGAATATGAAGACAAGGCTGGAAGACCGAGCAGGCCACCC
TCTCCCAAGCAGAATGTGAGGAAGAACCCTGGATTTTGGCCACTGTCCACCACGGCGCTCATCTCGAGG
ACAGACCTGCAAACTCTCCAGCCAAGCCAGCTGAAGAAGCTCAAAGCACAGACAGCAGTATGAAGAGAT
GGTGTGTCAGGCCAAGAAACGAGAGCTTAAAGAAGCCAGCGTCGGAGGAAGCAGCTGGAGGAAAGTTG
AAGTTGAAGAGAGCATTGGGAACGCTGTCTCACCTGGAACAATGAGATCTTGCCTAACTGGGAAACCA
TGTGGTGTCTAAGAAGTTTCGAGACTTGTGGTGGCAGGGAATCCCTCCGAGTGAAGAGGCAAAGTTTG
GAGCTTGGCCATTGGCAACGAGTTGAACATCACTCACGAGCTCTTTGACATCTGTCTTGTGCTCGGGCCAAG
GAGCGGTGGCGGTCCCTCAGCACTGGAGGCTCAGAAGTGGAGAATGAAGATGCTGGGTTTTTCAGCAGCAG
ACAGAGAGGCCAGTCTGGAGCTTATTAACCTGGACATTTCTAGAACATTTCTAATCTCTGCATTTTCCA
GCAAGGTGGTCCGTATCATGACATGTTGCACAGTATTTTGGGCGCTTATACTTGTACCGGCCGGATGTG
GGTTATGTTTCAAGGCATGTCTTTATAGCTGCAGTGTGATCTTGAACCTCGATACTGCAGATGCCTTTA
TTGCTTTTTCTAATCTTTTGAATAAACCTGTCAAATGGCGTTTTTCCGAGTGGACCATGGCCTTATGTT
GACTTACTTTGCTGCATTTGAGGTATTCTTTGAAGAAAACCTGCCGAAATTTGTTTGGCATTTCAGAAA
AACAACTTAACTGCAGATATCTACCTAATCGATTGGATTTTACCTTGTACAGTAAGTCCCTGCCCTTGG
ACTTGGCCTGCCGCATCTGGGACGTATTCTGTGAGATGGAGAGGAATTTCTCTTCCGCACAGCTCTGGG
CATCTGAAGCTTTTTGAAGATATCTGACCAGGATGGACTTCATCCACAGTGCCAGTTTCTGACCAGG
CTGCTGAGGACCTGCCCGCAGATGAAGTGTGTCAGCTATCTCCACAGTCCAGATGCAAAGCCGGAACA
AGAAATGGGCACAGGTGCTGTCTGCACTGCAGAAAACAGCCGCGAGATGGAGAAAGGAAGCCATCCCT
CAGACAC

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence: >MG209982 representing BC042515
 Red=Cloning site Green=Tags(s)


MGILDGRQGDSLQDLQHLSIKAAPRSLVPEYGPSLKLGALEDRHSLQSVDSGIPTLEIGNPEPVPCSVV
 HVKRKQSESEIVPERAFQSACPLPSCPTSPAPTCSEREQVVRKSSTFPRTGYDSVKLYSPTSKALSRSDNV
 SVCSVSSLGTELSTLVSINEDILDLMVTSNSSAIVTLENDDDPQFTDVTLSSINETSDDLHQDCVAETE
 EGRKLLKLLQPFSHFFTRNLLARKQNRALDRQDLGWKLFQKVPPLRETAQKDSKKTQKEYEDKAGRPSRPP
 SPKQNVNRKNLDFEPLSTTALILEDRPANLPAKPAEEAQKHRQQYEMVLQAKKRELKEAQRRRKQLEERC
 KVEESIGNAVLTWNNEILPNWETMWCSSKVRDLWWQGIPPSVRGVWVSLAIGNELNITHELFDICLAKAK
 ERWRSLSTGGSEVENEDAGFSAADREASLELIKLDISRTFPNLCIFQQGGPYHMLHSILGAYTCYRPDV
 GYYVQGMFIAAVLILNLDTADAFIAFSNLLNKPCQMAFFRVDHGLMLTYFAAFEVFFEENLPKLFQAFKK
 NNLTADIYLDWIFTLYSKSLPLDLACRIWVDFCRDGEEFLFRTALGILKLFEDILTRMDFIHSAQFLTR
 LPEDLPADEVFAAISTVQMQSRNKKWAQVLSALQKDSREMEKGSPLRH

TRTRPLE – GFP Tag – V

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:


ACCN: BC042515

ORF Size: 2039 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [BC042515](#), [AAH42515](#)

RefSeq Size: 4580 bp

RefSeq ORF: 2039 bp
Locus ID: 100855
Cytogenetics: 5 19.15 cM
Gene Summary: Negative regulator of starvation-induced autophagosome formation.[UniProtKB/Swiss-Prot Function]