

Product datasheet for **MG210041**

Ddr1 (BC065998) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ddr1 (BC065998) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Ddr1
Synonyms:	6030432F18, CD167a, PTK3A
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>MG210041 representing BC065998
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGGGACAGGGACCCTCTCATCTCTACTGCTGCTACTCTTGGTGACAATTGGAGATGCTGACATGA
AGGGACATTTTGACCCTGCCAAGTGCCGCTATGCCCTGGGCATGCAGGACCGCACCATTCTGACAGCGA
TATCTCTGTGTCCAGCTCCTGGTCGACTCTACCCTGCCCGCCACAGCAGGCTGGAAAGCAGTGATGGA
GATGGGGCTTGGTGCCCTGCAGGGCCTGTGTTCCCAAAGAAGAGGAGTACTTGCAGGTGGACCTTCGTA
GGCTACACCTGGTGGCTCTGGTGGGCACCCAGGGCCGATGCTGGGGTCTGGGCAAAGAGTTCTCCCG
AAGCTATCGTTGCGTTACTCCGAGATGGCCGCGCTGGATGGACTGGAAGGACCGCTGGGACAGGAG
GTGATTCGGGTAACGAGGATCCCGGGGAGTAGTGTGAAGGACCTGGGCCCCCATGGTGGCCCGGC
TGGTCCGCTTCTACCCAGGGCTGACCGGGTATGAGTGTCTGTCTTCGGGTGGAGCTCTATGGTCCT
CTGGCGGATGGACTCCTGTCATATACAGCCCCGTGGGGCAGACCATGCAGTTATCTGAGGTGATGGTA
CATCTCAATGATTCCACTTACGATGGATATACTGCTGGAGGGCTGCAATATGGCGGTCTGGGCCAGCTGG
CAGATGGCGTGGTGGCCTGGATGATTTCAAGCAGAGCCAGGAGCTGCGGGCTGGCCAGGCTATGACTA
TGTGGGATGGAGCAATCAGAGCTTCCCCACGGCTACGTGGAGATGGAGTTTGTGTTGATCGGTTGAGG
ACCTTCCAGACCATGCAGGTCCACTGTAACAACATGCACACTCTGGGAGCCCCCTACCAGGCGGGGTGG
AATGCCGTTTAAAAGGGTCCCGCCATGGCCTGGGAAGGAGAGCCTGTCCGCCATGCTCTGGGAGGCAG
CCTTGGAGACCCAGAGCCCGGCCATCTCAGTCCCCGGTGGCCACGTGGGCCGCTTCTGCAGTGC
AGATTCCTCTTTGAGGTCTTGGTACTCTTCAGTGAGATCTTTTCATCTCAGATGTGGTGAACGACT
CCTCTGACACCTTCCCACAGCCCCCTGGTGGCACCTGGCCCGCTCCCACCACTTCCAGCAGCTTGGTA
GCTGGAGCCCCGGGTCAACAGCCAGTGGCCAAGGCGGAGGGGAGCCCAACTGCCATCCTCATTGGCTGC
CTGGTGGCCATCATCCTGCTGCTGCTTCTCATCATCGCGCTGATGCTCTGGAGGCTGCACTGGCGCGGC
TGCTCAGCAAGGCTGAGCGCCGCTGTTGGAGGAGGAGCTGACGGTTCACCTTTCTGTCCCTGGGACAC
CATCCTCATCAACAACCGCCAGGACCCCGAGAGCCACCCCTTACCAGGAGCCCCGGCCTCGGGGACT
CCACCCATTCTGCACCTGCGTCCCCAACGGCTCTGCCTGCAGTGGGACTATATGGAGCCCGAGAAGC
CGGGTCCCCGCTTCTACCCACCTCCCAGAACAGCGTCCCCATTATGCCGAAGCTGACATTGTCAC
CCTGCAGGGCTCACTGGGGCAACACCTACGCTGTGCCTGCACTGCCCCAGGGCGGTGGGGATGGG
CCCCCAGAGTGGATTTCCCTCGGTACGGCTCCGCTTCAAGGAGAAGCTTGGCAGGGCCAATTTGGGG
AGGTACACCTGTGTGAAGTAGAGGACCCGCGAGCCCTTCCCAATGGGAGCCATCCAGGGGAAGCTGGAC
TCGAAAACCAGAGAAGCCTATGGCACCGCACCCACCGTCTCCATACTTGCCATTCCCAAGAGGCAAGT
TCTCTGCGGTGGGCTGGGCTGCCGAGGAGCAGATACACACTCCCTTCCGTTTCTTCTTCCATCT
GCCAAGTGCCACCCTCCAGCTGGCCCTGTGGTTGGGATCCTGCCTGACTTCTCCAGCTTCTTCTTGT
CACACTGGACGCCCGCAGCTGAGAAT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG210041 representing BC065998
 Red=Cloning site Green=Tags(s)


MGTGTLSSLLLLLLLLVTIGDADMKGHFDPKCRYALGMQDRTIPDSDISVSSSWSDSTAARHSRLESSDG
 DGAWCPAGPVFPKEEYLVQDLRRLHLVALVGTQGRHAGGLGKEFSRSYRLRYSRDGRRWMDWDRWGQE
 VISGNEDPGGVVLKDLGPPMVARLVRFYPRADRVMSVCLRVELYGCLWRDGLLSYAPVQGTMLSEVMV
 HLNSTYDGYTAGGLQYGGGLQGLADGVVGLDDFRQSQELRGWPGYDYVGWSNQSFPTGYVEMEFDRLR
 TFQTMQVHCNNMHTLGARLPGGVECRFKRGPAMAWEGEPVRHALGGSLGDPRARAI SVPLGGHVGRFLQC
 RFLFAGPWLLFSEISFISDVVNDSSDTFPPAPWWPPGPPPTNFSLELEPRGQQPVAKAEGSPTAILIGC
 LVAIIIIIIIIALMLWRLHWRLLSKAERRVLEEELTVHLSVPGDTILINNRPGPREPPPYQEPRPRGT
 PPHSAPCVPNGSACSGDYMEPEKPGAPLLPPPQNSVPHYAEADIVTLQGVTTGGNTYAVPALPPGAVGDG
 PPRVDFPRSRLRFKEKLGEGQFGEVHLCEVEDPAALSQWEPSSRGSWTRKPEKPMAPHPTVSI LAHSQEAV
 SLRVGWACPRSRYTHSLPFLPSICQVPPSQLALWLGSLTSSSFLLVTHWTPPAEN


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



Plasmid Map:


ACCN: BC065998

ORF Size: 2060 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [BC065998](#), [AAH65998](#)

RefSeq Size: 2582 bp

RefSeq ORF: 2060 bp

Locus ID: 12305

Cytogenetics: 17 18.7 cM

Gene Summary: Tyrosine kinase that functions as cell surface receptor for fibrillar collagen and regulates cell attachment to the extracellular matrix, remodeling of the extracellular matrix, cell migration, differentiation, survival and cell proliferation. Collagen binding triggers a signaling pathway that involves SRC and leads to the activation of MAP kinases. Regulates remodeling of the extracellular matrix by up-regulation of the matrix metalloproteinases MMP2, MMP7 and MMP9, and thereby facilitates cell migration and wound healing, but also tumor cell invasion. Promotes smooth muscle cell migration, and thereby contributes to arterial wound healing. Phosphorylates PTPN11 (By similarity). Required for normal blastocyst implantation during pregnancy, for normal mammary gland differentiation and normal lactation. Required for normal ear morphology and normal hearing.[UniProtKB/Swiss-Prot Function]