

Product datasheet for **MG210185**

Abcf3 (NM_013852) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Abcf3 (NM_013852) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Abcf3
Synonyms:	AI326318; AU016058; BB119416
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>MG210185 representing NM_013852
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGCGACTTGGCTGATATCCTGCGAAGCGAGTTCCTGAAATTGACGGGCAGGCTTCGACTACGTGA
CCGGTGTCTTGACAGTGGCAGCGCGGATTTTCAGTCTGTGGATGACCTGGTGGAAAGCTGTCGGGAACT
ATTGCAAGAAGTGTCCGGGGACAGCAAGGATGACGCTGGCATTAGGGCTGTCTGTGACGCGCATGTACAAC
ACTCTACGCTGGCTGAACCGCAGAACCAGGAAACAGCCAGGTGCTACTGGACGCCCCATCCAGTTGT
CGAAGATAATGGAAAATGACTGTGACACCAAATTCCAGGACTGCTAAAGAGGGAACAGTCTCTCGAC
AGTAAATGCGAAGAACTAGAGAAGGCTGAAGCTCGACTGAAGGCAAAGCAGGAGAAGCGTTTCAGAGAAG
GAAACTGAAGACCAGCAACCCTCTAGTCTGGAGGAGCGTCAGCTAGCCAGGCAGGCAGCAGAAAAG
AGAGTCGATTGGAATCATCTGGCAAGAACAAGTCTATGATGTGCGAATTGAGAATTTGATGTCTCTTT
TGGTGACAGGGTACTGCTGGCTGGAGCAGATGAACTTGGCATGGGGCCGACGCTACGGCTGGTGGGA
CGCAATGGTCTGGGGAAGACGACATTGCTAAAGATGCTGGCTACCCGGAGCTTGGTGTCCAGCCATA
TTCCCTACTGCATGTAGAGCAGGAAGTTGCTGGAGATGACACTCCTGCCTTGACAGAGTGTGCTGAAAAG
TGACACCGTACGAGAGGACCTGCTACGGCAGGAACGGGAGCTCAGCCTCAGGATTGCCGCTGGCAGGGCA
GAAGGTTGAGAAGCTGCACAGCTGGCAGAAATCTATGGGAAGCTGGAGGAGATTGAGGCAGATAAGGCAC
CTGCCAGAGCATCAGTATTCTTGCGGGCTTGGTTTTACTCCTAAAATGCAGCAGCAGCCTACCAGGGA
GTTCTCAGGTGGTGGAGGATGAGATTGGCCTTAGCCGAGCTCTTTGCCAGGCCAGACCTTCTGCTG
TTAGACGAACCCACAAACATGCTGGATGTAAGGGCCATCTGTGGTTGGAGAATTACCTGCAGACGTGGC
CCTCCACTATCCTAGTTGTCTCTCATGACCGAACTTCTGAATGCCATTGCCACAGACATCATCCACTT
GCATAGCCAACGGCTAGACGGTTACCGGGGGACTTTGAAACCTTCATCAAGAGCAAGCAGGAGCGCCTT
CTCAACCAGCAGCGAGTATGAGGCCAGCAGAGTATCGCAACACATCCAGTTTTTCATCGACCGTT
TTCGATACAATGCAAACAGACCTCTCAAGTGCAGAGTAACTGAAGATGTTGGAGAAGCTACCTGAACT
AAAACCGGTGGACAAGGAATCCGAGGTGGTGTGAAGTCCCCGATGGATTTGAGAAGTTCTCACCGCCC
ATTCTGCAGTAGATGAGGTGGATTTCTACTATGATCCTAAGCACAGCATCTTCAGTCGCTGTCTGTCT
CTGCTGATCTTGAGTCCCGCATCTGTGTGGTGGGGAGAATGGGGCTGGGAAATCTACTATGCTGAAGCT
GCTCATGGGGACCTGTCTCTGTTCGGGTATCAGGCATGCCACAGGAATCTGAAGATAGGCTATTTT
AGCCAGCACACGTGGAACAGCTAGACTTGAATGTCAGTGTGGAACCTCGGCTCGAAAATTTCTG
GGTACCTGAAGAGGAATATCGTCACCAATTGGGCCGGTATGGCATCTCTGGAGAAGTGGCCATCGGTCC
TGTTGCGAGCTTATCTGGGGCCAGAAAAGCCGGTAGCCTTTGCTCAGATGACCATGCCTTGCCCAAC
TTCTATATTCTGGATGAACCCACAAACCCTGGACATGGAGACCATTGAAGCTTTGGGCCAGGCTCTCA
ACAACCTCAGGGGTGGTGAATTTCTGGTATCTCATGATGAGCGCTTATCCGACTGGTGTCAAGGAGCT
GTGGGTGTGTGAGAACGGCAGTGTCACTCGGTTGAGGGGGCTTTGACCAGTACCGGGCCCTCTCCAA
GAACAGTCCCGCGGAGGGCTTCTCT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG210185 representing NM_013852
 Red=Cloning site Green=Tags(s)


MATCADILRSEFPEIDGQVFDYVTGVLHSGSADFESVDDLVEAVGELLQEVSGDSKDDAGIRAVCQRMYN
 TLRLAEPQNGNSQVLLDAPIQLSKIMENYDCDTKLPGLLKREQSSTVNAKKLEKAEARLKAKQEKRSK
 ETLKTSNPLVLEEASASQAGSRKESRLESSGKNKSYDVRIENFDVSFGDRVLLAGADVNLAWGRRYGLVG
 RNGLGKTTLLKMLATRSRVRPAHISLLHVEQEVAGDDTPALQSVLES DTVREDLLRQERELSLRIAAGRA
 EGSEAAQLAEIYGKLEEEIADKAPARASVILAGLGFTPKMQQPTRFSGGWRMLALARALFARPDLLL
 LDEPTNMLDVRAILWLENYLQTPSTILVSHDRNFLNAIATDIIHLHSQRLDGYRGDFETFIKSKQERL
 LNQQREYEAQQYRQHIQVFIIDRFYRANRASQVQSKLKMLEKLPVDPKSEVVLKFPDGFKEKFSPP
 ILQLDEVDFFYDPKHSIFSRLSVSADLESRICVVGENGAGKSTMLKLLMGDLSVVRGIRHAHRNLKIGYF
 SQHVEQLDLNVSARELLARKFPGLPEEYRHLGRYGISGELAMRPVAVSLGGQKSRVAFAMQTMPCPN
 FYILDEPTNHLDMETIEALGQALNNFRGGVILVSHDERFIRLVCKELWVCENGSVTRVEGGFDQYRALLQ
 EQFRREGFL

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:



ACCN: NM_013852

ORF Size: 2127 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_013852.2](#), [NP_038880.1](#)

RefSeq Size: 3288 bp

RefSeq ORF: 2130 bp

Locus ID: 27406
UniProt ID: [Q8K268](#)
Cytogenetics: 16 12.46 cM
Gene Summary: Displays an antiviral effect against flaviviruses such as west Nile virus (WNV) in the presence of OAS1B.[UniProtKB/Swiss-Prot Function]