

Product datasheet for **MG210236**

Agbl5 (BC057349) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Agbl5 (BC057349) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Agbl5
Synonyms:	4930455N08; 9430057O19Rik; CCP5
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MG210236 representing BC057349
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGAACATGAACAAGCAAAGCAAACGTATTCCAGGGCATGGCCCCCTTTGTGCGCACACTGCCTTCCC
 GGCCACGCTGGGAACGCATTAGAGAGCGGCCACCTTTGAGATGACAGAGACACAGTTTGTGCTATCCTT
 TGTTACCGTTTTCTAGAGGGCCGGGAGCCACCCTTCTTCGCTTCTGCTACCCCTTCTCTACAGT
 GACTGCCAGGATCTGCTGAGCCAGCTAGACCAGCGTTTTTCAGAGAATTACTCTACCCATAGCAGTCCCTC
 TGGATAGTATTTATTACCACCGAGAGCTTCTCTGCTATTCTCTGGATGGACTTCGTGTAGATCTACTAAC
 AATCACTTCTGCCATGGGCTTCGAGATGATCGGGAGCCCGACTAGAGCAGCTATTTCTGATCTCGGC
 ACTCCCGACCATTCCGTTTACAGGCAAAGGATATTCTTCTAAGCAGTAGGGTACACCTGGAGAGA
 CTCATCTAGCTTTGTGTTCAATGGCTTCTGGACTTCATCCTTCGACCTGATGATCCCCGGGCACAAC
 TCTCCGTCGACTTTGTTTTAAGTTGATTCCCATGTTAAACCCTGATGGTGTGGTCCGGGGCCACTAC
 CGCACGGACTCCCGTGGAGTGAATCTGAACCGTCACTACTGAAGCCTGATGCCGTCTGCACCCGGCCA
 TCTATGGCGCAAAGCTGTGCTGCTTTATCATCAGTGCCTCGAGGCTCAACGCCAAGAGTCCCACGAA
 CCAGCAACCCACTCTCCATCTCCCTCCGGAAGCCCTCTCTGACCTCGAGAAAGCCAATAATCTCCAC
 AATGAAGCTCACCTCGGGCAGTCACCTGATGGGAAAACCTGCGACTTGGCCTGAGACAGAGCCAGCAG
 AAGAGAAGACTGACCCTGTGTGGCTTATGCCACAGCCAATTCCTGAACCTCGAGGAGCCAGCCCTGACAC
 CATCCCCCAAAGAGAGTGGTGTGCTACTATGTAGACTTGCACGGACATGCTTCCAAAAGGGGCTGC
 TTCATGTATGGCAACAGCTTTAGTGACGAGAGCACCCAGGTGAAAACATGCTATATCCAAAGCTCATCT
 CCTTGAATTCAGCCACTTTGACTTCCAGGGTGCAATTTCTCAGAGAAGAACATGTATGCCCGAGACCG
 TAGAGATGGCCAGTCCAAGAGGGAAGCGCCGTGTCGCAATCTACAAAGCCTCAGGGATAATCCACAGC
 TACACACTTGAATGCAACTACAACACTGGACGCTCAGTCAACAGCATCCCTGCTGCTGCCATGACAACG
 GGCGCGCCAGCCCCCTCCCTCCGGCTTCCCTCCAGATACACTGTGGAACATTCGAGCAGGTGGG
 GCGAGCTATGGCCATCGCAGCTCTGGACATGGCAGAGTGAATCCGTGGCCCCGAATTTGTGCTGTCGGAA
 CACAGCAGCCTCACCAACCTGCGAGCCTGGATGCTGAGGCACGTGCGCAACAGCAGAGGGCTGACCAGTG
 CTGGGAACATGGGCGGAGCAAGAAGCGAGGGGCTCGAACCCCTCCAAAAGCAACAACAGCTTGCCTGT
 TTCCTGCTCTGAAAATGCTCTGAGCCGGTTTCGAGTTTTAGCACCGGCACGAGTACTGGTGGCAGCAGC
 AGCAGCCAACAAAACCTCCACAGATGAAGAACTCTCCAGCTTCCCTTTTCATGGCAGTCGGACTGCGG
 GGCTACCAGGCCTGGGCTCTAGCACCCAAAAGGTACGCCATCGGGTGTGGGCCCTGTGAGAGAGCCCCG
 ATGCTCAGACAGAAGGGCGGACAGCAGCCACTGAACCATCGCTCTACCAAGCAGCCTGGCCCCATCC
 CCAACTCTGGCTAGTTCTGGCCAACTCCTCACGAAACATGGGCTCCTGCCTCTACCCAATTCACTCA
 GCTTATCAGGGAGCAGCTGCTCATTCTCGTCTTCAGGGGACAAGCCAGAGGCTGTCATGGTGATTGGGAA
 AAGCCTACTAGGGGCTGGTGTGCTCGGATCCCTTGCATTAGGACTCGATTGCAGACTCCATTGTCCTGCCA
 GATGCACCCCAACCACAGCCTACTAGCTTTGGCTTCTACCCGTCT

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence: >MG210236 representing BC057349
 Red=Cloning site Green=Tags(s)



MNMNKQSKLYSQGMAPFVRTLPSRPRWERIRERPTFEMTETQFVLSFVHRFVEGRGATFFAFCYPFSYS
 DCQDLLSQLDQRFSENYSTHSSPLDSIYYHRELLCYSLDGLRVDLLTITSCHGLRDDREPRLEQLFPDLG
 TPRPFRFTGKRIFFLSSRVHPGETPSSFFVNGFLDFILRPDDPRAQTLRRLFVFKLIPMLNPDGVVRGHY
 RTDSRGVNLNRQYLKPDVAVLHPAIYGAVALLYHHVHSRLNAKSPTNQQTLLHLPPEAPLSDLEKANLH
 NEAHLGQSPDGENPATWPETEPAEKTDPVWLMPOPIPELEEPAPDTIPPKESGVAYYVDLHGASKRGC
 FMYGNSFSDESTQVENMLYPKLSLNSAHDFQGCNFKSEKNMYARRDRDQSQKEGSGRVAIYKASGIIHS
 YTLECNNTGRSVNSIPAACHDNGRASPPPPAFPSRYTVELFEQVGRAMAIAALDMAECNPWPRIVLS
 HSSLTNLRAWMLRHVRNSRGLTSAGNMGASKKRGARTPPKSNNSLPVSCSENALSRVRSFSTGTSTGGSS
 SSQQNSPQMKNPSFPFHGSRTAGLPGLGSSTQKVSHRVLGPVREPRCSDRRRRQQPLNHRSTTSSLAPS
 PTLASSGPTSSRNMGSCLLPNSLSLGSSCSFSSSGDKPEAVMVIKSLLGAGARIPICIRTLQTPLSLP
 DAPQPQTSFGFPTRP

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:



ACCN: BC057349

ORF Size: 2150 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [BC057349](#), [AAH57349](#)

RefSeq Size: 3332 bp

RefSeq ORF: 2150 bp

Locus ID: 231093

Cytogenetics: 5 B1

Gene Summary: Metalloprotease that mediates protein deglutamylation. Specifically catalyzes the deglutamylation of the branching point glutamate side chains generated by post-translational glutamylation in proteins such as tubulins (PubMed:20519502, PubMed:21074048). In contrast, it is not able to act as a long-chain deglutamylase that shortens long polyglutamate chains, a process catalyzed by AGTPBP1/CCP1, AGL2/CCP2, AGL3/CCP3, AGL1/CCP4 and AGL4/CCP6 (PubMed:25103237). Mediates deglutamylation of CGAS, regulating the antiviral activity of CGAS (PubMed:26829768).[UniProtKB/Swiss-Prot Function]