

Product datasheet for **MG210239**

Sec14l1 (BC005766) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Sec14l1 (BC005766) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Sec14l1
Synonyms:	1200017E04Rik; 2810012L19Rik; Naa-35
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>MG210239 representing BC005766
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGTGCAAAAATACCAGTCCCCAGTGC GG GTGTATAAACATCCCTTTGAGTTAATTATGGCTGCCTATG
 AGAGGAGGTTTCCTACATGTCCTCTGATTCCAATGTTCTGGACAGTGACACTGTGAGCGAGTTCAAGAG
 TGAGGATGGGGCTCTTCATGTCATCGAGAGGCGCTGCAAGCTGGATATAGATGCGCCAGGCTGTTGAAG
 AAGATCGCAGGAGTCGATTACGTTTATTTGTCCAGAAAACTCCCTGAATTCTCGAGATCGCACTTTGC
 ACATCGAGGCCACAATGAGACATTTTCTAACCGGTCATCATCCACGAACACTGCTGTTACACGGTTCA
 CCCTGAGAATGAAGACTGGACCTGTTTCGAACAGTCTGCAAGTTTAGATATCAAATCCTTCTTTGGTTTT
 GAAAGTACAGTGGAAAAATTGCCATGAAACATTATAAAGCAACATTAAGGAAGGAGATCATTG
 AATACTACCTACGGCAGCTGGAGGAGGAAGGCATCACCTTCGTGCCCCGCTGGACCCACCCCGTGGG
 ACCCTCAGAAACATGTTTCATCAAGCAAGAATCAAGTCACATCTGCAGCTGTCCTTGTCCCAGATGCTGCT
 GCCGTCATGGAGGGGCTGAGTGGGAAAACTCAGTAGTCCAGGAACAGCGTCTGAGCCCGTGGTGGGAA
 CCCCTGACGATAAATTAGATGCCGACTACATCAAGAGATACTTGGGTGACCTGACTCCCTCCAAGAGAG
 CTGCCTCATTGCTGCGCCAGTGGCTCCAGGAGACACAAAGGCAAGATCCCAAGGACGAGCATATC
 CTTGCGTTCTGCGTGCCGGGATTTTAAATTGACAAGGCCAGAGAGATCATGTGCCAGTCTCTAACCT
 GGCGGAAGCAGCACCAGTGGATTACATCCTGGACACTTGGACTCCGCCCCAGGTCTTCTGGACTACTA
 CGCTGGGGGCTGGCATCACCATGACAAAGATGGACGGCCGCTGTATGTGCTTAGGCTGGGCCAGATGGAC
 ACCAAAGGCTTGGTGGAGCTCTCGGAGAGGAGGCCCTACTGCGATACGTTCTTTCCATAAACGAAGAAG
 GCCTGAGACGATGCGAGGAAAACCAAGGTCTTCGGCCGGCCAATCAGCTCGTGGACCTGCTTGGTGGA
 CTTTGAAGGACTGAACATGCGCCATCTATGGAGACCAGGGGTCAAAGCCTTGCTGCGCATCATCGAGGTG
 GTGGAGGCCAACTACCCGAGACGCTGGGCCGCTCCTCATCCTCCGAGCCCCAGGGTCTTCCCTGTCC
 TCTGGACGCTGGTTAGTCCATTTATTGATGACAACACCAGAAGGAAGTTCCTCATTTATGCAGGAAATGA
 CTACCAGGGCCCAGGAGGCTGCTGGATTACATCGATAAAGAGATTATCCCGATTTCTGAGTGGGGAG
 TGCATGTGTGACGTGCCAGAGGGTGGATTGGTCCCAAATCCCTGTACAGGACTGCAGAGGAGCTGGAAA
 ATGAAGACCTGAAGCTCTGGACGGAGACCATCTACCAGTCTGCCAGCGTGTCAAAGGAGCTCCACATGA
 GATTCTCATTAGATTGTGGACGCTCCTCAGTGATCACCTGGGATTTTGTGTGTGCAAGGGGACATC
 GTCTTTAACATTTATCACTCCAAGAGGTCTCCCAGCCACCAAGAAGGACTCGCTAGGGGCTCACAGCA
 TCACTTCCCAGGAGGCAACAACGTGCAGCTGATAGACAAAGTCTGGCAGCTGGGCCGGGACTATAGCAT
 GGTGGAGTCTCCTCTAATCTGCAAAGAAGGAGAGAGCGTGCAGGGGTCCCATGTGACCAGGTGGCCAGGC
 TTCTACATCCTACAGTGGAAATTCACACCATGCCAGCATGTGCTGCCACCAACCTACCCGAGTGGACG
 ACGTGCTGGCCTCCCTGCAGGTCTCCTCCACAAGTGCAAAGTCATGTACTACACAGAAGTATCGGCTC
 TGAAGATTTAGAGGCTCCATGACCAGCCTGGAATCTAGCCACAGCGGGTTCTCACAGCTCAGTGCCGCC
 ACCACCTCCTCCAGCCAGTCGCAATCCAGCTCCATGATTTCAGTGAA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA


Protein Sequence: >MG210239 representing BC005766
 Red=Cloning site Green=Tags(s)

MVQKYQSPVRVYKHPFELIMAA YERRFPTCPLIPMFVSDTVSEFKSEDGALHVIERRCKLDIDAPRLK
 KIAGVDYVYFVQKNSLNRDRTLHIEAHNETFSNRV IIEHCCTVHPENEDWTCFEQSASLDIKSFFGF
 ESTVEKIAMKHYSNIKKGKEIIEYYLRQLEEGITFVPRWTPPPVGPSETCSSKNQVTSAAVLVPDAA
 AVMEGLSGENLSSPGTASEPVVGTPDDKLDADYIKRYLGDLTPLQESCLIRLRQWLQETHKGKIPKDEHI
 LRFLRARDFNIDKAREIMCQSLTWRKQHVDYILDWTWPPQVLLDYYAGGWHHDKDGRPLYVLRGQMD
 TKGLVRALGEEALLRYVLSINEEGLRRCENTKVFGRPISWTCCLVDLEGLNMRHLWRPGVKALLRIIEV
 VEANYPETLGRLLILRAPRVFVPLWTLVSPFIDDNTRRKFLIYAGNDYQGPGLLDYIDKEIIPDFLSGE
 CMCDVPEGGLVPKSLYRTAELENEDLKLWTETIYQSASVFKGAPHEILIQIVDASSVITWDFDVCKGDI
 VFNIYHSKRSPQPPKDSLGAHSITSPGGNNVQLIDKVWLGRDYSMVESPLICKEGESVQGSHTVRWPG
 FYILQWKFTMPACAATNLPRVDDVLASLQVSSHCKVMYYTEVIGSEDFRGSMTSLESSHSGFSQLSAA
 TTSSSQSSSMISSE

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: BC005766

ORF Size: 2148 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: BC005766, AAH05766

RefSeq Size: 2628 bp


RefSeq ORF: 2150 bp

Locus ID: 74136

Cytogenetics: 11 E2

Gene Summary: May play a role in innate immunity by inhibiting the antiviral RIG-I signaling pathway. In this pathway, functions as a negative regulator of DDX58/RIG-I, the cytoplasmic sensor of viral nucleic acids. Prevents the interaction of DDX58 with MAVS/IPS1, an important step in signal propagation. May also regulate the SLC18A3 and SLC5A7 cholinergic transporters. [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MG210239