

## Product datasheet for **MG210351**

### Cop1 (NM\_011931) Mouse Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Cop1 (NM_011931) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Cop1
Synonyms:	A1316802; C80879; Cop1; Rfwd2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide  
Sequence:

>MG210351 representing NM\_011931  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGTCAGGTAGCCGCCAGGCCGGCTCGGGTTCGCTGGGACGAGCCCCGGCTCCTCGGCCGCTCCTCGG  
TGACTTCCGCCTCCTCGTCCTTATCCTCGTCTCGGTCGCCGCCGTCCTGGCGGCCCTCGGCCGCGACGCT  
GGTGTCCGGCGCGTGGCTCCGGCCGCGGGCTCCGGCGGCCCTCGCGGCCCGGGCGGCCTGTGCTGGT  
GCGGCCGCGGTGTCTGGCAGCGCGAGCGCCGGCGGGCGGTGTCCGGGCCAGTCCCGGCTCAGCTGCG  
CGGCCAGGCCAGCGCCGGCTTGGAGGGAGCAGCTCCAGCCTCGGCAGTAGCAGTAGGAAGCGACCTCT  
GCTCGTCCGCTTTGTAACGGGCTCCTCAACTCCTACGAGGACAAAAGCAACGACTTCGTCTGTCCATC  
TGCTTTGATATGATTGAGGAGGCATACATGACAAAATGTGGCCACAGCTTTTGTACAAGTGATTATC  
AGAGTTTGGAGGACAATAATAGATGTCCCAAGTGAATTATGTTGTGGACAATATTGACCATCTCTATCC  
TAATTTCTTGGTGAATGAACTATTCTCAAACAGAAGCAAAGATTTGAGGAAAAGAGTTCAAATGGAC  
CACTCAGTGAGTAGCACCAATGGGCATAGGTGGCAAAATTTCAAGATCTACTAGGAACTGATCAAGATA  
ACCTTGATTTGGCCAATGTCAACCTCATGTTGGAATTATTAGTGCAAGAAGAAAACAACCTGGAAGCAGA  
ATCACATGCAGCTCAGCTACAGATCCTTATGGAATTCCTCAAGTTGCAAGGAGAAAATAAGAGAGAGCAA  
TTGGAGCAGATCCAGAAGGAACTAAGTGTGTTTGGAAAGAGGATATTAAGAGTGGAAAGAAATGAGTGGCC  
TCTACTCTCCTGTGAGTGAAGATAGCACAGTGCCTCAATTTGAAGCTCCTTCTCCATCACACAGTAGCAT  
TATTGATTCTACAGAATACAGCCAACCTCCAGGTTTCAAGTGGGACCTCTCAGACAAAAGAACAGCCTTGG  
TATAACAGCACATTAGCATCAAGACGAAAGCGACTCACTGCTCATTGAAAGACTTAGAGCAGTGTATT  
TTTCTACAAGGATGCTCTGATCTCAGATGACAGTGAAGTGAAGTGAAGTGAAGTGAAGTGAAGTGAAGT  
CTTGTCCAAGTTTACTCGATACAACTCAGTAAGACCGTTGGCCACATTGTCTTATGCTAGTGATCTCTAT  
AATGGTCCAGCATAGTTTCTAGTATTGAGTTTACCAGGATTGTGACTATTTTGAATTGCTGGAGTTA  
CAAAGAAGATTAAGTCTATGAGTATGGCACAGTCATCCAGGATGCAGTGGATATTACCTACCCTGAGAA  
TGAAATGACCTGCAATTCAAAAATCAGCTGTATCAGTTGGAGTAGTTACCATAAAGAACCTGCTAGCCAGC  
AGTGATTATGAAGGCACTGTTATACTATGGGATGGATTACAGGACAGAGGTCAAAGGTCTATCAGGAGC  
ATGAAAAAGGTGTTGGAGTGTGACTTTAATCTGATGGATCCTAACTTCTGGCTTCAAGTTCTGATGA  
TGCAAAAGTGAAGCTGTGGTCTACCAATTTAGACAATTTCTGTGGCGAGCATTGAGGCAAAGGCTAACGTG  
TGCTGTGTGAAGTTCAGCCCTCCTCCAGGTACCATCTGGCTTTCGGCTGTGCAGATCACTGTGCCACT  
ACTACGACCTTCGTAACACTAAACAGCCAATAATGGTATTCAAGGGACACCGAAAAGCAGTGTCTTATGC  
CAAGTTTGAAGTGGTGAAGAAATGTCTCAGCCTCAACAGACAGCCAACCTAAAGCTGTGGAATGTGGGA  
AAACCATACTGTCTACGTTCTTCAAGGTCATATTAATGAAAAAACTTTGTAGGCCTTGTCTCAATG  
GAGATTATAGCGTGTGGAAGTGAACAACCTCCCTCTACCTGTATTATAAAGGACTTTCTAAAACCTT  
GCTAACTTTAAGTTTGTATACCGTTAAGAGTGTATTAGACAAAGATCGGAAAGAAGATGACACAAATGAA  
TTTGTGCTGCTGTGTGTTGGCGGGCACTATCAGACGGGAGTCCAATGTGCTGATTGCTGCTAACAGCC  
AGGGTACAATTAAGGTGCTAGAATTGGTA

**ACGCGT**ACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence: >MG210351 representing NM\_011931  
 Red=Cloning site Green=Tags(s)

```


MSGSRQAGSGSAGTSPGSSAASSVTSASSSLSSSPSPPSVAASAATLVSGGVAPAAGSGGLGGPGRPVLV
AAAVSGSASAGGAVSAGQSRLSCAARPSAGVGGSSSSLGSSSRKRPLLVPLCNGLLNSYEDKSNDFVCP
CFDMIEEAYMTKCGHSFCYKCIHQSLLEDNNRCPKCNVVDNIDHLYPNFLVNELILKQKQRFEEKRFKLD
HSVSSTNGHRWQIFQDLLGTDQDNLDLANVNLMELELVQKKKQLEAESHAACLQILMEFLKVARRNKREQ
LEQIQKELSVLEEDIKRVEEMSGLYSPVSEDSTVPQFEAPSPSHSSIIDSTEYSQPPGFGSQTSTKQPW
YNSTLASRRKRLTAHFEDLEQCYFSTRMSRISDDSRASQLDEFQECLSKFTRYNSVRPLATLSYASDLY
NGSSIVSSIEFDRDCDYFAIAGVTKKIKVYEGTVIQDAVDIHYPENEMTCNSKISCISWSSYHKNNLLAS
SDYEGTVILWDGFTGQRSKVYQEHEKRCWSVDFNLMDPKLLASGSDDAKVKLWSTNLDNSVASIEAKANV
CCVKFSPSSRYHLAFGCADHCVHYDRLNRTKQPIIMVFKGHRKAVSYAKFVSGEEIVSASTDSQLKLVNNG
KPYCLRSFKGHINEKNFVGLASNGDYIACGSENNSLYLYYKGLSKTLLTFKFDTVKSVLDKDRKEDDTE
FVSAVCWRALSDGESNLIAANSQGTIKVLELV
  
```

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:



ACCN: NM\_011931

ORF Size: 2199 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_011931.4](#)

RefSeq Size: 5056 bp

RefSeq ORF: 2202 bp

Locus ID: 26374

UniProt ID: [Q9R1A8](#)

Cytogenetics: 1 H1

**Gene Summary:** E3 ubiquitin-protein ligase that mediates ubiquitination and subsequent proteasomal degradation of target proteins. E3 ubiquitin ligases accept ubiquitin from an E2 ubiquitin-conjugating enzyme in the form of a thioester and then directly transfers the ubiquitin to targeted substrates. Involved in JUN ubiquitination and degradation. Directly involved in p53 (TP53) ubiquitination and degradation, thereby abolishing p53-dependent transcription and apoptosis. Ubiquitinates p53 independently of MDM2 or RCHY1. Probably mediates E3 ubiquitin ligase activity by functioning as the essential RING domain subunit of larger E3 complexes. In contrast, it does not constitute the catalytic RING subunit in the DCX DET1-COP1 complex that negatively regulates JUN, the ubiquitin ligase activity being mediated by RBX1. Involved in 14-3-3 protein sigma/SFN ubiquitination and proteasomal degradation, leading to AKT activation and promotion of cell survival. Ubiquitinates MTA1 leading to its proteasomal degradation. Upon binding to TRIB1, ubiquitinates CEBPA, which lacks a canonical COP1-binding motif.[UniProtKB/Swiss-Prot Function]