

Product datasheet for **MG210403**

Zbtb20 (NM_019778) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Zbtb20 (NM_019778) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Zbtb20
Synonyms:	1300017A20Rik; 7330412A13Rik; A930017C21Rik; D16Wsu73e; DPZF; HOF; ODA-8S; Oda8; Zfp288
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>MG210403 representing NM_019778
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCCGCATCGCC

ATGCTAGAACGGAAGAAACCCAAGACAGCTGAAAACAGAAAGGCATCTGAGGAGAATGAGATTACTCAGC
CGGGCGGATCCAGCGCCAAGCCGGCCCTTCCCTGCCTGAACTTTGAAGCTGTTTTGTCTCCAGCCCAGC
CCTCATCCACTCGACACATTCACTGACAACTCTCACGCTCACACCGGGTATCTGATTGTGACATCAGT
TGCAAGGGGATGACCGAGCGCATTACAGCATCAACCTTCACAACCTTCAGCAATTCGGTGTCTGAGACCC
TCAACGAGCAGCGCAACCGTGGCCACTTCTGTGACGTGACGGTTCGCATCCACGGGAGCATGCTGCGCGC
ACATCGTGTGCTGGCAGCCGGCAGCCCTTCTTCCAAGACAAGCTGTGCTGGGCTACAGCGACATC
GAAATCCCGTCCGTGGTGTCCGTACAATCGGTGCAAAAGCTCATTGACTTCATGTACAGCGGTGTGCTGA
GAGTCTCACAGTCGGAAGCTCTGCAGATCCTCACAGCCGCCAGCATCCTGCAGATCAAAACAGTCATAGA
TGAGTGCCTCGCATCGTGTACAGAACGTGGGCGATGTGTTCCAGGCATCCAGGATTCGGCCAGGAC
ACACCAAGAGGCACACCAGAGTCAGGCACATCTGGCCAGAGCAGTGACACGGAATCAGGCTACCTGCAGA
GCCACCCACAGCATAGTGTGGACCGAATCTACTCCGCACTCTACGCCCTGCTCCATGCAGAATGGCAGCGG
CGAGCGCTCCTTCTACAGTGGTGCAGTGGTCAGCCACCACGAAACAGCTCTCGGCCCTGCCCGTGACCAC
CACATGGAAGACCCTAGCTGGATCACACGCATTATGAGCGCTCCAGCAAAATGGAGCGCTACCTGTCCA
CCACCCCTGAGACCACGCACTGCCGGAAGCAGCCCGGCCCTGTGCGTATCCAGACCCTGGTGGGTAAACAT
CCACATCAAGCAGGAAATGGAAGATGACTATGACTACTATGGGCAGCAAAGGGTGCAGATCCTAGAACGC
AATGAATCCGAGGAGTGACACAGAAGACTGACCAAGCAGAGGGCACTGAGAGCGAGCCAAAGGTGAAA
GCTTTGATTCTGGGGTCAGCTCCTCCATCGGCACCGAACCTGACTCAGTGGAGCAACAGTTTGGGCGAG
AGCCCCAAGGGACGGTACGGCAGAACCAGCCCAACCTGAGCAGGCAGCAGAAGCCCCAGCTGAGAGCAGT
GCCAGCCAAACAGCTAGAACCAGGTGCCCTCCTCCTGAGAGAAGCAACGAGTCAGAGATGGACAACA
CAGTCATCACTGTGAGTAACAGCTCCGATAAGGGCGTCTACAGCAGCCTTCAGTCAACACATCCATCGG
GCAGCCATTGCCAAGTACCCAGCTCTATTTACGCCAGACAGAAACCCTCACAGCAACCTGAGGATGCCT
CTGACCTTGACCAGCAACACACAGGTATTGGCACCGCTGGCAACACCTATCTGCCAGCCCTTTCACTA
CCCAACCCGCGGGCAGTGGCCCCAAGCCTTTTCTTTTCTTTCAGCCTGCCGAGCCCTGACAGGCCAGCAGC
CCAGTTTGTGACAGTGTCCAGCCCGGTCTGTCCACCTTACTGCACAGCTGCCAGCGCCACAGCCCTG
GCCTCATCTGCAGGCCACAGCACAGCCAGTGGGCAAGGGCAGAAAAAGCCTTATGAGTGCCTCTCTGCA
ACAAGACTTTCACAGCCAAACAGAACTACGTCAAGCACATGTTCCGTACATACAGGTGAGAAGCCCCACCA
GTGCAGCATCTGCTGGCGCTCCTTCTCCTTGAAGGATTACCTTATCAAGCACATGGTGACGCACACCGGC
GTGAGAGCGTACCAGTGTAGCATCTGCAACAAGCGCTTCAACCAGAAAGAGTTCCCTCAACGTGCACATGC
GCCTGCACCGCGGGGAGAAGTCCATGAGTGTACATCTGCAAAAAGAAGTTCTCCACAAAGACCCTGCT
GGAGCGACAGTGGCCCTGCACAGTGCCAGCAACGGGACCCCTCCGGCAGGCACGCCCCAGGTGCCCGC
GCGGGTCCGCCAGGCGTGGTGGCTGCACAGAGGGGACCCTTACGTCTGCTCCGTCTGCCAGCAAAGT
TTGACCAAATCGAGCAGTTCACGACCACATGAGGATGCATGTGTCTGACGGA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG210403 representing NM_019778
 Red=Cloning site Green=Tags(s)

MLERKKPKTAENQKASEENEITQPGSSAKPALPCLNFEAVLSPAPALIHSTHSLTNSHAHTGSSDCDIS
 CKGMTERIHSINLHNFNSVLETLNEQRNRGHFCDVTVRIHGSMLRAHRCVLAAGSPFFQDKLLLGYS
 DI EIPSVSVQSVQKLIIDFMYSGVLRVSQSEALQILTAASILQIKTVIDECTRIVSQNVGDVFP
 GIQDQSGQDTPRGTPESGTSGQSSDTESGYLQSHPQHSVDRIYSALYACSMQNGSGERSFYSGAVV
 SHHETALGLPRDHMEDPSWITRIHERSQQMERYLSTTPETTHCRKQPRPVRIQTLVGNIIHIKQ
 EMEDDYDYYGQQRVQILERNESECTEDTDQAEGETESEPKGESFDSGVSSISIGTEPDSVEQQF
 GAAAPRDGQAEPAPQEQAEEAPAESSAQPNQLEPGASSPERSNESEMDNTVITVNSSSDKGVL
 QQPSVNTSIGQPLPSTQLYLRQTETLTSNLRMPLT LTSNTQVIGTAGNTYLPALFTTQPAGS
 GPKPFLFSLPQPLTGQQTQFVTVSQPLSTFTAQLPAPQPLASSAGHSTASGQGDKKPYE
 CTLCKNTFTAKQNYVKHMFVHTGEKPHQCSICWRSFSLKDYLKHMVTHGTGVRAYQCSICN
 KRFTQKSSLNVHMRLHRGEKSYECYICKKKFSHKTLLEHVALHSASNGTPPAGTPPGARAG
 PPGVVACTEGTTYVCSVCPAKFDQIEQFNDHMRMHVSDG


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:


ACCN: NM_019778

ORF Size: 3153 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_019778.2](#), [NP_062752.2](#)

RefSeq Size: 2990 bp

RefSeq ORF: 2226 bp

Locus ID: 56490

UniProt ID: [Q8K0L9](#)

Cytogenetics: 16 28.44 cM

Gene Summary: May be a transcription factor that may be involved in hematopoiesis, oncogenesis, and immune responses (By similarity). Plays a role in postnatal myogenesis, may be involved in the regulation of satellite cells self-renewal (PubMed:27446912).[UniProtKB/Swiss-Prot Function]