

## Product datasheet for **MG210426**

### Med25 (NM\_029365) Mouse Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Med25 (NM_029365) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Med25
Synonyms:	2610034E13Rik; 2610529E18Rik; ESTM2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)



Protein Sequence: >MG210426 representing NM\_029365  
 Red=Cloning site Green=Tags(s)

MVPGSEGPARGGLVADVVFVIEGTANLGPYFEELRKHYLLPAIEYFNGGPPAETDFGGDYGGTQYSLVV  
 FNTVDCAPEYVQCHAPTSSAYEFVTWLDGIKFMGGGGESCSLIAEGLSTALQLFDDFKKMREQIGQTHR  
 VCLLICNSPPYLLPAVESTTYSGCTTESLVQKIGERGIHF SIVSPRKL PALRLLFEKAAPPALLEPLQQP  
 ADVSQDPRHMLVVRGLVLPVGGSSSTSGSLQTKQAVPLPPAPASAATLSAAPPQALPPVPPQYQVPGNL SA  
 AQVAAQNAVEAAKSQKAGLGPRFSPINPLQQAAPGVGPPFSQAPAPPLAPVPPGAPKPPPASQPSLVSTV  
 APGPVLAAPAQPGAPSLAGTVTPGGVNGPSAAQLGGPALGGQSVSNKLLAWSGVLEWQEKPKPASVDAN  
 TKLTRSLPCQVYVNHGENLKTEQWPQKLIMQLIPQQLLTTLGPLFRNSRMVQFHF TNKDLESKGLYRIM  
 GNGFAGCVHFPHTAPCEVRVLM LLYSSKKKIFMGLIPYDQSGFVNGIRQVITNHKQVQQQKLEQQRGMGA  
 QQAPPVLGP ILEEQARPPQNL LQLRAPQPQPGAVGASAAATGQPQPGATQAPT GAPQGGPAAAGPPPS  
 GPILRPQNPGANPQLRLLLLNPAPPQTGVPPPQASLHHLQPPGAPTL LPPHQSMGQPQLGPQLLHPPPAQ  
 SWPTQLPQRAPLPQMLLSGGPRGPVQPGLQPSVMEDDILMDLI


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



**Plasmid Map:**


**ACCN:** NM\_029365

**ORF Size:** 2235 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_029365.3](#)

**RefSeq Size:** 2326 bp

**RefSeq ORF:** 2238 bp

**Locus ID:** 75613

UniProt ID: [Q8VCB2](#)

Cytogenetics: 7 28.95 cM

**Gene Summary:** Component of the Mediator complex, a coactivator involved in the regulated transcription of nearly all RNA polymerase II-dependent genes. Mediator functions as a bridge to convey information from gene-specific regulatory proteins to the basal RNA polymerase II transcription machinery. Mediator is recruited to promoters by direct interactions with regulatory proteins and serves as a scaffold for the assembly of a functional preinitiation complex with RNA polymerase II and the general transcription factors. Required for RARA/RXRA-mediated transcription (By similarity).[UniProtKB/Swiss-Prot Function]