

Product datasheet for **MG210524**

Abcb9 (NM_019875) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Abcb9 (NM_019875) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Abcb9
Synonyms:	mKIAA1520; TAPL
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MG210524 representing NM_019875
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGAGGCTGTGGAAGCGGTAGTAGTGACGCTGGCCTTCGTGAGCACGGACGTCGGGGTACCACCGCCA
 TCTACGCCTTCAGCCACCTGGACCGTAGCCTGCTGGAGGACATCCGACACTTTAACATCTTCGACTCGGT
 GCTGGACCTCTGGGCTGCCTGCCTGTACCGCAGCTGCCTGCTGCTGGGGCTACCATAGGCGTGGCCAAG
 AACAGCGCTCTGGGGCCAAGGCGCTGCGGGCTTCGTGGCTCGTCATCACCTCGTGTGCCTCTTTGTGG
 GCATCTATGCCATGGCCAACTGCTACTCTTCTCAGAAGTGCAGAGGCCATCCGGGACCATGTTCTG
 GCGCTCTTCGTGTGGACCTACATCTCGCTGGCTGCCTCCTTCTGCTGTGGGGCTGCTGGCCACTGTG
 CGGCCGAGCAGAGGCCCTGGAGCCAGGAATGAGGGCTCCACGGGGAGGGTGGGGCCCTGCTGAGC
 AGGCATCGGGAGCCACACTACAGAAGCTGTTGTCCTACACAAAGCCTGACGTGGCCTTCTAGTAGCTGC
 CTCCTTTTCTCATCGTAGCAGCTCTGGGTGAGACCTTCTGCCCTACTACACTGGCCGGGCCATTGAC
 AGCATTGTGATCCAGAAAAGCATGGATCAGTTCACCACAGCCGTTGTCTGCTGTTGCCTGCTGGCCATCG
 GCAGCTCATTGGCCGAGGTATTCGGGGCGGTATTTTACCCTCGTATTTGCCAGACTGAACATTCGCCCT
 CCGAACTGTCTTCCGCTCACTGGTGTACAGGAGACAAGCTTCTTTGATGAGAACCACAGGGGAC
 CTCATCTCCCGCTCACCTCTGATACCACCATGGTCAGCGACCTGGTCTCCAGAACATCAACATCTTCC
 TGAGGAACACGGTCAAAGTACGGGGCGTGGTGGTCTTCATGTTACGCCTCTCCTGGCAGCTTCCCTGGT
 CACCTTATGGGCTTCCCTATCATGATGGTGTCCAACATCTACGGCAAGTACTACAAGAGGCTCTCC
 AAGGAGTCCAGAGCGCCTGGCCAGAGCCAGCACAGCCGAGGAAACCATCAGCGCATGAAGACGG
 TGCCGAGCTTTGCCAACGAGGAGGAGGCGGAGGTGTTCTCGGGAAGCTTCAGCAGGTGACAAGCT
 CAACAGGAAGGAAGCAGCGCCTACATGTCTAGCTCTGGGGCAGCGGGCTCACACTCCTGGTGGTCCAG
 GTCAGTATTCTACTACGGGGCCACCTAGTCATCTCGGGCAGATGAGCAGCGCAACCTCATCGCCT
 TCATCATCTACGAATTTGCTCTGGGAGACTGCATGGAGTCCGTGGGCTCCGTCTACAGTGGCCTGATGCA
 GGGAGTGGGGCCCGGAGAAGGTGTTGAGTTCATCGACCGGCAGCCACCATGGTGCATGACGGAAGC
 TTGGCCCCGACCACCTCGAGGGCAGGGTGGACTTTGAGAATGTGACCTTACCTACCGCACTCGGCCCC
 ACACACAGGTCTACAGAATGTCTCTTCAGCCTGTCCCCGGCAAGGTGACAGCTCTGGTGGGGCCGTC
 GGGCAGCGGAAGAGCTCGTGTGTGAACATCCTGGAGAATTCTACCCTCTGCAGGGCGGCCGTGATTG
 TTGGATGGCAAGCCATTGGCGCTATGACCACAAGTACCTGCATCGTGTGATCTACTGGTGGCCAGG
 AGCCCGTGTGTCGCCGCTCCATCACAGACAACATCTCCTACGGCCTGCCTACCGTGCCTTCGAGAT
 GGTGGTGGAGGCTGCACAGAAGGCAATGCCATGGCTTCATCATGGAATGCAGGACGGGTACAGCACA
 GAGACCGGGGAGAAGGGAGCCAGCTGTGAGTGGGAGGAGCAGCGGGTGGCCATGGCCCGAGCACTTG
 TGCGGAACCCACCTGTGCTTATCCTGGACGAAGCCACAGTGCCTGGATGCCGAGAGCGAATACCTGAT
 TCAGCAGGCCATCCACGGCAACCTGCAGAGACACAGGTGCTCATATTGCACACCGGCTGAGCACCGTG
 GAGCGGGCACCTCATCGTGGTGTGGACAAGGGCCGGTGGTACAGCAGGGCACACACCAGCAGCTGT
 TGGCACAGGGCGGCTCTATGCCAAGCTGGTGCAGCGTCAGATGTTGGGGCTCGAGCACCCCTTGGACTA
 CACAGCCAGCCACAAGGAGCCACCCAGCAACACTGAACACAAGGCC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG210524 representing NM_019875
 Red=Cloning site Green=Tags(s)


MRLWKAVVTLAFVSTDVGVTTAIYAFSHLDRSLLEDIRHFNIFDSVLDLWAACLYRSCLLLGATIGVAK
 NSALGPRRLRASWLVITLVCLFVGIYAMAKLLLFSEVRRPIRDPWFALFVWTYISLAASFLLWGLLATV
 RPDAAEALPGNEGFHGEAGAPAEQASGATLQKLLSYTKPDVAFVVAASFLLIVAALGETFLPYTGRAID
 SIVIQKSMQDQFTTAVVVVCLLAIGSSLAAGIRGGIFTLVFARLNIRLRNCLFRSLVSQETSFFDENRTGD
 LISRLTSDTTMVSDDLVSQINIFLRNTVKVTGVVVMFSLSWQLSLVTFMGFPIIMMVSNIYGKYYKRLS
 KEVQSALARASTTAEETISAMKTVRSFANEEEEAEVFLRKLQVYKLNKREAAAAYMSYVWGSGLTLLVVQ
 VSILYYGGHLVISGQMSSGNLIAFIIYEFVLGDCMESVGSVYSGLMQGVGAAEKVFEFIDRQPTMVHDGS
 LAPDHLEGRVDFENVFTFYRTRPHTQVLQNVFSLSLSPGKVTALVGPSSGKSSCVNILENFYPLQGGRVL
 LDGKPIGAYDCHKYLHRVISLSVQEPVLFARSITDNISYGLPTVPFEMVVEAAQKANAHGFIMELQDGYST
 ETGEKGAQLSGGQKQRVAMARALVRNPPVLIIDEATSALDAESEYLIQQAIHGNLQRHTVLIIAHRLSTV
 ERAHLIVVLDKGRVVQQGTHQQLLAQGGLYAKLVQRQMLGLEHPLDYTASHKEPPSNTTEHKA


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:


ACCN: NM_019875

ORF Size: 2286 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_019875.2](#), [NP_063928.2](#)

RefSeq Size: 3334 bp

RefSeq ORF: 2289 bp

Locus ID: 56325

UniProt ID: [Q9JJ59](#)

Cytogenetics: 5 F

Gene Summary: The membrane-associated protein encoded by this gene is a member of the superfamily of ATP-binding cassette (ABC) transporters. ABC proteins transport various molecules across extra- and intra-cellular membranes. ABC genes are divided into seven distinct subfamilies (ABC1, MDR/TAP, MRP, ALD, OABP, GCN20, White). This protein is a member of the MDR/TAP subfamily. Members of the MDR/TAP subfamily are involved in multidrug resistance as well as antigen presentation. The function of this half-transporter has not yet been determined; however, it may be associated with lysosome activity. [provided by RefSeq, Jul 2008]