

## Product datasheet for **MG210533**

### Satb1 (NM\_009122) Mouse Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Satb1 (NM_009122) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Satb1
Synonyms:	2610306G12Rik; AW413156
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide Sequence:**

>MG210533 representing NM\_009122  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**GCGATCGCC**

ATGGATCATTGAAACGAGCAACTCAGGGGAAAGAACATTGAGAAATGTCTAACAAATGTGAGTGATCCGA  
 AGGGTCCACCCGCAAGATTGCCCGCCTGGAGCAGAACGGGAGCCCTCTAGGAAGAGGAAGGCTTGGGAG  
 CACAGGTGAAAGATGCAGGGAGTGCCTTTAAACACTCGGGCCATCTCATGAAAACCAACCTTAGGAAA  
 GGAACCATGTTACCAGTTTTCTGCGTGGTGAACATTATGAAAACGCCATTGAGTATGATTGCAAGGAGG  
 AGCACGCGGAATTTGATTGGTGAGAAAGGATATGCTTTTCAACCAGCTGATAGAGATGGCGTTGCTGTC  
 TCTAGGCTATTACACAGCTCTGCTGCCAAGCCAAAGGGCTCATCCAGGTTGGGAAGTGGAAATCCAGTT  
 CCACTGTCGATGTGACAGATGCCCTGATGCCACGGTGCAGACATGCTTCAAGATGTGTATCATGTGG  
 TCACCCTCAAAATTCAGTTACACAGTTGCCCTAACTAGAAGACTTGCCTCCTGAACAATGGTCGCACAC  
 CACAGTAAGGAATGCTCTGAAGGACTTACTGAAAGATATGAACCAGAGTTCGTTGGCCAAGGAGTGCCCC  
 CTTTCACAGAGCATGATCTCCTCCATTGTGAACAGCACGTAATGCAAAATGTCTCAGCAGCAAAATGTC  
 AAGAATTTGGAAGGTGGTACAAACATTTCAAGAAGACAAAGGATATGATGGTTGAGATGGATAGTCTGTC  
 TGAATATCCCAGCAAGGTGCCAACACGTCATTTTGGCCAGCAGCCTGTCCCAGGAAACACAGCTGAG  
 CAGCCTCCATCCCTGCGCAGCTCTCTCACGGCAGTCAGCCCTCTGTCCGGACCCCTCTTCCGAACCTGC  
 ACCCTGGGCTTGTGTCAACACCGATCAGTCCTCAGCTGGTCAACCAACAGCTGGTATGGCTCAGTTGCT  
 GAACCAGCAGTATGCAGTGAACAGACTCTTAGCCAGCAGTCTTAAACCAACAGTACTTGAACCACCT  
 CCCCTGTGAGTGGTCTATGAACAAGCCTTTGGAGCAGCAAGTTCCACAAACACGGAGGTCTCTTCTG  
 AAATCTACCAGTGGGTGCGGGATGAACTGAAACGAGCCGGAATCTCACAGGCAGTATTTGCACGCGTGGC  
 TTTTAACCGAACTCAGGGATTGCTTTCTGAAATCCTCCGAAAGGAAGAGGACCCCAAGACTGCATCCAG  
 TCTCTGCTGGTAAACCTTCGGGCTATGCAGAATTTCTTACAGTTGCCGGAAGCCGAAAGAGACCGGATAT  
 ACCAGGATGAGAGGGAAAGGAGCTTGAACGCAGCCTCAGCCATGGGTCTGCCCGCTGCTGAGCACACC  
 ACCCAGCCGCCCTCCCAGGTGAAAACAGCTACCCTTGCCACTGAGAGAAATGGGAAGCCAGAGAACAAT  
 ACTATGAACATTAATGCCTCCATTTATGATGAGATTGAGCAGGAAATGAAGCGTGCTAAAGTGTCCCAAG  
 CACTGTTTGCAAAGTTGCCGCCACAAAAGCCAGGGATGGCTGTGTGAGCTGTTGCGCTGGAAAGAAGA  
 TCCTTCTCCAGAAAACAGGACCCTGTGGGAGAACCTGTCGATGATCCGAAGATTCCTCAGTCTGCCTCAG  
 CCTGAGCGGATGCCATCTATGAGCAGGAGAGCAATGCTGTGCATCACCATGGCGACAGACCTCCCACA  
 TCATCCAGTTCCAGCAGAACAGATTGAGCAGCAGCAACAGCAGCAGCAACAGCAGCAGCAGCAGCA  
 GCCACCGCCGCCACCGCCGACAGCCACAGCCACAGCCCGAGCCAGGCCCCAGACTCCCCACGGCAGCCC  
 ACCGTGGCCTCCTCCGCGGAGTCCGATGAGGAAAACCGGCAGAAAGACCAGGCCACGAACCAAAATTTCCG  
 TGGAAAGCCTGGGGATCCTTACAGATTTTCAAGATGTGGCCTGTACCCAGATGAAGAGGCTATCCA  
 GACTCTGTCTGCACAGCTGGACCTCCCGAAGTACACCATCATCAAGTCTTTTCAAGACCAGCGTACTAC  
 CTTAAGCACCATGGCAAGCTGAAGGACAACCTCCGGCTTGGAGGTGGATGTGGCCGAGTCAAAGATGAGG  
 AGTTGCTTAAGGATTTGGAAGAGAGCGTCCAGGATAAAAACGCCAACACCTTTTCTCAGTGAACCTAGA  
 GGAAGAGCTGTCGGTGGAAAGGAGCACAGACGTTAATGCCGACTTGAAAGAC

AG**CGGACCC**ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA


**Protein Sequence:** >MG210533 representing NM\_009122  
Red=Cloning site Green=Tags(s)

```
MDHLNEATQGKEHSEMSNNVSDPKGPPAKIARLEQNGSPLGRGRLGSTGGKMQGVPLKHSGLMKTNLRK
GTMLPVFCVVEHYENAI EYDCKEEHAEFVLVRKDMLFNQLIEMALLSLGYSHSSAAQAKGLIQVGKWNPV
PLSYVTDAPDATVADMLQDVYHVVTLKIQLHSCPKLEDLPPEQWSHTTVRNALKDLLKDMNQSSLAKECP
LSQSMISSIVNSTYYANVSAAKCQEFGRWYKHFKKTKDMMVEMDSLSEL SQQGАНHVNFQQPVPGNТАЕ
QPPSPAQLSHGSQPSVRTPLPNLHPGLVSTPISPQLVNQQLVMAQLLNQQYAVNRLLAQQSLNQQYLНP
PPVSRSMNKPLEQQVSTNTEVSSEIYQWVRDELKRAGISQAVFARVAFNRTQGLLSEILRKEEDPKTASQ
SLLVNLRAMQNFLQLPEAERDRIYQDERERSLNAASAMGPAPLLSTPPSRPPQVKATLATERNKGPENN
TMNINASIYDEIQQEMKRAKVSQALFAKVAATKSQGWLCCELLRWKEDPSPENRTLWENLSMIRRFSLPQ
PERDAIYEQESNAVHHHGDRPPHIHVPAEQIQQQQQQQQQQQQQPPPPPPQPQPQAGPRLPPRQP
TVASSAESDEENRQKTRPRTKISVEALGILQSF IQDVGLYPDEEAIQTL SAQLDLPKYTI IKFFQNRYY
LKHHGKLDKNSGLEVDVAEYKDEELLKDLEESVQDKNANTLFSVKLEEELSVEGSTDVNADLKD
```

SGPTRRRLE - GFP Tag - V

**Restriction Sites:** SgfI-RsrII

**Cloning Scheme:**



**ACCN:** NM\_009122

**ORF Size:** 2292 bp


**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_009122.2</a> , <a href="#">NP_033148.2</a>
<b>RefSeq Size:</b>	6228 bp
<b>RefSeq ORF:</b>	2295 bp
<b>Locus ID:</b>	20230
<b>UniProt ID:</b>	<a href="#">Q60611</a>
<b>Cytogenetics:</b>	17 C
<b>Gene Summary:</b>	Required for the switching of fetal globin species, and beta- and gamma-globin genes regulation during erythroid differentiation. Plays a role in chromatin organization and nuclear architecture during apoptosis (By similarity). Crucial silencing factor contributing to the initiation of X inactivation mediated by Xist RNA that occurs during embryogenesis and in lymphoma. Binds to DNA at special AT-rich sequences, the consensus SATB1-binding sequence (CSBS), at nuclear matrix- or scaffold-associated regions. Thought to recognize the sugar-phosphate structure of double-stranded DNA. Transcriptional repressor controlling nuclear and viral gene expression in a phosphorylated and acetylated status-dependent manner, by binding to matrix attachment regions (MARs) of DNA and inducing a local chromatin-loop remodeling. Acts as a docking site for several chromatin remodeling enzymes and also by recruiting corepressors (HDACs) or coactivators (HATs) directly to promoters and enhancers. Modulates genes that are essential in the maturation of the immune T-cell CD8SP from thymocytes.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MG210533