

Product datasheet for **MG210579**

6430527G18Rik (BC063253) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	6430527G18Rik (BC063253) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	6430527G18Rik
Synonyms:	6430527G18Rik; Eap1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MG210579 representing BC063253
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGTCGGCGGCGCAGGTGTCCTCGTCCCGAAGACAGTCTTGCTACCTGTGCGACCTGCCTCGCATGCCCT
 GGGCCATGATCTGGGACTTCTCCGAACCGTATGCCGTGGTTGCGTCAATTACGAAGGCGCCGACCGCAT
 AGAATTCGTGATCGAGACAGCACGCCAGCTGAAGCGGGCTCACGGTTGCTTCCAAGACGGTCGCTCCCCG
 GGGCCCGCCCGCCAGTCGGGGTCAAGACGGTGGCCTTGTCCGCCAAGGAAGCGGCAGCAGCAGCTGCCG
 CCGCACAGCAACAGCAGCAGCAACAACAGCAGCAACAGCAGCAGCTCAACCACGTCGATGGTTCCACCAA
 ACCTGCAGTGTGGCCGCCCGTCCGGTCTGGAGCGTTATGGCCTGAGCGCAGCGCCGCCCGCGGTGCC
 GCCCGCGTGGTGAACAACGACGCCGGTTCGAGTATCCGCCTCCTCCGGTGGTCTGGGGAGCAGCA
 GCCACGCCCGCGCTGCCAACGGCTGGGAGGTCTAACGGCTTCCCCAAGCCGGCACCAGAGGAAGG
 CCCCCCGAGCTGAACCGCCAGAGCCCCAACTCATCTCAGCAGCAACATCGGTGGCTTCTCGCGTGGG
 ACGCATAGTGGTGGTACCGGGCTGCCAACCCGGGTGGTGGCGGGGACCTCAGTAACCGTGGCCC
 CCAACCTATTGCCGACAGCGCTGCTAACCGGCCGGCCAGCGCCGCGGTGCTGCCTCCGCCCCACGGTT
 GGGGGGAAGCCGAGGGCCCCGACTCCAGCACCTCCGGGTGCTCCCGGGGGCCCCGCATGTCTCGGGGT
 CCCCCGGTGTGTAGCCACTGTATCCTCCGCTCCATCATCGACTTCGTCGACGGTGGCAGAGGTGGGGC
 TGGGTGCTGTGGCAAGAGGCTGGATCGGTGTGAGCAGCGACCGAGGAGCGTAAAGGAGAAGCA
 GCGCAACGCTGAGGCCTGGCAGAGTGTGAGCAGAGCCTGCGCAACCGCGCGGAGGAGTGGGCCAACAA
 CCCAAGATGGTCCGCGACAGCTGTCTCAGCTGGCTGGCTGCACGCCCTACGAGGTGCGCTTCAAGAAG
 ACCACTCGCTTTGGGTGCGCTTTGGCTTCGACGCGGTCTTAAGCCGGGCATGGACTATGAGTGA
 GCTGTTCAATTGAGTATCCACGGGCTCCGGCAATGTGTACTCCAGCGCGTCTGGGGTGGCCAAACAGATG
 TACCAGGATTGCATGAAGGACTTTGGCCGGGTGTGCTCTGGCTTCAAGTACCTGGAGTACGAGAAGA
 AGCACGGCTGTGGCGACTGGCGCTACTTGGGACCTACTGCCTGAGGCTGTGCGCTTCTCAAGGAGGG
 CGTGCCCCGGCGGACATGCTGCCCCAGCCCTACCTAGATGCTAGCTGTCCCATGCTACCCACGGCGCTG
 GTGAGTCTCAGCCGCGCCCCAGCGCACCTCCAGGAAGTGGAGCCTTCCACCCGCGCGCCACCGGGCC
 GGGGTGCAGCGTCCAGCCTGCGAAAAGGAAGGCATCCCCGGAGCCTCCAGACTCAGCTGAAAGCGCACT
 GAAGCTCGGTGAGGAACAGCAGAGGCAGCAGTGGATGGCAAACCAGAGCGAGGCGCTGAAGCTCACGATG
 TCCGCGGGGGCTTTGCGGCGCCGGGACTCAGCAGGGGACCACCTCCGCCCCCTCCACCTCTGGGAC
 CTATTCCAATCGGACCACCCCGCGGAGTCAGCCCCCAGAACGGTCCGTCCCCCATGGCTGCCCTCAT
 GTCGGTAGCAGACACTTTGGGCACAGCTCATTGCCCAAGGACGGCAGTTCAAGTGCCTCGACACTGCA
 TCTGCTAGGCGAAACAGCAGCAGCCGGTGTACCCCGCCTCCGTGCCGGGCCAGCGTGCCTGGCATCAC
 GCAATGGGGACCTGAATTTACAGGTGGCGCCCCCGCCCTAGCGCCACCCGGGCATGGACCAAGTACA
 CCCCCAAAACATCCGGATTCCTATGGCCAACAGCGGACCCCTCTGCTGTACCATTTGCCACGAACGT
 TTGGAGGACACGCATTTTCGTTTCAGTGTCCGTCTGTCCCTAGCCACAAATTTGTTCCCTTGTCTAGAG
 AGAGTATCAAGGCCAGGGGGCCACAGGTGAGGTGATTGCCCCAGCGGAGAGAAATGCCCCCTAGTCGG
 GTCTAATGTACCTTGGCCTTCATGCAGGGCGAAATTGCAACTATCTTAGCTGGGGATGTTAAAGTGAAA
 AAAGAGAGAGACCT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA


Protein Sequence: >MG210579 representing BC063253
 Red=Cloning site Green=Tags(s)

```
MSAAQVSSRRQSCYLCDLPRMPWAMIWDFSEPVCRCGVNYEGADRIEFVIETARQLKRAHGCFQDGRSP
GPPPPVGVKTVALSAKEAAAAAAAAAQQQQQQQQQQQLNHVDGSTKPAVLAAPSGLERYGLSAAAAAA
AAAAVEQSRFEYPPPPVSLGSSSHAARLPNGLGGPNGFPKPAPEEGPELNRQSPNSSAATSVASRRG
THSGLVTGLPNPGGGGGPQLTVPPNLLPQTLNNGPASA AVLPPPHGLGSGRGPPTPAPPGAPGPA CLGG
PPGV SATVSSAPSSTSTVAEVGVGAAGKRPGSVSSTDQERELKEKQRNAEALAE LSESLRNRAEEWANK
PKMVRD TLLTLAGCTPYEVRFKKDHSLLGRVFAFDVAVSKPGMDYELKLFIEYPTGSGNVYSSASGVAKQM
YQDCMKDFGRGLSSGFKYLEYEKKHSGDWRLLDLLPEAVRFFKEGVP GADMLPQPYLDASCPMLPTAL
VLSRAPSAPPGTGALPPAAPTGRGAASSLRKRKASPEPPDSAESALKL GEEQQRQWMANQSEALKLTM
SAGGFAAPGHSAGGPPPPPPPLGPHSNRTTPESAPQNGPSPMAALMSVADT LGTAHSPKDGSSVHSTTA
SARRNSSPVSPASVPGQRRLASRNGDLNLQVAPPPSAHPGMDQVHPQNI PDSPMANSGLCCTICHER
LEDTHFVQCPSVPSHKFCFPCSRESIKAQGATGEVYCPSGEKCLVGSNVPWAFMQGEIATILAGDVKVK
KERDP
```

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: BC063253


ORF Size: 2325 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

- Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
- Note:** Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
- RefSeq:** [BC063253](#), [AAH63253](#)
- RefSeq Size:** 3205 bp
- RefSeq ORF:** 2327 bp
- Locus ID:** 238330
- Cytogenetics:** 12 D2
- Gene Summary:** Probable E3 ubiquitin protein ligase involved in the proteasome-mediated ubiquitin-dependent degradation of target proteins. Through the degradation of CTNNB1, functions downstream of FOXF2 to negatively regulate the Wnt signaling pathway. Probably plays a role in the development of the central nervous system and in neuronal maintenance (By similarity). Also acts as a transcriptional regulator of genes controlling female reproductive function. May play a role in gene transcription by transactivating GNRH1 promoter and repressing PENK promoter (By similarity).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MG210579