

Product datasheet for **MG210595**

Kirrel3 (BC063072) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Kirrel3 (BC063072) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Kirrel3
Synonyms:	SST4, mKirre, NEPH2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:




ACCN:	BC063072
ORF Size:	2333 bp



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	BC063072 , AAH63072
RefSeq Size:	3585 bp
RefSeq ORF:	2333 bp
Locus ID:	67703
Cytogenetics:	9 A4
Gene Summary:	Synaptic adhesion molecule required for the formation of target-specific synapses (PubMed:23637329, PubMed:26575286). Required for formation of target-specific synapses at hippocampal mossy fiber synapses. Required for formation of mossy fiber filopodia, the synaptic structures connecting dentate granule and GABA neurons. Probably acts as a homophilic adhesion molecule that promotes trans-cellular interactions and stabilize mossy fiber filopodia contact and subsequent synapse formation (PubMed:26575286). Required for the coalescence of vomeronasal sensory neuron axons (PubMed:23637329). May be involved in the hematopoietic supportive capacity of stroma cells; the secreted extracellular domain is directly responsible for supporting hematopoietic stem cells (PubMed:12665856). [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MG210595