

Product datasheet for **MG210686**

Usp10 (NM_009462) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Usp10 (NM_009462) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Usp10
Synonyms:	2610014N07Rik; mKIAA0190; UBPO; Uchrp
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MG210686 representing NM_009462
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGCCCTCCACAACCCACAGTATATCTTTGGCGATTTCAGCCCTGATGAATTCAATCAGTTTTTTGTGA
 CTCCCGCTCTTCTGTTGAGCTCCCTCCATACAGTGGGACTCTGTGTAGCATACAGGCTGAAGATGAAC
 GCCAGATGGACAGGAGCACCAGAGGATCGAGTTTGGTGTAGATGAAGTCATCGAACCTAGTGAGGGGCTG
 CCGCAACCCCTAGCTACAGTATTTCAAGCACCTTGAACCCCGAGGACCCGAATTTATCCTTGGTTGTA
 CGACTTCTAAAAAGATCCCTGATGCCGTTGAAAAAGATGAGACCTACAGCTCCATTGACCAGTACCCCGC
 CTCGGCCTTGGCTCTGAAAAGCAACTCTAACGCAGAGGCTGAGACCCTGGAGAACGACAGTGGTGCCGGC
 GGCCTTGGTCAGAGGGAACGAAAAAGAAGAAGCGGCCCTGGGACTACAGTTATTTGAAAGATG
 GCGGTGAGGACAGTGCTTCCCAGCCACCCTGGTCAACGGCCATGCCACCTCAGTGGGCACAAGCGCGCA
 GGCTGTGGAGGACGCCGAGTTTCATGGACGTGCTTCTCCAGTCATGCCAGGACTTGTGACAGCCCTCAG
 AATCCTGTGGACTTCATCAGTGGCCCTGTCCCTGACAGTCCGTTCCCAGGACACTAGGAGGCGATGCCA
 GGACTGCAGGGCTGTGTGAGGGCTGCCATGAAGCTGACTTTGAGCAGCCCTGCCTCCCTGCAGACAGCCT
 GCTAAGGACAGCTGGGACGCAGCCCTACGTTGGCACCGATACTACTGAGAACTTCGCAGTTGCTAATGGA
 AAAACTTTGAATCCCGGGCGAGGACACTGCTGCCAGTGGGGCAGAAGTGCACACCGACGAGGGTGACG
 ACCTGGACCCTGCGAAGCCGGAGAGCCAGTCACCTCCTGCCGAAAGTGCCTCTCTGCCTCTGGCGCCAT
 CCCCATTAGCCAGCCTGCAAAGTCTGGGCTAGTCTCTTTTCATGATTCTAAGCCCTTCTCCTCCTACCC
 ATGGCGTATGTGAAACTAAGTGTCCCCCCTGTCCCGTCCCCCTGGCCTCTGAAAAACAGATGGAAG
 TCAAAGAAGGGCTTGTCCAGTGTGAGAGGATCCTGTAGCCATAAAGATTGCAGAGTTATTGGAGACTGT
 AACCTAATACATAAGCCAGTGTATTGCAACCCCGTGGGCTGATCAATAAAGGAAACTGGTGTACATT
 AATGCCACCCTGCAGGCAATTGGTGGCTTCCCTCCGATGTATCACCTGATGAAGTTCATTCTCTGTACT
 CAAAAGTGCAGCGCCTTGCACATCAACGCCCATGATAGATAGCTTTGTTGCGCTCATGAATGAGTTCAC
 TAACATGCCAGTGCCTCCAAACCCCGCAAGCTCTCGGGGATAAAATCGTGAGAGATATCCGCCAGGG
 GCTGCCTTTGAACCCACATACATATACAGACTCCTGACAGTCATCAAGTCGAGCCTGTCTGAAAAGGGCC
 GGCAGGAGGACGCAGAGGAGTATCTAGGTTTCATCCTCAATGGACTCCATGAAGAGATGCTGAGCCTGAA
 AAAGTACTCTCACCCACGCACGAAAAACTCTGTTTCCAATGGGCCAGAAAGCGACTTGATAGAGGAC
 GAGGAGCTGGAAGACACGGCAAGGGCAGCGAGGACGAGTGGGAGCAAGTGGTCCCAAGAATAAGACGT
 CCATACCCGGCAGGCGGACTTCGTCCAGACGCCCATCACTGGCATTTCGGGGACACATCAGGTCTGT
 GGTCTACCAGCAGAGTTCCAAGGAGTCAGCCACTCTGCAACTGTTTTTCACGCTGCAAGTTGGACATCCAG
 TCTGACAAGATACGCACAGTCCAGGATGCGCTGGAGAGCTTGGTGGCCAGGGAGTCTGTCCAAGGTTACA
 CCACAAAACCAAGCAGGAGGTTGAGGTGAGCCGAGAGTACTCTGGAGAAGCTGCCCCCTGTCCTCGT
 GCTCCACCTGAAGCGCTTTGTCTACGAGAAGACAGGCGGATGCCAGAAGCTGGTCAAGAACATTGACTAC
 CCTGTGGACTTGGAGATCAGCAGAGAATACTTTCTCCAGGCATCAAAAATAAAATTTTAAATGCCAGA
 GAACCTATAGGCTGTTTGCAGTGGTCTACCACCATGGCAACAGCGCTACAGGCGGCCACTACACTACGGA
 CGTCTTCCAGATTGGGCTTAATGGCTGGCTACGAATCGATGACCAAACGGTCAAGGTTATTAACCATAC
 CAGGTGGTGAACCGCCTGCCGACCGCACAGCCTACCTCCTATATTACCGCCGCTGGACCTGCTG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG210686 representing NM_009462
 Red=Cloning site Green=Tags(s)

```
MALHNPQYIFGDFSPDEFNQFFVTPRSSVELPPYSGLTCSIQAEDELPGQEHQRIIEFGVDEVIEPSEGL
PRTPSYSISSTLNPQAEFILGCTTSKKIPDAVEKDETYSSIDQYPASALALESNSNAEAETLENDGAG
GLGQREKKKKKRPPGYYSYLKDGGEDSASPATLVNQHATSVGTSGEAVEDAEFMDVLPVMPRTCDSPQ
NPVDFISGPVPDSPFPRTLGGDARTAGLCEGCHEADFEQPCLPADSLLRTAGTQPYPVGTDTTENFAVANG
KILESPGEDTAASGAELHTDEGADLDPKAPESQSPPAESALSASGAIPISQPAKSWASLFHDSKPSSSSP
MAYVETKCSPPVPSPLASEKQMEVKEGLVPVSEDPVAIKIAELLETVTLIHKPVSLQPRGLINKGNWCYI
NATLQALVACPPMYHLMKFIPLYSKVQRPCTSTPMIDSFVRLMNEFTNMPVPPKPRQALGDKIVRDIRPG
AAFPTYIYRLLTVIKSSLSEKGRQEDAEYLGFI LNGLHEEMLSLKLLSPTHEKHSVSNRPSDLIED
EEL EDTGKSEDEWEVQVGPKNKTSITRQADFVQTPITGIFGGHIRSVVYQQSSKESATLQLFFTLQLDIQ
SDKIRTVQDALESLVARESVQGYTTKTKQEVESRRVTLEKLPVVLVHLKRFVYEKTGGCQKLVKNIDY
PVDLEISRELLSPGIKKNFKCQRTYRLFVAVYHHGNSATGGHYTTDFVQIGLNGWLRIDDQTVKVINQY
QVVKPPADRTAYLLLYRRVDLL
```

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_009462


ORF Size: 2379 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_009462.1</u> , <u>NP_033488.1</u>
RefSeq Size:	3236 bp
RefSeq ORF:	2382 bp
Locus ID:	22224
UniProt ID:	<u>P52479</u>
Cytogenetics:	8 E1
Gene Summary:	<p>Hydrolase that can remove conjugated ubiquitin from target proteins such as p53/TP53, BECN1, SNX3 and CFTR. Acts as an essential regulator of p53/TP53 stability: in unstressed cells, specifically deubiquitinates p53/TP53 in the cytoplasm, leading to counteract MDM2 action and stabilize p53/TP53. Following DNA damage, translocates to the nucleus and deubiquitinates p53/TP53, leading to regulate the p53/TP53-dependent DNA damage response. Component of a regulatory loop that controls autophagy and p53/TP53 levels: mediates deubiquitination of BECN1, a key regulator of autophagy, leading to stabilize the PIK3C3/VPS34-containing complexes. In turn, PIK3C3/VPS34-containing complexes regulate USP10 stability, suggesting the existence of a regulatory system by which PIK3C3/VPS34-containing complexes regulate p53/TP53 protein levels via USP10 and USP13. Does not deubiquitinate MDM2. Deubiquitinates CFTR in early endosomes, enhancing its endocytic recycling. Involved in a TANK-dependent negative feedback response to attenuate NF-kappaB activation via deubiquitinating IKBKG or TRAF6 in response to interleukin-1-beta (IL1B) stimulation or upon DNA damage. Deubiquitinates TBX21 leading to its stabilization. [UniProtKB/Swiss-Prot Function]</p>

Product images:



Circular map for MG210686