

## Product datasheet for **MG210869**

### **Bahcc1 (BC060615) Mouse Tagged ORF Clone**

#### **Product data:**

|                           |                                          |
|---------------------------|------------------------------------------|
| Product Type:             | Expression Plasmids                      |
| Product Name:             | Bahcc1 (BC060615) Mouse Tagged ORF Clone |
| Tag:                      | TurboGFP                                 |
| Symbol:                   | Bahcc1                                   |
| Synonyms:                 | MGC79150, B930044J06                     |
| Mammalian Cell Selection: | Neomycin                                 |
| Vector:                   | pCMV6-AC-GFP (PS100010)                  |
| E. coli Selection:        | Ampicillin (100 ug/mL)                   |



[View online »](#)

## ORF Nucleotide Sequence:

>MG210869 representing BC060615  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCCCGCATCGCC

ATGTTGCAGCCAGTCCTGAGGCGGAAGAATGGGCTTGTCCATCGCTCTGTGAGCCGCAACGCCAAGG  
 CCATCCTGGGAAAGAGCAGGAAGCTGACGAAGGTGAAGAGAGAGGCTGTGAGCAAGCAGGGCCAGGGGCG  
 GGCGGTGAGCCGGCTGCTGGAGAGCTTTGCTGTGGAGGATGACTTTGAGTTTGACGAGGATGACACCAGC  
 TTCTCGGATGAGGAAGAAGAAGAGGAAGAAGCCGGTGTACAGCTCAGCGCAGAGCAGAGCGCCGCCCTGG  
 CACGGTCTGCACCATCCACAAGGAAGACCTGCAGGACGAGTGGCCGCTCTGATCCCGAAGGAGGACAG  
 CCTGCTGTATGCGGGCAGTGTGAGGACCTGCAACCCCCGACATCTACAGCATCGTCATTGAGGGGGAG  
 AGAGGCAACCGTCAACGTATCTACTCGTGGAACAGCTGCTGAGGAGGCGGTTCTTGATGTCCAGCCAC  
 AGTCCAGTCGGTACCTTCCACCTGGTACCCGGGTCTGTGCTACTGGAGTCAGAAGTCCCGGTGTCTATA  
 TCCAGGCAATGTGTTTCGAGGTGCTTCAAGCGATGAAGAAGACCTGGACTCAGTGTGGTAGAGTTCGAT  
 GACGGAGACACAGGCCACATAGCTGTCTTAACATCAGGCTGCTGCCGCCGACTTCAAGATCCAGTGTA  
 CAGAGCCCTCGCCAGCCCTGTTGGTGTCCAGCAGCTGCCGGAGGACCAAAAAAGCAGCCAATGAGGGCCA  
 CCCGCCAGCGAAGCCCCACACCCAGCCTGTCCCTAAAGTGCCGGATGGCCCTGAAACTTCTAAGACC  
 CCTGGGAAGAAATCTGGCAGCAAGACAAGCTGGCAAGTGCACCTCCTAACCTCAGGTGCCAAGTCCC  
 CCACGGGGGCTCAGACCACTTCTAGGCCGACAGGACGCCCCCTGTGAGCTGGTCAGCAGTGGCTCA  
 GACCAAGCGGAAAGCGGTGGCTGCGGCAGCGCGGCGGAGGTGGCAAGGGCCGGGGTACTACAGAAC  
 CTCTTCCAGCTCAACGGAAGCACCAAGAAGCTGCGGGCCCGGACACCCTGTTTCCCATGCATAGCATGG  
 CTACCCCGGTGTTTGGTAACAGTTCGAGCTGACTCCTTACAGCAGCTGGCGAGTTCCTACACACCCTT  
 CCTTGGAGGGGCTGGGGCAGGCCTTCCCGAGGAGGCCCAAGCTGCTTAGGGCCAAGAAGGCTGAACGG  
 GCTGAAGCAGAAAAGGCCGGGAGGCGGCGGCGAGGCGGCGAGTTTCTGGTTAAGCTGGACCATGAGGGTG  
 TGACCTCCCCAAAAACAAAACTGCAAGGCCCTGCTCATGAGCGACAAGGACTTCGACCCAAGCTGGG  
 GCGACCTCTTAGCAACCCAGTTACGCACACCCAGCCCTCATTGGCAAGGACAAGAAGGGTCGGGCGCCT  
 GTGCATCCGCTACCATGGGGCTGGCTCTGCGCAAGTACCCACTGCCCTGTGATAGTGACTGCCCCAGCT  
 CCTACTCCGATGAGGACGAGGATGGGCCGGGGCTGGCCACCGGTGTTCCCTCCCGATTCTCACCCGCCT  
 ATCCATGTCGTCCTCTCTGGCTCGTCCACGTCTTCTTCTCAGGCTCCGTGTCCACGTCCAGCCTC  
 TGTTCTCAGACAATGAGGACTCATCTTATAGCTCAGATGACGAGGACCCTGCCCTGCTGCTGCAGACCT  
 GTCTACCCGCCCTGTGCCCCCCTCTGGCCCCACCCGAAGCCCTGCGCTCTAAGGGTAGCAGCCCCCA  
 CGCGCACACCCACGCCAGCGCTGCTTCTGTCCAGGGCTGGGGTGGCTGGTGCAGGTGCAGGGGCCAGC  
 CCCAGTGGTAGCAATCCAAGTTCAAGCGCAAGGAGGCCCTGAGCTTCTCCAAAGCCAAAGAGCTTTCTC  
 GGAGGCAACGGCTGCCCTCCGTAGAAAACCGGCCAAAGATCTCAGCCTTCTGCCCTGCCCGGCAGCTCTG  
 GAAGTGGTCAGGAAACCCACGCAGAGACGAGGCATGAAGGGGAAAGCCAGGAACTATTCTACAAGGCC  
 ATTGTCCGAGGCAAGGAGACACTGCGCATCGGGGACTGTGCGGTCTTCTTTGCGCCGGCCGGCCCAACC  
 TGCCCTACATAGCCGCATCGAGAGCTTGTGGGAGTCTGGGGCAGCAACATGGTGGTGAAGGTCAAATG  
 GTTCTACACCCAGAGGAGACCAAGCTGGGGAAGCGGCAGAGTGACGGGAAGAACGCATGTACCAAGTCC  
 TGCCACGAAGATGAGAATGACGTGACAGACCATCTCCACAAATGTCAGGTGGTGGGCCGAGAGCAGTATG  
 AACAGATGATGCGGGGCCGCAAGTACCAGGACAGCAGGACCTGTACTACTTGGCAGGCACCTATGATCC  
 TACTACTGGGCGCCTGGTGACAGCTGACGGCGTCCCCGTCTGTGC

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

## Protein Sequence:

>MG210869 representing BC060615  
 Red=Cloning site Green=Tags(s)


```
MLQPVLRRKNGALSIALSARNAKAILGKSRKLTkVvKREAVSKQGQGRAVSRLLESFAVEDDFEFDEDDTS
FSDEEEEEEEAGVQLSAEQSAALARSCTIHKEDLQDGLPVLIPKEDSLLYAGSVRTLQPPDIYSIVIEGE
RGNRQRIYSLEQLLQEAVLDVQPQSSRYLPPGTRVCAYWSQKSRCLYPGNVVRGASSDEEDLDSVLVEFD
DGDGTGHIAVSNIRLLPPDFKIQCTEPSPALLVSSSCRRTKKAANEGHPPSEAPTPSLSPKVDPGPETSKT
PGKKSQKDKAGKVDLLTSGAKSPTGASDHFLGRRGSPLLWSAVAQTKRKAVAAAAAAGGKGPVQLQN
LFQLNGSTKKLRARDTLFPMHSMATPVFGNSFRADSFSSSLASSYTPFLGGAGAGLPGGAHKLLRAKKAER
AEAEGAGRRRAGGEFLVKLDHEGVTSPKNKNCKALLMSDKDFGPKLGRPLSNPSYAHPALIGKDKKGRAP
VHPLPMGLALRKYPLPCDSDCPSSYSDEDEDGPGLATGVPSRFLTRLMSSSSSSGSSTSSSSGSVSTSSL
CSSDNEDSSYSDDDEDPALLLQTCLTRPVALLAPPEALRSKGSSPHAHTHAQRCFLSRAGVAGAGAGAS
PSGSKSKFKRKEALSFSKAKELSRQRLPSVENRPKISAFLLPARQLWKWSGNPTQRRGMKGKARKLFYKA
IVRGKETLRIGDCAVFLSAGRPNLPIYIGRIESLWESWGSNMVVKVWVFYHPEETKLGKRQSDGKNALYQS
CHEDENDVQTI SHKCQVVGREQYEQMMRGRKYQDQQLYYLAGTYDPTTGRLVTADGVPVLC
```

TRTRPLE – GFP Tag – V

## Restriction Sites:

SgfI-MluI

Cloning Scheme:



|                 |                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|-----------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| ACCN:           | BC060615                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| ORF Size:       | 2496 bp                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a> |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.                                                                                                                                                                                                                                                                                                               |
| Components:     | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).                                                                                                                                                                                                                                                               |

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**Note:** Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

**RefSeq:** [BC060615](#), [AAH60615](#)


**RefSeq Size:** 5004 bp

**RefSeq ORF:** 2498 bp

**Locus ID:** 268515

**Cytogenetics:** 11 E2

## Product images:



Circular map for MG210869