

## Product datasheet for **MG211401**

### Ankrd27 (BC065093) Mouse Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Ankrd27 (BC065093) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Ankrd27
Synonyms:	MGC25907, Varp
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide Sequence:**

>MG211401 representing BC065093  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**


ATGGCTCTATATGATGAAGATCTCTGAAAAACCTTTCTACCTGGCCCTTCAGAAGTGGCGCCCGACT  
 TGTGCAGCAAGGTGGCCAGATCCATGGCATTGTCTAGTGCCTTGCAGAGGAAGCTGCCGGCAGTGT  
 GCAGGCCCTCTGCCAGTTCGAGTCTATGTTTTGGTACCCACGGAAGGACACTTTCAGACCTTAGATGGA  
 AAGGCTGTCGTCATTGAAGGAAACAGGATTAAGCTAGGAGCAGGTTTTGCTTGCCTTCTCTGTGCCCA  
 TCCTCTTTGAGGAGACTTCTACAACGAGAAAGAGGAGAGTTTCAGCATTCTCTGCATTGCTCATCTTT  
 GGAGAGGAGAGAGACTTCAGAAGAACCTTCAGCGCCTGCAGATCTTTCTCCCTGAAAACCATCGAAGAT  
 GTGAGAGAATTTTTGGGAAGACTCAGAGAAATTTGACAAGAACATTGCCTCTTTCCACCGGACGTTCC  
 GAGAATGTGAAAGAAAGAGCCTCCGCCACCACATAGACTCCGTAATGCTCTCTACACCAATGCCTCCA  
 GCAGTCTCTCGAGACTCTCACCTGAAGTTCTTGCAAAGCAGGAGGCCAGATGAACCTGATGAAGCAG  
 GCCGTAGAGATGTATGTCCATCACGATATTTACGACCTGATCTTTAAATACGTGGGGACCATGGAGCGCA  
 GCGAGGATGCCGCCTTTAACAAAATCACGAGAAGCCTTCAAGATCTCCAGCAGAAAAGACATCGGCCTGAA  
 ACCCGAGTTACGCTTCAACATCCCTCGGGCAAAGAGAGAGCTGGGTGAGTTGAACAAGTGTACGTGCCA  
 CAGCAGAAGTGTGTGCTGAGGAAGTGGTCCAGCTCATGACACAATCTCCAGCCAGAGAGTGAACCT  
 TGGAGACCATGTGTGCCGATGATCTCTCTCGTCTGTTATATCTGCTCGTGAAGACAGAGATCCCTAA  
 CTGGATGGCAAATTTGAGCTACATCAAAAACCTCAGATTTAGCAGCTCGGCCAAAGATGAGCTAGGATAC  
 TGCTGACCTCAGTCGAGGCTCGATTGAATACATCCGGCAAGGAAGTCTCTCCACGAAGACGCCTGACG  
 CTGAGGGATTTGGCGACAGACTATTTCTCAAGCAGCGGATGAACCTGCTGTCCAGATGACATCAACACC  
 CATCGACTGCCTGTTCAAGCACATCGCCTCTGGTAACCAAAAAGGAAGTGGAAAGACTTCTGAGCCAAGAT  
 GACCAGGACAAAAGACCCATGCAAAAAGATGTGCCACCCACTTTGTTCTGTGAGGACTGTGAGAAGCTTA  
 TCTCTGGGAGGCTGAATGACCCGTGAGTTGCTACTCCTTTCTCCAGAGATGACAGAGGTCAAACACCCCT  
 CCATGTGGCTGCCCTCTGTGGTCAGGCATCCCTCATTGACTTCTGGTCTCCAAAGGTGCCGTAGTGAAC  
 GCCACAGACTATCATGGTCCACCCCTTACATCTGGCATGTGAGAAAGGCTTCCAAAGTGTGACGCTGC  
 TGCTGCTGCACTACAAGGCCAGCACTGAGGTTCCAGGACAACAACGCCAACACCCCGCTGCACCTGGCCTG  
 CACCTATGGGCAGGAGGATTGCGTGAAGGCCCTGGTCTACTATGATGTCCAGGCTGCAGATTGGATATT  
 GGTAATGAGAAAGGAGACACAGCCCTGCACATTGCTGCACGCTGGGGTTACGAGGGCATCATAGAGACAC  
 TGCTGCAGAATGGAGCCCCACAGCAGTCCAGAACAGACTGAAGGAAACACCACTCAAGTGTGCTCTGAA  
 CTCGAAGATTCTGTCTATCATGGAAGCTCATCACCTATCCTCTGACAGGAGGCCAAGGCCTTCTGAGGTT  
 CCTGCACAGTCCCTACACGCTCTGTGGATTCCATCAGCCAGGGCTCCTCCACCTCCAGCTTTTCTTCCA  
 TATCGGTGAGCTTACAGCAAGAAGAGGTCAAAAAGGACTACAGGGAGAAACTGGCCAGGATTTCTGCCAA  
 TGGGCTTAGTGTGAATGTGACCAACCAGGACGGCTTCTCTCCCTGCACATGGCTGCCCTGCATGGCCGC  
 ACAGACCTCGTCCCCCTCTGCTGAAGCACGGTGTACTCAGGTGCCAGAAACAAGCCAAGCCGTCC  
 CGCTCCACCTGGCCTGCCAGCAGGGCACTTTCAGGTGGCAAAATGTCTTTAGATTCCAATGCAAAACC  
 TAATAAAAAGGACCTCAGTGGGAACACACCCCTCATTGTGCTGCTGCTGCCGGCCACCATGAAGTCGCT  
 GCGCTGTTGCTACAGCACGGGGCTCCATCAATGCTTGCAACAACAAGGCAACACGGCCCTGCATGAGG  
 CTGTGATGGGAAGGCACACCCTGGTAGTGGAGCTGCTTCTGTTCTACGGGGCGTCGGTTGACATCTGAA  
 CAAGAGGCAGTACACGGCGGCGACTGTGCCGAGCAGGACTCAAAAATAATGGAGTTGCTGCAGGTAGTG  
 CCAGGCTGTGTTGCGTCTCTAGACAGGCTCGAGGAGGCAGACCACGAGGGTTACGTGACTGTTGAGATCA  
 GGAGAAAATGGAACCCCAAAATGTACAACCTGCCAGAGGAACCTCTTAGGAGACAGTTTTGCCTCATTAA  
 CCCGGGGGTGCGTTCCAGGAGAGGACATCCAGGGAGACTATGGGAAGGGATAGAAGTGTGCCGATCTT  
 GCAGGACGTTCTTTCAGGAGCCAGAGAAGCAAAGAGTCAACGGAAAGCAGAGTGACCTGTCAGACCTGA  
 GCAGGTATCAGACCTCTGAGGAAGGAAACAAGGGGCTCCCGGAGAGGCCCGTGTGACAGACAGGCTGTCC  
 AGGACACCGCCGATGGTCCGCAGACACACAGTCAATGACGCAGCCATACTCCAGTCCCAGAGGTGACT  
 GTCACCTCACCCACAGGGCAGTGTCTCAGTCT

**ACGCGTACGCGGCCGCTCGAG** - GFP Tag - GTTTAA



<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">BC065093</a> , <a href="#">AAH65093</a>
<b>RefSeq Size:</b>	4105 bp
<b>RefSeq ORF:</b>	2981 bp
<b>Locus ID:</b>	245886
<b>Cytogenetics:</b>	7 B2
<b>Gene Summary:</b>	May be a guanine exchange factor (GEF) for Rab21, Rab32 and Rab38 and regulate endosome dynamics (By similarity). May regulate the participation of VAMP7 in membrane fusion events; in vitro inhibits VAMP7-mediated SNARE complex formation by trapping VAMP7 in a closed, fusogenically inactive conformation (By similarity). Involved in peripheral melanosomal distribution of TYRP1 in melanocytes; the function, which probably is implicating vesicle-trafficking, includes cooperation with Rab32, Rab38 and VAMP7 (PubMed:19403694, PubMed:21187289). Involved in the regulation of neurite growth; the function seems to require its GEF activity, probably towards Rab21, and VAMP7 but not Rab32/38 (PubMed:19745841, PubMed:22171327). Proposed to be involved in Golgi sorting of VAMP7 and transport of VAMP7 vesicles to the cell surface; the function seems to implicate kinesin heavy chain isoform 5 proteins, GOLGA4, RAB21 and MACF1. Required for the colocalization of VAMP7 and Rab21, probably on TGN sites (By similarity). Involved in GLUT1 endosome-to-plasma membrane trafficking; the function is dependent of association with VPS29 (By similarity). Regulates the proper trafficking of melanogenic enzymes TYR, TYRP1 and DCT/TYRP2 to melanosomes in melanocytes (PubMed:26620560).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MG211401