

Product datasheet for **MG212190**

Klk9 (NM_028660) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Klk9 (NM_028660) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Klk9
Synonyms: 1200016C12Rik; AI324041
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG212190 representing NM_028660
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAAGCTGGGACTCACCTTGGTGCTGTTCTCTCTGCTGGCAGGCCACTGTGGAGCTGACACCCGTGCCG
TCGGAGCCCGGAATGCGTACGCAACTCTCAGCCCTGGCAGGCTGGCCTTTCTACCTCACCCGGCAGCT
CTGTGGAGCCACGCTCATCAATGATCAATGGCTGCTCACAGCTGCTCACTGCCGAAAGCCGTACCTGTGG
GTCCGCCTCGGGAGCATCATCTGTGGCGGTGGGAGGCCCTGAGCAGCTGCTTCTGGTCACAGACTTCT
TCCCCATCCTGGTTCAATCCAGATCTCTGCCAATGACCACAATGATGATATTATGCTGATCCGGCT
ACCTCGGAAGGTCCGGCTGACTCCAGCTGTGCAACCCCTCAACCTCACCGAGAGCCGGCCACCTGTGGGC
ACCCAGTGCCTTATCTCTGGCTGGGTCAGTGTGTCAAGCTCCAAATTACAGTATCCAATGACACTGCAGT
GTGCCAACATTAGCATTCTGGACAACAAGCTCTGCCGCTGGGCATACCCGGGCCACATCTCTGAGAAGAT
GCTGTGCGCTGGCCTGTGGGAGGTTGGCCGAGGCTCTTGCCAGGGTGACTCTGGGGGACCCCTGGTGTGC
GAAGGAACCTTGGCAGGTATAGTGTCTGGAGGATCTGAGCCTTGTCCAGACCCCGGAGACCCGCTGTCT
ACACTAACGTCTTTGACTACCTGGAGTGGATTGAAAGCACCATGGAAAAAACC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG212190 representing NM_028660
Red=Cloning site Green=Tags(s)


MKLGTLVLVFLSLAGHCGADTRAVGARECVRNSQPWQAGLFYLTRQLCGATLINDQWLLTAAHCRKPYLW
 VRLGEHHLWRWEGPEQLLLVTDFFPHPGFNPDLSANDHNDIMLIRLPRKVRLTPAVQPLNLTESRPPVG
 TQCLISGWSVSSSKLQYPMTLQCANISILDNKLCRWAYPGHISEKMLCAGLWEGGRGSCQGDSSGGPLVC
 EGTLAGIVSGGSEPCSRPRRPAVYTNVFDYLEWIESTMEKN

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:



ACCN: NM_028660

ORF Size: 753 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_028660.3 , NP_082936.2
RefSeq Size:	1475 bp
RefSeq ORF:	756 bp
Locus ID:	101533
Cytogenetics:	7 28.26 cM