

## Product datasheet for **MG212252**

### C130026I21Rik (NM\_001037909) Mouse Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	C130026I21Rik (NM_001037909) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	C130026I21Rik
Synonyms:	4930565N07Rik; OTTMUSG00000029174
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG212252 representing NM_001037909 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCAGGAGGCAACAATGAACTCAGCAGCAGGACAATCCAGAGGACCAGAATGAAGAGGAGTCAGATG  
ACTATCAGCTCATGTTCAAGCACTTCAAAGAAAACAAGGTGGAGATCGCAAGTCCATCACCAAGCCGTT  
TCCTTTCCTCATGAGCCTTCGAGACCGAGACTTCATCTCGGAGCAGAAGTTTCAGGAATATCAAGAAACA  
TGTAAGAACCTGGTCCCAGTGGAGAGAGTCGTGTACGACGTCCTCAGCAATGTACAGAAGAAGTTTAGCC  
GGGATCTTCTGAAAGTGATATTCAGCAAAACACATCTGAAGGCCTACCCGACTTAAAGGAGACTCTAAA  
ACACTTCTTCTAAATGCATCCAAAACGAATGATGAACAAGCAGAGGAGATGCTTAGCCTCCACAAATGC  
AATGGAGGAGAGGGAAGTAGTTCTTGTGAACAAACGTGTGATGAACAAGAGCCCCAGGATGACCTGCCCT  
CATCCCTGAGACAAGAAGCAGCAAGAAGAAAAAACGTCCGAACTGGTCAAATTCAAAAGAAGACTGCA  
GAAAAAAGCCCGTCAAGATGAGATGATGGGAGTGCCCTCACCTGGACATGGAGTTCAAGAGAAGCTC  
AAGGCGGTGAGCAGGAGGACTTTGTGAAAAGATGACTCATCTACGAACGTGAAGGAGGTGACCAAGACAC  
TGAGAGCAAGGATGAGGTGTGGCACCTGCCATTGGAAAGGCTCAGCCTGGGCAAATCTCTCAGAGCTGGT  
GCAGGGTTCGAGGGATACA

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >MG212252 representing NM\_001037909  
 Red=Cloning site Green=Tags(s)

MAGGNNELSSRTIPEDQNEEESDDYQLMFKHFKENKVEIASAITKPFPLMSLRDRDFISEQKFQEYQET  
 CKNLVPVERVVYDVL SNVQKKFSRDLLKVI FSKTHLKAYPDLKETLKHF FLNASKTNDEQAEMLSLPQC  
 NGGEGSSSCEQT CDEQEPQDDL PSSLRQEAARRKKRPNWSNSKRRLQKKKPRQDEMMPVSPGHGVQEKL  
 KAVSRRTLWKDDSSSTNVKEVTKTLRARMRCGTCHWKGSWANLSELVQGSRDT

TRTRPLE - GFP Tag - V


**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**

Cloning sites used for ORF Shutting:



**Plasmid Map:**



**ACCN:** NM\_001037909

**ORF Size:** 789 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001037909.3</a> , <a href="#">NP_001032998.2</a>
<b>RefSeq Size:</b>	1887 bp
<b>RefSeq ORF:</b>	792 bp
<b>Locus ID:</b>	620078
<b>Cytogenetics:</b>	1 C5