

## Product datasheet for **MG212404**

### **Pibf1 (NM\_029320) Mouse Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	Pibf1 (NM_029320) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Pibf1
Synonyms:	1700017E21Rik; 4930513H15Rik; 4933438D16Rik; 4933439E17Rik; D14ErtD581e
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide Sequence:**

>MG212404 representing NM\_029320  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGTCCTCGCAAAATTGCCAAGGAACCAAAAAAGTAAATATCTCTAGTTCTCTGGAGTCTGAAGATTA  
 GTTTGAAACAACCATTCATACAGATGATGTCTCATCATCAGAGGAGCGAGAAGGTAAGTCAAAATCAC  
 CAGGCAGTTAATCGAAAGAAAAGAGATACTTCATAATATTCAGTTACTGAAAATCGAGCTATCCCAGAAA  
 AACATGATGATCGACAACCTTGAAGATGGATTATCTTACAAAGATTGAGGAGCTAGAGGAAAACTTAATG  
 ACGCCCTTACCAGAAGCAGCTGCTAACTTTGCGATTAGACAATCAGTTGACTATTCAACAGAAAGATGC  
 CAAAAATATCAAGAATAATGAAACAAGAAATGGAAACCATTTATTGCGACAGAAGCAACTGGAAGAA  
 ACAAAACCATCAGCTGAGAGAAAAGGCTGGAGATGTTCCGCCGAAATCTGCGAGACTTTGAGCTGACAGAAG  
 AGCAGTATGTGAAGCTAAAATCTTTTCTGAAGATCAACTCTCTATTCTGAATATGTATCTATTCGCTT  
 CTATGAGCTCGTGAACCCATTAAGAAAGGAAGTCTGTGAGCTACAGGTGAAGAAGAGTGAAGTCTCTGAA  
 GAAGTGAAGTACAAGTAAAGGCCAACTGAAGCAGCTGACGGAGACATATGAAGAAGATCGAAGAAAACACG  
 CTGAACCTCTAATTCGATGTCAACGTTTACCTTGAATTAGCAGACACAAAACAGTTAGTTTCAGCAAGG  
 TGATTACCGTCAAGAGAACTATGACAAAGTGAAGAGTGAACGCGATGCTCTGGAACAGGAGCTACTCGAG  
 CTTAGAAGAAAACACGAAGTACTTGAAGCCTCTCACATAGCTCAAGCTAAGGAAAGGAATGAATTATCAA  
 AGGAGGTCAGCAGCCTGCAGCAGACAGTCAACCTGCTGCAGAAGGATAAAGACTACCTCAATCGCCAAAA  
 CATGGAAGTCAAGTACGCTGTGCCATGAGGAGGATCGGCTGGAAGGCTGCAAGTTCAACTGGAAGAC  
 ACCAAAAAGGCTAGAGAAGAGATGTATGAGAAATATGTACGTCAGAGACCATTATAAAACAGAATATG  
 AAAATAAATCATGATGAAGTGAACAAATCAAATGAAAATCTAGAAATGATCAGCTTTCGAAG  
 TGCTCTAGGGAAATGTATGAACGAGAAAACAGAAATCTCCGTGAAGCAAGGGATAATGCACTCGCTGAA  
 AAGAACCAGCAGTGGCAGCGAAAAGGACGCTCTGGAAAGCATGAGCAGCTCCTAGACAGGTACAGAG  
 AACTCCAGCTCAGTACAGAGAGCAAGGTATCTGAGTTTCTCCATCAGAGCAAGTTGAAGTCTTTGAAAAG  
 TGAGCGTGTCAACTCCTGCAAGAGGAACTGCAAGAAATCTCACGCAGTGCCAGTTGGAGTGTAAAAA  
 TATCAGAAGAAATGGAGGTTTTAACTAAAGAATTTATAGTCTCAAACCTTCTTCTGAAAAACGCATTA  
 CCGAACTCGAGGCACAGAACTCAGAGCATCAGGCAAGGTTAGACATTTACGAGAAGTTGAAAAGGAGCT  
 TGATGAGATAAATGCAAAACAGCAGAAATGAAAATGAAGATGAGGCTGAAAGAATTCTTTATTCCTAT  
 GGTATGGTGCTAATGTTCCCAACAGCTAAAAGACGACTAAAGCAAAGTGCCACTTGGCAAGAAGAG  
 TTCTTCAGTTAGAAAAACAAAATTCATTAATTTTAAAAGATCTGGACCATCAAAGAACCAAGTAAGACA  
 GCTTTCACAAGAGCTTGACAGAGCCAATCACTGTTGAATCAGACTCAGCAGCCCTACAGATACCTCATA  
 GAGTCTGTGCGACAGAGAGATGCCAAGATTGACTCACTGATGAAGTCTACAGCTCAACTTGAAGAAAGATG  
 TCAGCAACTTAAATAAAGAGAAGTCAAGCCTGCTGCAGACGAAGAACCAGATGGCACTGGATCTGGAGCA  
 GCTCCTCAGTCACCGGAGGAATTTGCAGCTATGAAGCAGATCATCATTAAATATGTGTAGTAAACATTCT  
 GAGAACAACCTATTTCTACGAAAATGGAATCAAAAAGTGTGACAGAAAATCAAGCAAAGACTTTGAATA  
 TGCCAAGAGAACATGAAGAGAATATATTTATACCCAAGCCAACGCTCTTTACTAAAAAGGAAGCACAAAGA  
 GTGGCCCAAGAGTCAGAAGATGAAGACC

**ACGCGTACGCGGCCGCTCGAG** – GFP Tag – GTTTAA

**Protein Sequence:** >MG212404 representing NM\_029320  
Red=Cloning site Green=Tags(s)

```
MSRKIAKEPKKVNISSSLESEDISLETTIHTDDVSSSEEREGKVKITRQLIERKEILHNIQLLKIELSQK
NMMIDNLKMDYLTKIEELEEKLNDAHQKQLLTRLDNQLTIQQKDAKKYQELMKQEMETILLRQKQLEE
TNHQLREKAGDVRRNLRDFELTEEYVVKLSFPEDQLSIPYVVSIRFYELVNPLRKEVCLEQVKKSELSE
ELSTSKGQLKQLTETYEDRRNNAELLIRCQRLTLELADTKQLVQQGDYRQENYDKVKSERDALEQDVLE
LRRKHEVLEASHIAQAKERNELSKEVSSLQQTVTLLQKDKDYLNRQNMELSVRCAHEEDRLERLQVQLED
TKKAREEMYEKYVTSRDHYKTEYENKLHDELEQIKLKTNLEIDQLRSASREMYERENRNLREARDNALAE
KNRAVAAEKDALGKHEQLLDYRELQLSTESKVSEFLHQSKLKSFESEVQLLQEETARNLTQCQLECEK
YQKKLEVLTKEFYSLQTSSEKRITELEAQNSEHQARLDIYEKLEKELDEIMQTAEIENEDEAERILYSY
GYGANVPTTAKRRLKQSVHLARRVLQLEKQNSLILKDLHDHQNQVRQLSQELDRANSLNQTQQPYRYLI
ESVRQRDAKIDSLMKSTAQLEKDVSNLNKEKSALLQTKNQMALDLEQLLSHREEFAAMKQIIINMCSKHS
ENNLFLTKMESKSVTENQAKTLNMPREHEENIFIPKPTLFTKKEAQEWPKSQKMKT
```

TRTRPLE - GFP Tag - V

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_029320

**ORF Size:** 2268 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_029320.3](#), [NP\\_083596.2](#)


**RefSeq Size:** 3271 bp

**RefSeq ORF:** 2271 bp

**Locus ID:** 52023

**Cytogenetics:** 14 49.27 cM

**Product images:**



Circular map for MG212404