

Product datasheet for **MG212529**

Gm436 (NM_001085504) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Gm436 (NM_001085504) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Gm436
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG212529 representing NM_001085504 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCAGTTCTGTGCTGCTACTGCTAGTTGGGCTGCCTACATTCTTACTGGGAGTCTTTATATGGGCTG
TGGTAAAGCATTTTATCACAGCAGAAATTCCTTCTTTACAGCATCCAGTCAAGTTCAGGTGTATGCA
TTGTCTGATGCTTTCATAATGGCCTTGGGTAATATATGCAGAAGCTGAGAATCTGCTCTCTGCCCAGC
TTCATTATATTTATGCATGATTTAGTCCCATAAGGAAGGACCCTACCTTAGTGGTGACAAACATGCATT
TGGGGACAATCCCTGTGAGACTGTTCCAGCCCAAGGCAGTGTCCCTCCAACTCCGGAGGGGCATCATCTT
CTACCATGGAGGGCTGCAGCTTTTGGTAGCCTGGACATTTACCATAACCTGTGCAGTTTCTGGTCCGG
GAGACAGACTCAGTGTTGCTTTCAGTGGGTACCGCAAGCTTCCCTAAGCATCATCACCCCTGCTGCTCTTA
ATGATTGCCTGAGTGCCACCATCTGTTTCTGAAGGCCCTGAAACTTATGGAGTGGACCCATCTCGGGT
GGTGTCTGTGGAGATAGTTTTGGAGGTTGGGTTGTAGCCTCCATCTCCAGACCTTGGTCAGTACCCT
AGTCTCCCCAGATCCGGGCTCAAGTCTTGATTTATCCAATTTGCAGGTCATTAATTTCTTTTACCAT
CTCTACAACAGAATAAAAAATCCCACTCTTTCCAAAGACTTGTGATTATATTTGTGTGTAACACTT
GGCCATTGATCTCTTTGGGAGATTCCATGCTTACAGGTGCTGTATACCTCTGGACAAGTGAAGAAA
TACAGAAAATGGCTCAGCTCTGACAATATCCAGAAGGCTCTGGAGCCAAGACACAACTGAGTTTC
TTGGACATTTAATGAGGCTGCCTATCTAGAAACAAAACACATTTTAAATCCAGAAATTTACCTATCAT
CGCAGATGATAAGACCATTGCTCAGCTTCTGAGGCATTCTAGTGAGCTGTGAGTATGACATAGTTCGT
GATGACACCTTGCTTTACAAGAAGCGCTTGAAGATCAGGGAGTTCCTGTGACCTGGCACCATGTGATGG
ATGGCTTTCATGGATGCGTCTTCTGTTGATAGAAAGTTTCTATCTTTCCCTGCTCCATAAAAGTAGT
AAACACTGTTGTCAGTTGCATAAAAAGCATA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG212529 representing NM_001085504
 Red=Cloning site Green=Tags(s)


MAVLCLLLLVGLPTFLLGVFIWAVVKHFITAEIPSSLQHPVKFRMCHCLMLFIMALGNILQKLRICSLPS
 FIIFMHDLVPIRKDPTLVVTNMHLGTIPVRLFQPKAVSSKLRGGIIFYHGGAAAFGLSDIYHNLCFLVR
 ETDSVLLSVGYRKLPKHHHPAALNDCLSATILFLKALETYGVDPSSRVVLCGDSFGGWVVASISQTLVSTP
 SLPQIRAQVLIYPILQVINFLPLSLQONKNIPLLSKDLLIIFVCNYLAIDL SWGDSMLTGAVIPLDKWKK
 YRKWLSNDNIPRRLWSQDTQPEFLGHFNAAAYLETKHIFNPEISPIIADDKTIAQLPEAFVLSCEYDIVR
 DDTLLYKKRLEDQGVPTWHHVMDGFHGCVLLFDRKFLSFPCSIKVVNTVVSCIKSI

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



Plasmid Map:



ACCN: NM_001085504

ORF Size:	1221 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001085504.1 , NP_001078973.1
RefSeq Size:	1224 bp
RefSeq ORF:	1224 bp
Locus ID:	230890
Cytogenetics:	4 E1