

## Product datasheet for **MG212689**

### Cd209f (NM\_026956) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Cd209f (NM\_026956) Mouse Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** Cd209f  
**Synonyms:** 1810029C22Rik; AI647041; SIGNR8  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >MG212689 representing NM\_026956  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGGCAGGAGAGCCAGAGTTCCATCAGTCAAAGAATCATGAGGAGGAGGTCACATCTGGAGAGCAGGGAT  
 TTGCTGAGCATGACCCCGAGGGTCTCATTTGCAGCTCCAAGAGTTTACAAGGGCATCTGACCCAAGTCCC  
 CTGGCTTCTGCTGCTTCTCATCTCTTTGGGCTCTTTTGGCTTATGTTGGCCACCCTGGTTCAAGTCTCC  
 AGGATTCGTGCATACTCACAGGGACAGACCCAGGATCAGCAGGGTAGCTCCAGCTTGACAAGGTTGCCG  
 TTCTCCGGAGCAGACACTCCGGCTTGAGCAGATCCAGCAGATTCAGCAGCAGCTGACTCAGTTCAA  
 CGCCTCGCTGGCTGGCCTGTGCCGGCCCTGTCCCTGGGACTGGGAGCTTCCAGGGAAGCTGTACCTC  
 TTCTCCAGAACCCTGGGCAGCTGGGAAACCTCTGCCTCCTCCTGCGAGGATCTCGGGGCCACCTGGTGA  
 TTGTCAACAGTGTTCAGAGCAGAGATTCATGAAATATTGGAATGTGAGGAAAAACCAACGCTCCTGGAT  
 CGGCCTCAGTGACCATATACATGAAGTTTCTGGCAGTGGGTGGACGGCAGCGCCCTCAAGTTCAGCTTT  
 TGGAAAGAGGGGAAACCTAACAATGACGGAGATGAGGACTGCGTGGAGCTGTTTCATGGATGATTGGAATG  
 ATAATAAGTGCAGTGAACAGAACCTTTGGGTCTGTGAGCAGCCCTCTGCTCCCTGCCCTCATCAC

**ACGCGTACGCGGCCGCTCGAG** - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >MG212689 representing NM\_026956  
Red=Cloning site Green=Tags(s)

MAGEPEFHQSKNHHEEVTSGEQGFAEHDPEGLICSSKSLQGHLTQVPWLLLLLISLGLFLLMLATLVQVS  
 RIRAYSQQQTQDQQGSSSLDKVAVPREQTHSGLEIQIQQQLTQFNASLAGLCRCPWDWELFQGSCYL  
 FSRTLGSWETSASSCEDLGAHLVIVNSVSEQRFMKYWNVRKNQRSWIGLSDHIHEGSWQVWDGSAKLFSE  
 WKEGEPNNDGDEDCVELFMDDWNDNKCTEQNFVWCEQPSAPCPHH

TRTRPLE - GFP Tag - V

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**

Cloning sites used for ORF Shutting:



**ACCN:** NM\_026956

**ORF Size:** 765 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_026956.2](#), [NP\\_081232.2](#)


**RefSeq Size:** 1109 bp

**RefSeq ORF:** 768 bp

**Locus ID:** 69142

**Cytogenetics:** 8 A1.1

### Product images:



Circular map for MG212689