

## Product datasheet for **MG212707**

### 1700029F12Rik (NM\_025585) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** 1700029F12Rik (NM\_025585) Mouse Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** 1700029F12Rik  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >MG212707 representing NM\_025585  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGAGCTCCAGGCAGGAATCTGACATTGGCCATGACTGGGATGTAAGCTCAGACAGGAATCTGACATCAG  
GGCATGATAAGGAAGACAGCTGTGGAGCTCAAGCGCACCCAGTGTGACAGTGAAGCCACCACTGACCCG  
GAGAACCAATGTCACATGGACGGCCAACACTGCCTGTAGCCACTGCTATACTGCGTGGACTCCATCCGCT  
GACACCTGCCATCATTACCACGTGGATTGTGAGTGTGCTGCAGGCCATGGCTGCTGCGGCCACGACTG  
CTATGAGACTCACAGACTTCTACCACGATTACCAGCCACTGTCCACAGAAGACCAGCTGTCCACAGAGGA  
ACCTGATGCTAATGTTTACCAAGTGGATCCTGATGGTAAAAGTCAGCTACACCCAAAACCTGACGACGCC  
CCAAGACAGCAGGAGGCAGCT

**ACGCGTACGCGGCCGCTCGAG** - GFP Tag - GTTTAA

**Protein Sequence:** >MG212707 representing NM\_025585  
Red=Cloning site Green=Tags(s)

MSSRQESDIGHDWDVSSDRNLTSGHDKEDSCGAQAAPSVTVKPPLTRRTNVTWTANTACSHCYTAWTPSA  
DTCHHSRGLSVLLQAMAAAATTAMRLTDFYHDYQPLSTEDQLSTEEPANVYQVDPDGESQLHPKTD  
DPRQEEAA

**TRTRPLE** - GFP Tag - V

**Restriction Sites:** Sgfl-MluI



[View online »](#)



<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_025585.3</a></u>
<b>RefSeq Size:</b>	1130 bp
<b>RefSeq ORF:</b>	444 bp
<b>Locus ID:</b>	66479
<b>Cytogenetics:</b>	13 D1