

Product datasheet for **MG212713**

1700017N19Rik (NM_001081246) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	1700017N19Rik (NM_001081246) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	1700017N19Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG212713 representing NM_001081246 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAATGCAGCAAACTGCAGCATTTTCATGCTTCTGGAAACCAACCTCTGGTTGTGTGAAGATCA
GTTGTATCTTTTACCACAGCAAACCCGAAATATCAATGGATTGTTTTGCCACCAAGTAGCAATACCAC
CTTACAGAAAGAAAGTCAGGAAGGAAGCCCATCCCAACCCAGAGTCAGGAATCTCTGAAACCCATGGAG
AACGTATCAAGACCCATTTCATCCATTAGTTTTAAAACTAACTTTGAAGAAGAGGAGGAGGAAG
AGGAACCAATGATGCTTCTAGTTTATGGACAAAGACCCCTGAAGAAATTGAAGAGAAAAGAGCAATAAA
GGAGATGTGTTATAAATCTGGTGAATACTATAGATTCCACGCTCCTCCAGACATCCCATCATCAAAAAGC
GTAGCTCCGACAGTGGAGAAAGAGTTAGAAAAGCCTTTGGAACTGGCAGTGAATTACAAGAAGGAGATG
GTCTTACAGTTCCAACAAAGTTTAGCCTATTTGAAAGGGCAGGAGAGGCAAAAACATCATTGGATAGGAA
ACCAAGGACTGACATTGCTGCTTTTAAAATGGAGGAGGTGACTGCTATGCTCCACAGAGGATCATATTT
CTTGGAGTGGATGAAAGTGAAGCATTAGCTGAAGACAAAGAAGCCACCTCAAAATGCTCAAATGTCAAAG
AAGAAAGCAAGAATGGTCTTACCCGAAGCATCCTCTAACTACTCGCCTCATGCCTACAACACATGTGTT
AAATGCCACTGAGAATATCAGCATGAAGTGCAGAGAGACCCCTCTTCAATGAATGATGTCCAGCCAGTG
AGGAAGCCTCATTTTAAAGGTGTGAAGAAAAGAAAATGGATTTATGATGAACCAAGAAGCTTTCTGGCT
CTGGAATGCGGAGAGCAGTGCACATCCCAAATCCCAAACATAAAATGAGTTACCATCATATAATAAAAA
CCGGAATGCCGAGAACGCTCCTATACGCATGCTCCAAGAGATGCAGTAAGAAGTGTCACTGAACGCC
CCTCCCGAAGCCGGCCACGAGCCGGTCTACAACAAAGCTGATGTCAACAAGGAACCCAACTCAACC
TCTGTCCAGATAAACATATGTCAACATCGTATAATGGTTCAGCCTGGCGAAGAAGAATTCCTTTTTCAA
AACTTATTCAAAGACGAAAAGGTATATCCAGAGCCAAGAAGAAATGGGAGTAAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)


Protein Sequence: >MG212713 representing NM_001081246
 Red=Cloning site Green=Tags(s)

```
MEMQQNCSISCFWETQPLGCVKISCFYHSKPRNINGLFLPPSSNTTLQKESQEGSPPTQSQESLKPME
NVSRIHHPLVLKTNFEIEEEEEEPNDASSLWTKTPEEIEEKRAIKEMCYKSGEYYRFHAPPDIPSSKS
VAPTVEKELEKPLETGSELQEGDGLTVPTKFSLFRAGEAKTSLDRKPRTDIAAFENGGDCYAPQRIIF
LGVDESEALAEDKEATSKCSNVKEESKNGLHPKHPLTTRLMPHTHVLNATENISMKCRETPSSMNDVQPV
RKPHFKGVKRRKIYDEPKNFPGSGMRRRAVHIIPNPKHKMSYHHHNKRNNAENASYTHAPRDAVRTVTLNA
PPRSRPTSRSYNKADVNEPKLNLCPDKHMSTSYNGSAWRRRIPFSKTYSKTEKVYPEPRRNGSK
```

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001081246

ORF Size: 1245 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001081246.1](#), [NP_001074715.1](#)


RefSeq Size: 1690 bp

RefSeq ORF: 1248 bp

Locus ID: 66605

Cytogenetics: 10 D1

Product images:



Circular map for MG212713