

Product datasheet for **MG212899**

Olfr107 (NM_146511) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	Olfr107
Synonyms:	MOR156-3
Mammalian Cell	Neomycin
Selection:	
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >MG212899 representing NM_146511
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**

ATGAACTGTAGTAAGACTCCTGGCTTTATCCTGTTGGGACTTTCCAGTGACCCAGAGAAATGGCAGCCCC
TCTTTAATATTTTCTCTGTCTCTACTTGCTGGGCCTTTTAGGGAACCTGCTACTTCTGTTAGCTATTGG
CACTGATGTCCACCTCCACACCCCATGTATTTCTCCTCAGTCAGCTCTCACTAGTGGATCTTTGCTTC
ATCACCACAACAGCCCCATAAATGCTGGAGGCTCTGTGGACTGGAGATGGATCAATCTCATTCTCTGGAT
GCCTGACTCAGTTTTACTTCTTTGCTGTTTTGCTGACATGGATAACCTGCTTCTGGCAGTCATGGCTAT
TGACCGCTATGCTGCTATCTGCCACCCCTGTTCTACCCATTTCTAATGACTCCTTGTAGATGCGAAGTT
CTTGCCAGTGGGTCATGGGGAATAGCTCATTGTGTCTCTGTTCTATACCTTATTACTCTCTCAGTTTT
ATTATCATACCAATCAAGGGATACCTATTTTTCTGTGATTCTAGGCCTCTTACTGCTTTCTTGCTC
AGACACTCACCTCAGTGAGGGCCTGATGATGGCTTTGTCTGGAGTTTTGGGAATGAGCTCAGTTCTCTGT
CTTGTAAGTTCTTATGTTTGTATTTTTATGCTGTGGCTAGGGTTCATCAGCACAGGGGAAAAGGAAAG
CCCTGGCCACATGCAGTCCCACCTCTCTGTAGTCCTCTATTCTATAGCACAGTCTTTGCTACCTACCT
GAAGCCCCATCTACTTCTCACTCCTCTGCAGAGGTGGTAGCTGCTGTTATGTATACCTTGTTACTCCC
ACTCTGAACCCCTTCATTATAGTCTTAGAAATAAGGATGTTAAGAGTTTCATTGAGAAAGATCCTGAACA
TGGACAAGTTTCAGGGC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



Protein Sequence:

>MG212899 representing NM_146511

Red=Cloning site Green=Tags(s)

MNCSKTPGFILLGLSSDPEKWQPLFNIFLCYLLGLLGNLLLLLAIGTDVHLHTPMYFFLSQLSLVDLCF
ITTTAPKMLEALWTGDSISFSGLTQFYFFAVFADMDNLLAVMAIDRYAAICHPLFYFPLMTPCRCEV
LASGSWGIAHCVSLFYTLSSQFYHYTNQGIPIHFFCDSRPLLLSCSDTHLSEGLMMALSGVLGMSSVLC
LVSSYGCIFYAVARVPSAQGRKALATCSSHLSVLLFYSTVFATYLPKPSTSHSSAEVVAVMYTLVTP
TLNPFYISLRNKDVKSSLRKILNMDKFQG


TRTRPLE - GFP Tag - V

Restriction Sites:

SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



CTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCCGATCGCCGCCGCCAGATCT

EcoRI BamHI KpnI RBS Kozac Consensus SgfI AscI

CAAGCTTAAGCTAGCTAGCGGACCG ACG CGT ACG CGG CCG CTC GAG ATG GAG AGC GAC --- --- ---

HindIII NheI RsrII MluI NotI XhoI GFP Tag

T R T R P L E M E S D - - -

--- GAA GAA AGA GTT TAA ACGGCCGCCCGGAGCT

- - E E R V Stop

ACCN:

NM_146511

ORF Size:

927 bp

OTI Disclaimer:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:

This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:

The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_146511.2](#), [NP_666722.2](#)

RefSeq Size: 1033 bp


RefSeq ORF: 930 bp

Locus ID: 258504

Cytogenetics: 17 B1

Gene Summary: Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]

Product images:



Circular map for MG212899