

Product datasheet for MG212927

Aadacl3 (NM_001085503) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Aadacl3 (NM_001085503) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Aadacl3
Synonyms:	Gm435
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG212927 representing NM_001085503 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCCGCATCGCC

ATGGTGGTTCTGGCCCTGACACTCCTTGTAGGTTCAAGTTCAGTTGCTGTGTTCTCACTGGGCAGTCTGCTGTGGG
TGGTCGGCAAACACTTCTGGACTGAACACATCCCTGAAGGCATCACTCACCCCTGGAGGCTTCGAATCCT
CAGTTGCTTATTTTCATCTGACAATGACATGGGGGATGATTTTAAAACTGGGACTCTGTTATGCACCC
CAGTTTCCAGTTCATGATCTGAAACCACTGAAGAGAGATCCTGATGTTGTGGTGAAGATCTTC
ACTTTGGGACAATCCCTGTGAAGCTTTACAAGCCCAAGAAACCTCCTCCATCCCAAGGCTTGAATCAT
TTTCTTCATGGTGGTGAACCATCATAGGGAGCTTGAGAACCCACAACAGCATATGCTTGCGTTTGTCC
AAGGAGTGTGACTCTGTTGTCGTTTCAGTTGGATACAGGAAATCTCCCATGTACAAGTACCCAGTGATGA
AAGATGACTGTGTAGTAGCCACTACCCACTTTCTGGAATCATTGGATGTGTATGGAGTGGATCCAGCTCG
AGTGGTGACCTGTGGTGACAGTGTAGGAGGAAGTGCAGCCACTGTGACTAGTCAAATGTTAGTGACCCGA
CCAGATCTCCCCGGATCAAGGCTCAGATTTTGATCTACCCCTTACTGCAATTGATAGATTTTGGGAGCC
CCTCTTATCAACAGAATAGAAACATTCCTGCTCTCCTGGGATCTGGCCTTCTACTGCTTCTGTTGCCA
CCTAGATGTGAACATTTCTGGAAAAGTGTGTAAGAATGGGATGCATTTGCCTCCGGACGTCTGGGAG
AAATACAGAAAAGTGGCTAGGTGCAGAAAACATCCAGAGAGGTTTAAGAATCGAGGCTACAAATCTATAC
CCTGGGGTCTGTGAATAATGATGCTTACCAAGAGATAAAACGTTCACTGAACTATACCTGCTCACCCCT
CATCTCAGAAGACAGTATAGTATCTCAGCTCCCAGAAACCTGCATTGTGAGCTGTGAGTATGACCTTCTC
CGGGACCATTTCTCTTATATAAGAAGAGGCTAGAAGTCTGGGGTGCCTGTACGTGGCATCATATGG
AAGATGGCTTCATGGAGTGCTTAGTGCCCTTGACTATGGATTATTGAGCTTCCCTTGCCTCGAGAAT
TATGGACCTGATTATTCAGTTCATTAGAAAGTTC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA


[View online »](#)


Protein Sequence: >MG212927 representing NM_001085503
 Red=Cloning site Green=Tags(s)

MVVLALTLLVGSVAVFSLGSLLWVVGKHFWEHIPEGITHPWRLRILSCLFHLTMTWGMIFEKLGCLYAP
 QFASFLHDLKPLKRPDPVVVDLHFGTIPVKLYKPKKPSSIPRLGIFFHGGGTIIIGSLRTHNSICLRLS
 KECDSVVVSVGYRKSPMYKYPVMKDDCVVATTHFLES LDVYGVDPARVVTGDSVGGTAATVTSQMLVHR
 PDLPRIKAQIL IYPLLQLIDFGSPSYQQNRNIPLLSWDLAFYCFCHLDVNI SWKSVVKNGMHLPPDVWE
 KYRKWLGAENIPERFKNRGYKSIPWGPVNNDAYQEIKRSLNNTCSPLISEDSIVSQLPETCIVSCEYDLL
 RDHSLLYKKRLEDLGVPVTVHHMEDGFHGVLSALDYGLLSFPCASRIMDLIIQFIRKF

TRTRPLE – GFP Tag – V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001085503

ORF Size: 1224 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001085503.2](#), [NP_001078972.1](#)

RefSeq Size: 3153 bp


RefSeq ORF: 1227 bp

Locus ID: 230883

UniProt ID: [A2A7Z8](#)

Cytogenetics: 4 E1

Product images:



Circular map for MG212927