

## Product datasheet for **MG212988**

### **Olfr97 (NM\_146512) Mouse Tagged ORF Clone**

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Olfr97 (NM_146512) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Olfr97
Synonyms:	MOR156-2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG212988 representing NM_146512 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGATGAAGTGCAGTCAGGCTCCTGGCTTCATCCTGTTGGGACTGTCCAGTAACTCAGAGAAATGGCAAC  
CCCTCTTAGCATCTTCCTTGTCTCTACTTGTGGCCCTTTAGGGAACCTGCTACTTCTGTTAGCTAT  
TGGCACTGATGTCCACCTCCACACCCCATGATTTCTCCTCAGTCAGCTCTCACTCGTGGATCTTTGC  
TTCATCACCACCACAGCCCTAAAATGCTGGAGACTCTGTGGACTGGAGATGGATCCATCTCGTTCTCTG  
GATGCCTGACTCAGTTGACTTCTTTGCTGTTTTGCTGACATGGATAACCTGCTTCTAGCAGTCATGGC  
TATTGACCGCTATGCTGCTATCTGCCACCCACTGCTGTACCCACTTCTAATGACTCCTTGTAGATGTAGA  
GTTCTGGTCAGTGGGTCATGGGAGTAGCTCATTGTGTCTCTGACCCATATTTGTTGCTCTCTCAGT  
TATATTTTCATACCAATCAAGAAATTCCTCATTTTTCTGTGACTTTGGCCCTCTTATTGCTTTCCTG  
CTCTGATGCTCACCTCAATGAGAGTCTAATGATGGCTTTGGCTGGAGTTTGGGAATCAGTGCACCTGCTC  
TGCATTGTAAGTCTTATGGTTGATTTTTATGCTGTGGCTAAGGTTCCATCAGCACAGGGAAGAAAGGA  
AAGCCTGGCTACATGCAGTCCACCTCTCTGTAGTCCCTCTATTCTATAGCACAGTCTTTGCTACCTA  
CCTGAAACCCCATCTAGTCTCGTTCCTCTGGGAGGTGGTAGCTGCTGTTATGTATACCTTGGTCACT  
CCCACTAAACCCCTTCAATTTATAGTCTGAGAAACAAGGATGTTAAGAGTTCATTGAGAAGGATCTTGA  
ACATGGTGAATCTCAGGAC

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >MG212988 representing NM\_146512  
Red=Cloning site Green=Tags(s)

MMNCSQAPGFILLGLSSNSEKWQPLFSIFLVLYLLGLLGNLLLLLAIGTDVHLHTPMYFFLSQLSLVDLC  
 FITTTAPKMLETLWTGDGSI SFGCLTQLYFFAVFADMDNLLLAVMAIDRYAAICHPLL YPLLMTPCRCR  
 VLVSGSWGVAHCVSLTHILLLSQLYFHTNQEIPHF FCFGLLLLSCSDAHLNESLMMALAGVLGISALL  
 CIVSSYGCIFYAVAKVPSAQGKRKALATCSSHLSVLLFYSTVFATY LKPPSSSRSSGEVVAAMVTLVT  
 PTLNPF IYSLRNKDKVSSLRRLNMVKSQD

TRTRPLE - GFP Tag - V


**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**

Cloning sites used for ORF Shuttling:



**Plasmid Map:**



**ACCN:** NM\_146512

**ORF Size:** 930 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_146512.1</a> , <a href="#">NP_666723.1</a>
<b>RefSeq Size:</b>	933 bp
<b>RefSeq ORF:</b>	933 bp
<b>Locus ID:</b>	258505
<b>Cytogenetics:</b>	17 B1
<b>Gene Summary:</b>	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]