

Product datasheet for **MG213124**

Tpgs2 (NM_001142698) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Tpgs2 (NM_001142698) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Tpgs2
Synonyms: 5730437P09Rik; 5730494M16Rik; AI666318; Pgs2
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG213124 representing NM_001142698
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAGGAGACAACGCCCCACTGCAAGCCGGCAGCAAGCCGCACCTTGAGAAGTTGACCCTAGGCGTCA
CTCGCATCCTGGAATCTTCCCAGGTGTGACCGAGGTGAGCATAATAGAAAAGCTGCCTGCTGAGCGTCA
TATGATTTCTTCATGGGAACAGAAAACAATTGTGTGATGCCTGAGGATGTGAGAACTTCTACTTGATG
ACCAATGGTTTTACATGACCTGGAGTGTGAAGCTAGATGAGCACATTATCCACTGGGAAGCATGGTGA
TTAACGGCATTTCAAAAGTACTGACTCAGCTCATCCAGTCTTCTGTGTACTCACTTCCTAATGCACCAACTCT
GGCAGACCTAGAGGACGATACTCAAGAAGGAAATGAAGACCACCAGCTAGAGAAGCCTCACTTTGACTGT
CGCAGTGAATATTTGAGCTGGATTCATGCAATGGGAATGGGAAAAGTTTGCCTCGTCTACAAAAATGGGA
AACCTGGGTTAGCTCATGACACTGAGATATGGTTCCTGGACAGAGCATTACTGGCACTTTCTCACGGA
TACCTTTACCGCCTATTACCGCCTGCTCATACCCACCTGGGCCTGCCCCAGTGGCAGTATGCCTTTACC
AGCTATGGCATTAGCCACAGGCCAAGGTAGCTTCTATCTATAGTCAAGAAGAAAGCATGAATTTCCAAG
CAACATCTGTTAGAAACCAATGG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG213124 representing NM_001142698
Red=Cloning site Green=Tags(s)

MEETTPPLQAGSKPHLEKLT LGVTRILESSPGVTEVSIIEKLP AERHMISSWEQKNNCVMPE DVRNFYLM
 TNGFHMTWSVKLDEHI IPLGSMVINGISKLTQLIQSSVYSLPNAPT LADLEDDTQEGNEDHQLEKPHFDC
 RSAIFELDSCNGNGKVC LVYKNGKPLAHDTEIWFLDRALYWHFLDTFTAYYRLLITHLGLPQWQYAF
 SYGISPQAKVASIYSQEESMNFQATSVRNQW

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:



ACCN: NM_001142698

ORF Size: 723 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001142698.1 , NP_001136170.1
RefSeq Size:	3607 bp
RefSeq ORF:	726 bp
Locus ID:	66648
UniProt ID:	Q66JT5
Cytogenetics:	18 A2