

## Product datasheet for **MG213429**

### Poteg (NM\_026256) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Poteg (NM\_026256) Mouse Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** Poteg  
**Synonyms:** 4921537P18Rik; 4930456F22Rik; Gsarp1; Gsarp2  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >MG213429 representing NM\_026256  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGTCCTTTTCCATGAAGAAGATCTTTGGGGCTAAGGAAAAGACACCCTTGGGTTTCTGTGATGCCCCAC  
AAATGACTGTCTCTGAGTTTTTCTCGGTGAAAATCCAACTATCCGAAATACCATACAACATACAGACC  
GTTGGGGCATATTCACAGGGTTGCTGCTGAGGGTATGCAGCCCGGATGGAGATTCTCTCACGCTTGA  
CAATGCAACGTGTACCACAGAGACCGCAAGGACAGGACAGCTTTGCACTTGCCTGTGTCTATGGCCGGC  
TACCAGTAGTGACCATTCTAGTGAAGAATAATTGTGAGATCGATGCCTTGGACAAGAACCATACCACACC  
TCTGATGAAGTCTGTCCAATGCTGGAAGCAGAAATGTGCGACCGTGTCTGCTGGAGCACGGCGCTGATCCG  
AACATTAGAGACAGCAGTGGCAACAGTGCCTTGCCTACGCTGTGTACAACGGGCATCAGGAGATGGCCT  
CCCTACTGCTTACGTACAACGCTGATATCGAACAAAAACCAAGGATGGCTTACACCACTGTTGTGCTGGC  
ACTCCGAGAAAAGCGGGTGGAGGTGGCAGAGTTTTTAGTAAGAATGGGAGCAGACATACAGTGGTAGAT  
GAGCTGCAAGAAAACACCTTGATATACGCTATAAGATGTGGATCAAAAGATCTGTCTGTATTACTCCTTG  
AGAAAGTATTGACTTCTTCTATAAGGATGTCTTCGGATGGACAGCATTACGCTACGCTATTGAAGGCCA  
CTGTACTTTTAGACAGACACTTCTGGATTTTGAAGAAAGTCTACATAGTAACAAAAAAGATAAGGAACCA  
GAACTACAAGTGGAGCCAGTACCTCCTGCATTTAGCAAAACAAATAGATGCTGGAATGACTCACTAG  
CAAGGATCTCCAGCTGTCTAGCCCTGAGACTCCTGTGCTTACAATGAAGGAAGACGAGAACAGTATGC  
CAGGGTAAACTGCTCCTCTGTGGAGTTACTTCTTGGCCTAAGTTGTTTTGTTTATTTC

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >MG213429 representing NM\_026256  
 Red=Cloning site Green=Tags(s)

MSFSMKKIFGAKEKTPLGFCDAPQMTVSEFFFGENPNYPKYHTTYRPLGHIHRVAAEGDAARMEILLTLG  
 QCNVYHRDRKDRDTALHFACVYGRLPVVTILVKNNCEIDALDKNHTTPLMKSVQCWKQKCATVLLLEHGADP  
 NIRDSSGNSALHYAVYNGHQEMASLLLQYNADIEQKTKDGFPTLLLLALREKRVEVAEFLVRMGADIHVVD  
 ELQRNTLIYAIIRCGSKDLSVLLLEKIDFFYKDVFGWTALRYAIEGHCTFRQTLLDFEESLHSNKKDKPE  
 ELQVEPSTSCILAKQIDAGNDSLARISSCPSPETPVLTKEDENSARVNCSSVELLPWPKFVFGLF

TRTRPLE - GFP Tag - V

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_026256

**ORF Size:** 1041 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_026256.2](#), [NP\\_080532.1](#)


**RefSeq Size:** 1440 bp

**RefSeq ORF:** 1044 bp

**Locus ID:** 70952

**Cytogenetics:** 8 A2

### Product images:



Circular map for MG213429