

Product datasheet for **MG213479**

Prps111 (NM_029294) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Prps111 (NM_029294) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Prps111
Synonyms:	1700011K15Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG213479 representing NM_029294 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCCAAACATCAAACCTTTCAGCGGAAGCTCCCACCAGGACTTATCCCAGAAGATTACTGAGCGCCTGG
GCCTCGAGCTCAGCAAGGTGGTCACTAAGAAGTTCAGCAACCAGGAGACCTGCGTGAAATTGGTGAGAG
TGTGCGTGGAGAAGATGTCTACATCGTTCAGAGCGGTTGTGGTGAGATCAACGACAGTCTAATGAACTT
TTGATCATGATCAATGCCTGCAAGATCGTTCAGCCAGCAGGGTGACTGCAGTCATCCCGTGCTTCCCTT
ATGCCCGCAGGATAAAAAGGATAAGAGCCGGGCTCCCATCTCTGCGAAGCTCGTAGCAAATATGCTGTC
TATAGCAGGTGCAGATCACATTATACCATGGACCTACATGCTTCTCAGATCCAGGGTTTTTTGACATC
CCGGTAGACAATTTGATGCAGAGCCAGCTGTGCTGAAGTGGATAAAGGAGAATATTTCTGACTGGAGGA
ACTGCACGATTGTCTCTCCTGATGCTGGCGGAGCCAAGAGGGTCACTTCCATCGCGGACCATTTGAATGT
GGATTTTGCCTTGATTCACAAAGAACGGAAGAAAGCCAATGAAGTGGACCGCATGGTGCTTGTGGGTGAT
GTGAAGGATCGGGTGGCTATTCTTGTAGATGACATGGCTGATACTTGTGGTACAATCTGTTACGCGGCGG
ACAAACTTCTCTCAGAGGGAGCTACCAGAGTTTATGCCATCTTGACTCATGGCATCTTTTCTGGCCAGC
TATTGCTCGCATCAATGGCGCATGCTTTGAAGCAGTAGTAGTACAAAACACCATACTCAAGAAGACAAG
ATGAAACAGTGTAGCAAAATCCAGGTGATTGACATCGCCATGATCCTTGTGAAGCCATCCGGAGAACC
ATAACGGAGAATCTGTTTCTACCTGTTTCTCAGTCACGTACCCCTTG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG213479 representing NM_029294
 Red=Cloning site Green=Tags(s)

MPNIKLFSGSSHQDL SQKITERLGLLELSKVVTKKFSNQETCVEIGESVRGEDVYIVQSGCGEINDSLMEL
 LIMINACKIASASRVTAVIPCFPYARQDKKDSRAPISAKLVANMLSIAGADHIIITMDLHASQIQGFFDI
 PVDNLYAEPVAVLKWIKENISDWRNCTIVSPDAGGAKRVTSIADHLNVDFALIHKERRKKANEVDRMVLVGD
 VKDRVAILVDDMADTCGTICAYAADKLLSEGATRVYAILTHGIFSGPAIARINGACFEAVVVTNTIPQEDK
 MKQCSKIQVIDIAMILAEAIRRTHNGESVSYLFSHVPL

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



Plasmid Map:



ACCN: NM_029294

ORF Size: 954 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_029294.2 , NP_083570.1
RefSeq Size:	1676 bp
RefSeq ORF:	957 bp
Locus ID:	75456
Cytogenetics:	12 A3