

Product datasheet for **MG213624**

Btnl10 (NM_138678) Mouse Tagged ORF Clone

Product data:

| | |
|---------------------------|---|
| Product Type: | Expression Plasmids |
| Product Name: | Btnl10 (NM_138678) Mouse Tagged ORF Clone |
| Tag: | TurboGFP |
| Symbol: | Btnl10 |
| Synonyms: | BUTR-1; Butr1 |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-AC-GFP (PS100010) |
| E. coli Selection: | Ampicillin (100 ug/mL) |
| ORF Nucleotide Sequence: | >MG213624 representing NM_138678 Red=Cloning site Blue=ORF Green=Tags(s) |

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCAAGGGCACACCCTGGAGATGCTACTCTGCCAAGCATCTTGGTCAGCTTCATCTTCTGCAACTCC
TGACTTCTGGGAACGGAAAATCAGATTTCTCGTCCTTGGACCTCCCCACCCACTTCTGGCCATAGTTGG
TCAAGACAAAGAGTTGCCCTGCAAAGTCACTCAACATCAGCGCTGAGGGCATGGAGCTGAGGTGGTAC
AGGGACAAGCCCTCGTCAGTCGTGCATGTGTACAAGAATGGGGAGGATGTGTATGACGAGCAGATGGTGG
AATATAAGGGAAGGACCTCATTCAACGGCAGCCATGTAGCCAGGGGAGAAGCTGCTGTGAAGATCCACAA
CGTCACCGTGTTTGACAACGGGACTTATCACTGTGTGTTAAGGAGTACACATCCCACAGCCAGGCCACC
CTGTGGCTGAAGGTGGCAGGGCGGGGCTCGTCCCCAGAATCCGTGTGACTGACACCCAGGACAAAGGCA
TCCGGGCCGAGTGACCTCAGCAGGCTGGTACCCGAGCCCAAGGTGGAGTGGCTGGACCTTAAAGGACA
ACCAGTGTCTGCTGAGAGTCACTTCTCAGTTTCAGCCAGCACTGGCCTCGTGGCCCTGCTCTCCATTGTG
ACTCCCAGGATACGGCTGTGGGGGCTGACTTGCTCCATCTCCAATCCCCTCTCCAGAGAAGGACA
CTGGATTCCTTGACCGCTAGTTAAGGTATCTGTGAGCGGAGCTCACACAGGGAACATTGGTCAGAGTGT
ACAGTCACATGGCAGCATCATCAAAAGTTCAGAGTCTTTCTCTGTAAAAGTTCCC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG213624 representing NM_138678
 Red=Cloning site Green=Tags(s)

MARAHPGDATLPSILVSFIFLQLLTSGNGKSDFLVLGPPHLLAIVGQDKELPCKLSLNISAEGMELRWY
 RDKPSSVVHVKNGEDVYDEQMVEYKGRTSFNNGSHVARGA AVKIHNVTVFDNGTYHCVFKEYTSHSQAT
 LWLKVAVRGRSSPRIRVTDTDKGI RAECT SAGWYEPKVEWLDLKGQPVSAESHFSV SASTGLVALLSIV
 TPQDTAVGGLTCSISNPLLPEKDTGFLAAVVKVSVSGAHTGNIGQSVQSHGSI IKSSEFSVKVP

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:



ACCN: NM_138678

ORF Size: 825 bp

| | |
|-------------------------------|---|
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | NM_138678.2 , NP_619619.2 |
| RefSeq Size: | 1028 bp |
| RefSeq ORF: | 828 bp |
| Locus ID: | 192194 |
| UniProt ID: | Q9JK39 |
| Cytogenetics: | 11 B1.3 |