

## Product datasheet for **MG213939**

### **Klrc2 (NM\_010653) Mouse Tagged ORF Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Klrc2 (NM\_010653) Mouse Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** Klrc2  
**Synonyms:** NKG2C  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >MG213939 representing NM\_010653  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGAGTCACCTGCTTGGAACTGAACAGGCCTCAGAACCCAGGAAAGCCACAAAGAAAACCCAGGCGCAGT  
TCCGTTTCAGTAATGCAGCAGGAGATCGCCTGCTGGAGTTTAGCCTTCAACATGCTTCCCAGAAGCATCT  
GAAGGCCAGCAGGCATGGCTACTGCAAAAATTTGCATCACCTCCAGAGAAGCTCATCGTGGAATCCTG  
GGCACCATCTGGTTTACCTTATTGATTGCTCTAGTAATCTCAACCAGAATTGTTTCTCCATATGCTGAAC  
TGAAGAAGCAGATCCACCCCTTCTGGTTAGAACCAGAAAGGCTCACACGCCTTGCTCTCAAAGCACA  
GCCTTGTCTCATTGTCCAAAGGAGTGGATATCATATTCTCACAATTGTTATTTTCATTGGTATGGAGAGA  
AAATCTTGGAAATGACAGTTTGGGGTCTGCCTTCCAAAACTGTAGTCTGCTTCACATAGATAGTGAAG  
AGGAGCAGGATTTTCTGCAATCTTTCACTGGTCTCATGGACTGGAATCCTTCGGAAGGGCAGAGGTCA  
GGCCTGGGACTGGAAAAAAGACTCAATTTCAAACCAAAGATAGCAGAGATTTCCCATGATGAATGTAAC  
TGTGCCATGATGTCAGCTTCTGGTCTCACAGCAGACAACGTACAACCTTACATCCATATCTTTGTAAGT  
GTAATTTCCCATC

**ACGCGTACGCGGCCGCTCGAG** - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >MG213939 representing NM\_010653  
Red=Cloning site Green=Tags(s)

MSHLLGTEQASEPRKATKKTQAQFRFSNAAGDRLLLEFSLQHASQKHLKASRHGYCKNFASPPEKLIAGIL  
 GTIWFTLLIALVISTRIVSPYAEKKQIHPFLVRTQKGSALLSKAQPCPHCPKEWISYSHNCYFIGMER  
 KSWNDSLGSCLSKNCSLLHIDSEEEQDFLQSLSLVSWTGILRKGRGQAWDWKKSIFKPKIAEISHDECN  
 CAMMSASGLTADNCTTLHPYLCKCKFPI

TRTRPLE - GFP Tag - V


**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**

Cloning sites used for ORF Shutting:



**Plasmid Map:**



**ACCN:** NM\_010653

**ORF Size:** 714 bp

|                               |   |
|-------------------------------|---|
| <b>OTI Disclaimer:</b>        | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>  |
| <b>OTI Annotation:</b>        | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.  |
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).  |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol> |
| <b>RefSeq:</b>                | <a href="#">NM_010653.4</a> , <a href="#">NP_034783.3</a>   |
| <b>RefSeq Size:</b>           | 717 bp  |
| <b>RefSeq ORF:</b>            | 717 bp  |
| <b>Locus ID:</b>              | 16642   |
| <b>Cytogenetics:</b>          | 6 F3  |