

Product datasheet for **MG214020**

Rarres1 (NM_001164763) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Rarres1 (NM_001164763) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Rarres1
Synonyms:	5430417P09Rik; AI662122
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG214020 representing NM_001164763 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGCAGCCCCGTCAACCCCGATGCCAGCGCTGTTGCTGAGCCTCTGGCTGCTGCCGAGCCTCGCACTGG
CCGCTGCTGCTCACTGAGCCCGCGGACCTGGAATACACCGAAGTCCCCGGCAGCCCGCGGCGTGGGGCT
GCCACGCGGGATCCTGCAGCTGGCGGCGCGCTGCGCTGCACTTCTTCAACTCCGCGCGGGCTCGCCC
AGCGCGCTACGGGTGCTGGCCACCGTGCAGGAGGGCCGCGGTGGTCAACCCTCAAGAAGGCTGCGAAG
TGGACTTGGTATTTAGCACAGAGCAATACAACCCGAGCAGGAAGGTGAACGTTTGGGCAATGTTCTGC
GAGGGTGTTTTTCAAGAATGAAAAACCCAGACCGCCGTTAATGTAACTGCGCAGGCTCTTTGATAAA
GTGAAGAGACAAGAGAAGGATTATCGGCTTTACAAGCAGATGAAGCAATTGAAGACTCCTCTGCACAGCA
TCAGCATACTGACAGTCATGGACACATTGATAATCACTGAGACCCCTGTGGGACTTGGCCTTCCTCGG
CAGCTCATACGTGATGTGGGAGAAGACAACCCAGGTCCTGCACTATTACTTGGTCCAGCTCAGCAGTGTG
GAGCGGCTGAAAACGGATGATGACAGCATCGACTTTGACTTCACTGTTCTGCTTCATGAATTCACACAC
AGGAAATCATTCCCTGCCGTATTCATTGGTCTGGTACCCAGGTAACCCCTTAAAGTCAATTACCACTG
TCAAGAGCAGCAGAGCCAGAAGAAGCCTCTGGAACATGGAAGCATCAGCTATGGCACCGCCAGAGTTC
GGTAATTC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG214020 representing NM_001164763
 Red=Cloning site Green=Tags(s)

MQPRQPPMPALLLSLWLLPSLALAAAVTEPADLEYTEVPRQPAAVGLPRGILQLAARAALHFFNFRAGSP
 SALRVLATVQEGRAWNPQEGCEVDLVFSTEQYNPEQEGERLGKCSARVFFKNEKPRPAVNVTCARLFDK
 VKRQEKYRLYKQMKQLKTPLHSISIPDSHGHIDNSLRPLWDLAFLGSSYVMWEKTTQVLHYYLVQLSSV
 ERLKTDSDSIDFDFTVLLHEFSTQEIIPCRIHLVWYPGKPLKVNHYHCQEQQSPEEASGTMEASAMAPPEF
 GNF

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001164763

ORF Size: 849 bp


OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

- Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
- RefSeq:** NM_001164763.1, NP_001158235.1
- RefSeq Size:** 1411 bp
- RefSeq ORF:** 852 bp
- Locus ID:** 109222
- Cytogenetics:** 3 E1

Product images:



Circular map for MG214020