

Product datasheet for **MG214073**

1700011L22Rik (NM_026315) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: 1700011L22Rik (NM_026315) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: 1700011L22Rik
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG214073 representing NM_026315
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGACACACTTCTTGACTTGACTCCAGAAATCTGCTCCCTTTAGCCCTCTCACTTCTACAGAGTTCCG
 AACTGATCAGACGCAAGGCCCAACAACCTGTGGCAGGATGAGACACGGTGGTCTGCTTCTCCATGACAAC
 CTACTCAGGGAGTTACCGGGAAAAGCAACTGGACGAGGCCACCTGCAACCGGCTGGCCAGAGAGTAGGC
 CAGCCGAGTTTGTAGTACAAACCAACACCATTGCCAGGCAGCTCTGCCTACAACACTCTACCCGGCCAGC
 CTGGTTCCAAGAGGCTGCAGATGGCAAAGGAAGACTCCCTGATATAGCAAGCCCTCCAGGGACAGCCC
 ACTGAACATCAAGCACAAAGTGGCACATCAGATTTGGGGTTCTGAAGTCCCTGTCCAACATTTCTAGCT
 TACAGAGGTCTCCGTATGTCCCCCTGTATGCAGCCAAAAAGCCAGGCTTTGAGCTCCTCATGAGTTACA
 GGAAACAGAGGGAAGGCCCTGCTGAAACGGCTTCAGCGACAATGGGATTATGAAAGCAAGTTGGGCTCCTC
 GGAGGACTCAGGGACGGATCGATTCTCCTCCAATACCAGTGGAAGCTCAGGGCGCAAATTC AAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG214073 representing NM_026315
 Red=Cloning site Green=Tags(s)

MTHFLYL TPEILLPFSPLTSTEFELIRRKAQQLWQDETRWSASSMTTYSGSYREKQLDEATCNRLAQRVG
 QPQFEYKPTPLPGSSAYNTLPGHAGSQEAADGKGRLPDIASPRDSPLNIKHKVAHQIWGSEVPCPTFLA
 YRGLRMSPCMQPKKPGFELLMSYRNRGKALLKRLQRQWDYESKLGSSSEDSGTRDFSSNTSGSSGRKFK

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_026315

ORF Size: 624 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_026315.1, NP_080591.1</u>
RefSeq Size:	978 bp
RefSeq ORF:	627 bp
Locus ID:	67687
UniProt ID:	<u>Q9DAF8</u>
Cytogenetics:	8 C1