

## Product datasheet for **MG214240**

### **Ythdc1 (NM\_177680) Mouse Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	Ythdc1 (NM_177680) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Ythdc1
Synonyms:	A730098D12Rik; C80342; mKIAA1966
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide Sequence:**

>MG214240 representing NM\_177680  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGCATCGCC**

ATGGCGGCCGACAGCCGGGAGGAGAAAGATGGGAACTTAATGTTTTGGATGATATTTGACTGAAGTAC  
 CAGAACAGGATGATGAATTGTATAATCCAGAGAGTGAACAAGATAAAAAAGGAAAAAGGATCAAAAAG  
 AAAAAAGTGAAGAAATGGAGTCTACTGACACCAAGCGACAAAAACCTTCCATCCATTCAAGACAATTGATT  
 TCTAAGCCACTAAGCTCATCTGTAAGCAATAATAAAGAATAGTTAGTACAAAAGGAAAGTCGGTTACAG  
 AATATAAAAAATGAGGAATATCAAAGATCTGAAAGGAACAAGCGTCTAGATGCTGATCGAAAAATTCGTCT  
 GTCAAGCAGTTCCCTCTAGAGAACCTTACAAGAGTCAACCAGAAAAAACTTGTCTACGGAAAAAGGATTCT  
 GAAAGAAGGGCCAAGTCTCTACACCAGATGGTCTGAGAGAATTGGGCTGAAGTTGATAGACGTGCAA  
 GCAGATCCAGCCAGTCTTCAAAGGAAGAGGTGAACTCTGAAGACTATGGCTCTGACCACGAGACAGGAAG  
 CAGTGGTCTTCTGAACAGGGCAACAACACTGAGAATGAGGAGGAAGGAGGGGAGGAAGATGTAGAGGAA  
 GATGAAGAGGTAGATGAAGATGCAGAAGATGATGAGGAAGTGGATGAAGATGCAGAGGAGGAGGAGGAAG  
 AGGATGAAGAGGAAGAAGAGGATGAGGACGAGGATGAAGAGGAAGAAGAATATGAACAGGATGAGAGAGA  
 TCAGAAGGAAGAAGGAAATGATTATGACACCCGAAGTGAAGGCGAGTATTCTGGTCCGAGTCTGTTTCC  
 TTCACAGATGGATCTGTGATCTGGTTCAGGAACAGATGGATCAGATGAGAAAAAGGAAAGGAAAGGA  
 GAGCTCGAGGCATATCACCCATTGTCTTTGATAGAAGTGGCAGTCTGCATCAGAGTCATATGCAGGTTT  
 AGAAAAAGCATGAGAAATTATCATCTCCGTTCTGCTGTCCGAAAAGATCAAACCAGTAACTCAAA  
 TATGTCCTTCAGGATGCAAGATTTTCTCATAAAGAGTAATAACCATGAGAATGTGTCTCTTGCCAAAG  
 CAAAGGGTGTATGGTCCACATTGCTGTAATGAGAAGAAATTAATCTTGCAATTTAGATCTGCAAGGAG  
 TGTATATTAATATTTTCTGTGAGAGAAAGTGGAAAAATTTCAAGGTTTTGCCAGATTGTCTTCAGAAATCA  
 CACCATGGCGGATCTCCTATACACTGGGTGCTTCCAGCAGGAATGAGTGCTAAAAATGCTTGAGGTTT  
 TAAAAATAGACTGGATTTGCAGGCGTGAATTACCTTTACTAAATCAGCTCATCTACCAATCCCTGGAA  
 TGAACATAAGCCAGTAAAGATTGGACGTGATGGACAGGAAATGAACTTGAATGTGGAAGTCAAGTCTTGT  
 CTTCTGTTTCCCCTGATGAAAGTATTGACTTGTATCAGCTCATTCAAAAATGCGTCACAAGAGAAGAA  
 TGCATTCTCAGCCTCGATCAAGAGGACGTCCATCCCGTCGAGAACCAGTCCGGGATGTGGGAAGGCGTCCG  
 ACCAGAAGATTATGATATTCATAACAGCAGAAAAGAAACCAAGGATTGACTATCCCCCTGAGTTTACCAG  
 AGACCAGGATATGTAAGGATCCCCGATACCAGGAAGTGGACAGTTTTACAAATCTTATCCCAACAGAC  
 GATTTTCAGGAGTTCGCCGAGATGTGTTTTAAATGGGTCTACAATGATTATGTGAGGGAATTTATAA  
 CATGGGACCACCGCTCCTTGGCAAGGAATGCCTCCTTACCCGGGAATAGAACAACCTCCACACCATCCT  
 TACTACCAGCACCATGCCCCACCTCCTCAAGCCACCCCTTACTCAGGACACCATCCGGTACCACATG  
 AAGCAAGATACAGAGATAAACGAGTACATGACTATGATATGAGGGTTCGATGATTTCTTCCGCCGACGCA  
 AGCTGTGGTCAGTGGTCCGAGAAGTAGACCTCGAGAAAGAGATCGGGAGCGAGAGCGAGACCCCCAGA  
 GATAACAGAAGAGATAGAGAGCGAGACAGAGGTCTGTGATCGAGAAAGAGAGAGAGAAAGAATATGTGATC  
 GGGACAGAGACCGAGGGGAGAGAGGTCGTTATCGAAGA

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA

**Protein Sequence:** >MG214240 representing NM\_177680  
 Red=Cloning site Green=Tags(s)


MAADSREEKDGELNVLDLILTEVPEQDDEL YNPESEQDKNEKKGSKRKSERMESTDTRKQKPSIHSRQLI  
 SKPLSSSVSNKRIVSTKGSVTEYKNEEYQRSENRKRLDADRIRLSSSSSREPYKSQPEKTCRLKRDS  
 ERRAKSPTPDGSEIRIGLEVDRRASRSSQSKEEVNSDYGSDHETGSSGSSEQGNNTENEEEGEEDVVEE  
 DEEVDEDAEDDEVEDAEAAAAAAAAEEEEDE  
 FTDGSVRSVSGSGTDGSDKKEKRRARGISPIVFDKSGSSASESYAGSEKKHEKLSSSVRAVRKDQTSKLLK  
 YVLQDARFFLIKSNHENVSLAKAGVWSTLPVNEKLNLAFRSARSVILIFSVRESGKFGQFARLSSSES  
 HHGSPHWWLPAGSAKMLGGVFKIDWICRRELPTKSAHLNPNWNEHKPVKIGRDGQEIIELECGTQLC  
 LLFPPDESIDL YQLIHKMRHKRRMHSQPRSRGRPSRREPVRDVGRRRPEDYDIHNSRKKPRIDYPPPEFHQ  
 RPGYVKDPRYQEVDSFTNLIPNRRFSGVRRDVF LNGSYNDYVREFHNMGPPPPWQGMPPYPGIEQPHPHP  
 YYQHHAPPPQAHPPYSGHHPVPHEARYRDKRVHDYDMRVDDFLRRTQAVVSGRRSRPRERDRERDRPR  
 DNRRDRERDRGRDRERERERICDRDRDRGERGRYRR


TRTRPLE - GFP Tag - V

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**

Cloning sites used for ORF Shutting:



**Plasmid Map:**


**ACCN:** NM\_177680

**ORF Size:** 2208 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_177680.4](#)

**RefSeq Size:** 3034 bp

**RefSeq ORF:** 2211 bp

**Locus ID:** 231386

UniProt ID: [E9Q5K9](#)

Cytogenetics: 5 E1

**Gene Summary:** Regulator of alternative splicing that specifically recognizes and binds N6-methyladenosine (m6A)-containing RNAs (By similarity). M6A is a modification present at internal sites of mRNAs and some non-coding RNAs and plays a role in the efficiency of mRNA splicing, processing and stability (By similarity). Acts as a key regulator of exon-inclusion or exon-skipping during alternative splicing via interaction with mRNA splicing factors SRSF3 and SRSF10 (By similarity). Specifically binds m6A-containing mRNAs and promotes recruitment of SRSF3 to its mRNA-binding elements adjacent to m6A sites, leading to exon-inclusion during alternative splicing (By similarity). In contrast, interaction with SRSF3 prevents interaction with SRSF10, a splicing factor that promotes exon skipping: this prevents SRSF10 from binding to its mRNA-binding sites close to m6A-containing regions, leading to inhibit exon skipping during alternative splicing (By similarity). May also regulate alternative splice site selection (By similarity). Also involved in nuclear export of m6A-containing mRNAs via interaction with SRSF3: interaction with SRSF3 facilitates m6A-containing mRNA-binding to both SRSF3 and NXF1, promoting mRNA nuclear export (By similarity). Also recognizes and binds m6A on other RNA molecules (By similarity). Involved in random X inactivation mediated by Xist RNA: recognizes and binds m6A-containing Xist and promotes transcription repression activity of Xist (By similarity). Involved in S-adenosyl-L-methionine homeostasis by regulating expression of MAT2A transcripts, probably by binding m6A-containing MAT2A mRNAs (PubMed:29262316).[UniProtKB/Swiss-Prot Function]