

Product datasheet for **MG214357**

Olfr1385 (NM_001011805) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	Olfr1385
Synonyms:	MOR256-41P; MOR256-42P; Olfr1549-ps1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >MG214357 representing NM_001011805
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**


ATGGATAGTTTCAATGCCACGTTAGAAGAGAGATTCTTTTTGGTGGGCTTTTTAGATTGGCCTCAGCTGG
AACTCATCCTTTTTGTTTTATTTCATTTTTCTACTCCCTAACTATCTTTGGCAACACCGCCATCATCGC
TCTCTCCAAATGGACCTTCGACTGCACACCCCATGTACTACTTCTCTCCACCTCTCCTTCTGGAC
CTCTGCTACACCACCAGCACTGTGCCACAGCTTCTATAAACTTTCATGGACTTGACAGGACCATCAGCT
ATGGAGGGTGTGTGGCCAGCTTTTCATATCTCTTGCCTGGGATCCACTGAGTGTGTGCTCTTAGTGGT
GATGGCTTTTGACCGCTATGCTGCTGTGTGTCCTCCACTGCACTACATGACCATCATGCACCCCTTCTC
TGTCAGGCATTGGCTCTTGCTTCATGGGTAGGAGGATTCTTGAACCTCTGATTCAAACAGGCCTCATGA
TGGCCATGCCACTCTGTGGACATCGACTGAATCACTTCTTCTGTGAGATGCCTGTCTTCTCAAGTTGGC
CTGCCAGGACACAGGAGGAACAGAGCCAAGATGTTTGTGGCCAGAGCCATAATCTTGGTCTTCCCTGCA
ACATTAATCTTGGCTCCTATGGACACATTGCTAAAGCACTAAAGGTCAAGTCAACAGCTGGACGCA
GAAAAGCTTTTGGGACATGTGGTCCCACCTTCTGGTGGTTTCCCTGTTTTATGGCTCAGCCATTACAC
ATACTTGCAACCAAGAGCAGTTATTCTGAGAGTGATGAAAATTTGTGGCCCTTTTTATACTATTGTC
ACCCCATGCTCAACCTCTGATCTATACCCTGAGGAACAAGGATGTGAAGGGGCTCTGTGGAAGGTGC
TAGGGAGAGGCACAGAC

ACCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001011805.1</u> , <u>NP_001011805.1</u>
RefSeq Size:	930 bp
RefSeq ORF:	930 bp
Locus ID:	258027
Cytogenetics:	11 B1.2
Gene Summary:	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]

Product images:



Circular map for MG214357