

Product datasheet for **MG214420**

Trif1 (NM_001039488) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Trif1 (NM_001039488) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Trif1
Synonyms:	2010012G17Rik; 4933421E11Rik; A1450568
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MG214420 representing NM_001039488
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGTCTAATAGTCTCCAGAGCGTCATTCTGAAAACCCAGAGGAAAAGTCAGGCTCGCGCTGTATTTTCAG
 GGTGCATGTACCAAGTAGTTCCTACGATTGGCAGCGATGGAAAAAACTTCTGCAATTACTTCCAATTTT
 TAAATCTTCTGGAAATCTTATCCAGTAGTTCAGTCTCCAGTCATGTCTCATGGTTTGAAAAGCGAATACT
 GAAAAACAGTTCAAGTTACTTTTCAAACCTCAGATTTCCAGCTCTCCACAAGTGCATCTGTTCAATTGC
 CTGTTTTTCAGCCAGCCAATACTACAAAATGTTTTTTTACAGGAGCAATAGATACAACAGGAAAAGATAG
 AGTTACTTCTGTGAGGACTGGAAATTTTACACCACCAGTTTCTAACATTCAGAATCATGGTGTGAAAATT
 CATAAACTCACCAGGCAAACATTTACTATTCTCCCTCAACACAAAATGATTCATCTCATTTTTATATTCA
 ATACCCCAAGTCTCTACCAAATGTGAATCTTCAATTTTGCCTTCTGGGAATCATTTAAAGATTCCAGC
 CCATGCTGAAGTGAATCTGTACTTGCATCATCTTGCCTCCTTTGGTTCAGCAAAAGATACTTGGAACT
 GCAACCACAAGTACCTCAGGAACAGTTGAGGCCTCCCAAATACCAACAGTTGTTTATGTACATCCTGTAA
 ATTCAGTGAAAATTTGTTGTACCAAGAAAACCTCAGACTATTTACCCAAACAGTGACATTCATACCTT
 GCAAAATCCACCAAATGTTGCTACAGAAAACACAATTAAGGTTGGTCAGCATCCTCAAGCTGCTCCAGTG
 AACTCGATTTTTTCAGGAGTATCTACAGCCTGGAATCCCATGTAATTATTCTGTAAAGTCTTCAAATAATG
 TGGCTACAAAGGTTTTTAAACACTTTTGTAGGTAGGAAAAATCTGGGAGATAATACAATAGATACACCTCT
 GTTGAATACCAACAGTTCTGGCAGGACACACTCCGTAAGTGAACCTATTAAGATAATGCTTTGATAATG
 TTTAATGGGAAAGTCTGGTTGAATGAAAAGGGGACATGCGGTCTGCCATCAAAAATGACCAGCAGAATT
 CTGTTTTCTTGATATCCACTAAAAGATTCATCACAGCTAGTAAGTTCAAGTATAGTCACAGAAATATC
 CAGAGAGATTTTTAAATAGTGTTTTGGTTAAAAGTAAATCTTTTTCAGTTGAAGACAAAATCACTTTCAAAT
 TCCAGCTTGCTTCCATGGCCAATCTAAGAGCAGAGAAGAAATGAAAAAGTGGAGAGACCATCTTTTTCTG
 TCACAAATCCACATACTATGAACCAGTCCACTCACTGCCTTAAACAGAGTAAGACTGTATTTATAAACCC
 AGTTTTTCCAGATGGATTTAGGACAGGACAGAATGCCCCAGAAAAGGAAATCTCGTCCAGAACATAGAG
 AAAATATGTTCTCTGTTGATGCAGCAACTGTTACTTCTCAACAGTGTGTTTTTCAGAGACCAAGAATCAC
 AGACGCAGTATGAGATGGCATCAATAGTAAAAAGGAAATTCAGAAAAAGGAAATAACAAGAAATATTC
 TCAAGGAAGCCATACTAATATTAAGGCATCATGTCTGAAGAATGATGCTGAATTTAAAAAGCTATTTGGT
 CTCACTAAAGATTTGAGAGTATGTCTTACTCGAATTCCTGACCATCTGAGCTCTGAAAAAGTTTTAATT
 CCTTTAACAGTTTGATGAAAAGTAGTCTTACAAGATGCAAACATTGTGGTGAAGAAAGAGAAAAAAA
 ACAGAGTTTCTCTAAGAAAAGAAAAGCGGAAACCATGAAGATGGGTAACAAAAAAGATAAAAATTGAG
 AATGCTGATGACACGGTCATGAGCATCATGAATGGAACCTGATGTTGCTAGTTCCCAACCTCTCAGTAGTA
 TTTTACCAACATCAGACATATCAACAATAACATTGTCACAAGCCATAGTACAACCAGGGAAGACAAGAG
 AACTGAGGCAGAACACTGTTCTCATGAGAAGCAGGAAAAAGGCACATTAAGTTCAAGCACATCTTTTGAA
 CAAAGTACTTTTTTAAATAAAAACCTTTATGGAAGATATTTTTCCAGTGACACCACCAGAAGTAAAGAAA
 CCATCAGAGATGAAAAAATAAGAAGACTTAAGCAAATTTAAGAGAAAAAAGAAGCAGCTCTTGAAGAATT
 GCGTAAGAAGATGTACCAAAAACAA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG214420 representing NM_001039488
Red=Cloning site Green=Tags(s)

```

MSNSLQSVILKTAEEKSGSRCISGCMYQVVPTIGSDGKKLLQLLPISKSSGNLIPVVQSPVMSHGLKANT
EKPVQVTFQTQISSSSTSASVQLPVFQPANTTKCFFTGAI DTTGKDRVTSVRTGNFPPVSNIQNHGVKI
HKLTRQFTTIPPSTQNDSSHFIFNTPSLLPNVNSSILPSGNHLKIPAHAEVKSVLASSLPPLVQQKILGT
ATTSTSGTVEASQIPTVVVYVHPVNSVKFVVTKKTQTIYPKPVTFNTLQIPPVATETQLKGGQHPQAAPV
NSIFQEYLQPGIPCIIPVKSSNNVATKVLNTFVGRKNLGDNTIDTPLLNTNSSGRTHSVSEPIKDNALIM
FNGKVVWLNKGTGCLPSKIDQQNSVSSDIPLKDSSQLVSSSIVTEISREILNSVLVKSFSQFKTKLSLSN
SQLASMANLRAEKNEKVERPSFSVTNPHTMNQSTHCLKQSKTVFINPVFPDGFRTGQNAPRKNLVQNIIE
KICSSVDAATVTSQQCVFRDQESQTQYEMASIVKKEIQEKGNKKYSQGSHTNIKASCLKNDAEFKLLFG
LTKDLRVCLTRIPDHLSSGKSFNSFNLMKSSSYKDANI VVKKEEKQSFSKRKAETMKMGNTKKIKIE
NADDTVMSIMNGTDVASSQPLSSILPTSDISQHNIVTSHSTTREDKRTEAEHCSHEKQEKGLSSSTSFE
QSTFLNKNFMEDIFPVTPELEETIRDEKIRRLKQILREKEAALEELRKKMYQKQ
    
```

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001039488

ORF Size: 2265 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001039488.1](#), [NP_001034577.1](#)

RefSeq Size: 3370 bp

RefSeq ORF: 2268 bp


Locus ID: 321000

UniProt ID: [Q8CDD9](#)

Cytogenetics: 3 F2.3

Gene Summary: Together with SMCHD1, involved in chromosome X inactivation in females by promoting the compaction of heterochromatin. Also able to repress the ligand-induced transcriptional activity of retinoic acid receptor alpha (RARA), possibly through direct recruitment of histone deacetylases.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MG214420