

## Product datasheet for **MG214505**

### Olf503 (NM\_001011527) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Tag:** TurboGFP  
**Symbol:** Olf503  
**Synonyms:** MOR34-8P; MOR34-9; MOR34-12; Olf548  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)

**ORF Nucleotide Sequence:** >MG214505 representing NM\_001011527  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCCGCGATCGCC

ATGGTCATGTCAGCGAGGAATAACCCAGATCTGACCCAGCTTCCTTTGTTCTGAATGGGATCCCAGGCC  
 TGGAAACACATGCACATCTGGATCTCCTTCCCTTTCTGTTCCATGTATGCAGTGGCTATGATGGGAACTG  
 TGGACTCCTCTATCTCATCTTCTTTGAGGATTCCTTACACAGACCCATGTAATACTTTTAGCAATGCTT  
 TCTCTAACAGACCTTGTCATGTGCTCTAGTACAATCCCAAAGCCCTGTGCATCTTCTGGTTTCATCTCA  
 AAGAAATTGGATTTGACGACTGTCTGTACAGATGTTCTTCATCCATACCTTCATGGGGATGGAGTCTGG  
 AGTGCTCATGCTCATGGCTCTGGACCACTATGTAGCCATTTGCTACCCTCTGCATTACTCCACCATTCTC  
 ACCAATCCTATCATTGCCAAGATTGGGTTAGCCACCTTCTGAGGGCGTTCTGCTAATTATTCCATTCA  
 TATTTCTCTCCAAGTGCCTACCCTACTTTGAGGCAATATAATAAACCATACTACTGTGACCACATGTC  
 TGAGCCAAGTTGCCCTGTGGCAATGTCAAGGTGAATGCCATTTATGGTCTAATGGTTGCTCTTTGATT  
 GGGGTTTTGACATCCTGTGTATCACAATCTCCTATATCATGATTCTGAGGGCAGTGGTCAGCTTATCAT  
 CAGCAGATGCTAGGCAGAAGCCCTCAGCACCTGCACTGCCACATCTGTGCCATTGTTTTCTCCTATAG  
 CCCAACCTTCTTGTCTTCTTTCCACCCTTTGGGGGACACACAATCCCTCCATCTTGCCACATCATT  
 GTGGCTAATTTATCTGCTTTTGCCTCCCACTATGAACCTGTTGCTATGGAGTGAACCAACCAACAGAG  
 TACAAGACTGTGCATAAGGATTTTTTCAGAGTCTAAGGATTCCAAAGCTCACGGTATA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



**Protein Sequence:** >MG214505 representing NM\_001011527  
 Red=Cloning site Green=Tags(s)

MVMSARNNPDLTPASFVLNGIPGLEHMHIIWISFPFCSMYAVAMMGNCGLLYLIFEDSLHRPMYYFLAML  
 SLTDLVMCSSTIPKALCIFWFHLKEIGFDDCLVQMFFIHTFMGMEGVLMALDHYVAICYPLHYSTIL  
 TNPIIAKIGLATFLRGVLLIIPFIFLSKCLPYFRGNIINHTYCDHMSVAKLPCGNVKVNAIYGLMVALLI  
 GGFIDILCITISYIMILRAVVSLSADARQKAFSTCTAHICAIVFSYSPTFLSFFSHHFGGHTIPPSCHII  
 VANIYLLLPTMNPVVYGVKTKQIQDCVIRIFSESKSKAHGI

TRTRPLE - GFP Tag - V

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**

Cloning sites used for ORF Shuttling:



**ACCN:** NM\_001011527

**ORF Size:** 969 bp


**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>Note:</b>	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
<b>RefSeq:</b>	<u><a href="#">NM_001011527.1</a></u> , <u><a href="#">NP_001011527.1</a></u>
<b>RefSeq Size:</b>	972 bp
<b>RefSeq ORF:</b>	972 bp
<b>Locus ID:</b>	259143
<b>Cytogenetics:</b>	7 E3
<b>Gene Summary:</b>	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]

Product images:



Circular map for MG214505